



Comfort-Tech E2601




Competition Comparison

Comfort-Tech

Jay Basic

ROHO Mosaic



Avg Pressure at IT's (mmHg)

62.53

89.06

65.20

Relative Humidity (%)

3.10

31.08

67.30

Postural Alignment

Yes

No

No

Testing was completed by Central Coast Testing, 290 Ridge Road, Ben Lomond, California 95005 Test cushions were acquired through commercial retail outlets. Where applicable the cushions were adjusted for fit per manufacturers instructions.