

# CUTIMED® SORBION®

## SUPER-ABSORBENT WOUND DRESSING

Cutimed® Sorbion® dressings have the excellent ability of quickly absorbing large amounts of exudate and locking it within their inner core, even under compression, therefore reducing the frequency of dressing changes. Cutimed® Sorbion® dressings also have the potential to remove chronicity-causing factors and thereby support wound healing; they provide gentle debridement, reduce the concentration and activity of MMPs and immobilize micro-organisms. The unique fiber matrix offers stability even when highly saturated, thus preventing the formation of lumps and pressure points and ensuring maximum wound contact.

### INDICATIONS

- For all moderately to highly exuding wounds such as venous, diabetic or pressure ulcers

### FEATURES AND BENEFITS

- High absorption capacity, long dressing wear time
- Fewer dressing changes (up to 4 days wear time)
- Securely locks in exudate, even under compression
- Removes wound debris and locks in wound germs
- Maintains a moist wound environment by gel formation
- No binders or adhesives; well tolerated
- No risk of leakage, high tear strength
- Able to apply on both sides
- Less odor formation
- Sterile
- Latex Free



Description		Dimensions	Code	Box Contents	Selling Unit
<b>Cutimed® Sorbion® Sachet S</b> The classic form suited for many wounds types and areas of application 		7.5 cm x 7.5 cm	7323200	10 dressings	Box
		10 cm x 10 cm	7323206	10 dressings	Box
		5 cm x 12 cm	7323203	10 dressings	Box
		10 cm x 20 cm	7323209	10 dressings	Box
		15 cm x 15 cm	7323212	10 dressings	Box
		20 cm x 20 cm	7323215	10 dressings	Box
		20 cm x 30 cm	7323218	10 dressings	Box
<b>Cutimed® Sorbion® Sachet Multi Star</b> Flexible form for applications on toes, heel, foot, sacrum, chest or wound cavities 		8 cm	7323701	10 dressings	Box
		14 cm	7323703	10 dressings	Box
<b>Cutimed® Sorbion® Sachet XL</b> For large difficult-to-dress wounds on lower leg, sacral region, thorax, thigh, lower part of abdomen 		25 cm x 45 cm	7324000	10 dressings	Box
<b>Cutimed® Sorbion® Sachet Border</b> With hypoallergene adhesive border 		10 cm x 10 cm	7323600	10 dressings	Box
		15 cm x 15 cm	7323602	10 dressings	Box
		15 cm x 25 cm	7323604	10 dressings	Box
<b>Cutimed® Sorbion® Sana</b> With particularly atraumatic surface for pain-sensitive patients 		8.5 cm x 8.5 cm	7323300	10 dressings	Box
		12 cm x 12 cm	7323303	10 dressings	Box
		12 cm x 22 cm	7323306	10 dressings	Box
		22 cm x 22 cm	7323309	10 dressings	Box
<b>Cutimed® Sorbion® Sana Multi Star</b> With particularly atraumatic surface for pain-sensitive patients 		11 cm	7323801	10 dressings	Box
		17 cm	7323803	10 dressings	Box

\* See legend on page 18