

6855

High Capacity Waist-High Stand-On Scale

Designed for obese and unsteady patients, this 600 lb capacity scale features a balanced, lightweight design. It's easily tipped back for transportation on its hard, large diameter, rubber wheels yet narrow enough to fit through 30" doors. The battery-powered, LCD display rotates for easy readability. The low-profile platform, convenient handrail and skid-proof mat add to patient assurance.



SPECIFICATIONS

Model:	6855
Capacity:	600 lb x 0.2 lb / 270 kg x 0.1 kg
Platform Size:	18"W x 14"D / 46 cm x 26 cm
Display:	<ul style="list-style-type: none">• 0.75" high digits, high-contrast LCD• Body Mass Index calculation• LOCK button holds weight after weighing
Power:	6 "AA" size batteries (not included) or optional AC adapter.
Shipping Weight:	70 lb / 32 kg

Options include height rod and AC adapter.



Large LCD and Body Mass Index calculation