


New Eclypse® Super Absorbent Dressings Available Now...

Eclypse Foot®

Super absorbent foot wrap dressing



Eclypse Foot® is a super absorbent dressing, designed to be applied securely around the foot area to manage wound exudate.

The dressing utilizes the high fluid handling capabilities of the Eclypse® range in a pre-shaped dressing designed to fit comfortably and effectively around the foot. It is easy and simple to apply and offers enhanced patient comfort.




Eclypse Boot®

Super absorbent leg wrap dressing



Additional sizes now available



Product	Overall Size	Calf circumference	Leg length	Stock code	Dressings per box	HCPCS code
Eclypse Foot®	13in x 18.9in	n/a	n/a	CR4219	5	A6198*
Eclypse Boot® Small	26.1in x 18.5in	≤ 15in	8in (Mid-calf to heel curve)	CR4237	5	A6198*
Eclypse Boot®	27.6in x 23.6in	≤ 21in	14-19in (posterior knee to heel curve)	CR4014	5	A6198
Eclypse Boot® XL	27.9in x 31.5in	≤ 24in	≥ 19in (posterior knee to heel curve)	CR4229	5	A6198*

*pending approval
≤ (less than or equal to)
≥ (more than or equal to)

Disclaimer

Advancis Medical USA LLC has made every effort to ensure the information in this leaflet is accurate. As HCPCS codes and associated coverage policies are frequently revised, we recommend you consult with a reimbursement specialist and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

The listing of a HCPCS code does not guarantee coverage or payment.

Advancis
Medical
www.advancismedical.com

New Eclypse® Super Absorbent Dressings Available Now...

Eclypse Border® and Eclypse Border Oval®

Super absorbent dressing with silicone contact layer and border



Eclypse Border® and Eclypse Border Oval® are super absorbent dressings with a soft silicone contact layer and border with added hypo-allergenic adhesive on the border to ensure the dressing remains in place.

The dressing has a rapid wicking face and a highly absorbent moisture locking system with an additional wicking layer to further maintain the integrity of the dressing and help to distribute exudate across the entire dressing.

Dressing features

- ✓ Increased absorbency capacity
- ✓ Atraumatic and hypo-allergenic
- ✓ Breathable, fluid-repellent backing
- ✓ Provides moist wound healing environment
- ✓ Additional wicking and moisture locking central core
- ✓ Unique large adhesive border
- ✓ Unique oval shape

Product	Size	Stock code	Dressings per box	HCPCS code
Eclypse Border®	5.9in x 5.9in	CR4296	10	A6197*
Eclypse Border®	7.9in x 11.8in	CR4297	10	A6198*
Eclypse Border Oval®	3.9in x 7.9in	CR4300	10	A6197*
Eclypse Border Oval®	5.9in x 7.9in	CR4299	10	A6197*
Eclypse Border Oval®	7.9in x 11.8in	CR4298	10	A6197*

*pending approval

