

Specifications and Dimensions

Item	Overall Height	Overall Width	Overall Depth	Width Between Arms	Arm Height From Floor	Arm Height From Seat	Seat Height	Seat Width	Seat Depth	Depth Fully Reclined In Trendelenburg	Width With Both Side Tables in Upright Position	Shipping Weight	Max. Weight Capacity
FR566G	50"	31"	43"	21"	27 ½"	8"	22"	20 ½"	20"	75"	51"	145 lb	350 lb
FR566DG	51 ½"	33"	38"	23"	28"	8"	22"	20 ½"	19"	76"	53"	135 lb	350 lb
FR577RG	47"	31"	35"	22"	25 ½"	8"	22"	21 ½"	21"	74"	54"	127 lb	290 lb
FR587W	50"	33"	36"	23"	26 ½"	8"	21 ½"	23 ½"	20 ½"	76"	53"	138 ½ lb	450 lb
FR588W	55"	38"	59"	25 ¼"	32"	8"	25 ½"	25 ¾"	24"	71"	58"	308 lb	700 lb
FR3302	44"	31"	38"	23"	19 ½" - 25 ½"	5"	19"	21 ½"	20"	-	-	109 lb	250 lb
574G	47"	25"	41"	20"	26"	9"	21"	19"	19"	-	-	110 lb	250 lb
565G	48"	28"	36"	19"	28"	8 ½"	20"	19"	19 ½"	-	-	80 lb	250 lb
565WG	48"	32"	37"	24"	28"	8 ½"	20"	24"	19 ½"	-	-	91 lb	400 lb
565TG	48"	28"	36"	19 ½"	23 ½" - 29 ½"	0" - 8 ½"	22"	19"	19 ½"	-	-	90 lb	250 lb
565DG	48"	28"	36"	19 ½"	23 ½" - 29 ½"	0" - 8 ½"	22"	19"	19 ½"	-	-	85 lb	250 lb
GF4404	42"	25 ¾"	26"	18"	31"	6"	26"	18"	18"	-	-	52 lb	350 lb

For the most up-to-date information and product specifications, please visit our website www.grahamfield.com.



GF Health Products, Inc. Atlanta, Georgia 30360
800.347.5678 www.grahamfield.com