


## What is a booster pad?

Booster pads have a flow-through design that allows the pad to fill to capacity first and then pass additional fluid to the host (primary) garment. Booster pads are placed inside any disposable undergarment with a moisture barrier, boosting absorbency to extend wear-time, reduce changes and help contain bowel incontinence. Strategically placed over the hip helps prevent leaking for side sleepers.

For additional placement ideas, see page 20.



ADDED

PROTECTION



TOP SELLER

## TopLiner® Booster Pads

Thin design allows for placement in target areas where additional absorption is needed.

Product Features:



ITEM #	SIZE	LENGTH/WIDTH	CAPACITY	PACKAGING	HCPC
2072	Mini	10.5" x 2.75"	5.7 oz/169 mL	8 bags of 25 (200)	T4535
2070	Regular	14" x 4"	10.9 oz/322 mL	8 bags of 25 (200)	T4535
2060	Super	15" x 4.25"	14.2 oz/420 mL	8 bags of 25 (200)	T4535

Absorbency: 14.2 oz/420 mL, size Super



ADDED

PROTECTION



TOP SELLER

## TopLiner® Booster Contours

Wide hourglass shape helps contain bowel. Soft texture and no adhesive strip makes cleanup easier.

Product Features:



ITEM #	SIZE	LENGTH/WIDTH	CAPACITY	PACKAGING	HCPC
3096	Contour	21.5" x 13.5"	13.6 oz/403 mL	10 bags of 12 (120)	T4535
3097	Super-Plus*	32" x 14"	27.5 oz/813 mL	8 bags of 12 (96)	T4535

Absorbency: 13.6 oz/403 mL, size Contour

\* Ideal for use with sizes XL and greater.

