

### 3-Position Recliner



FIGURE A



FIGURE B

- 3-positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels “pop off” for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy-to-move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy-to-clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117
- Standard

#### PRODUCTS

#### Description

D574-BR	Blue Ridge, 1/cs
D574-J	Jade, 1/cs
D574-R	Rosewood, 1/cs

#### SPECIFICATIONS

#### Description

Overall Dimensions	47" (H) X 25" (W) X 41" (D)
Width	20" between arms
Seat to Floor Arm Height	9" - 26"
Seat Dimensions	21" (H) x 19" (W) x 19" (D)
Carton Shipping Weight	110 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts

#### PARTS

#### Description

D574P-1058BR	Footrest Blue Ridge, 1/bx
D574P-1058J	Footrest Jade, 1/bx
D574P-1058R	Footrest Rose, 1/ea
D574P-1055BR	Legrest Blue Ridge, 1/bx
D574P-1055J	Legrest Jade, 1/bx
D574P-1055R	Legrest Rose, 1/bx
D574SEAT-BR	Seat Blue Ridge, 1/bx
D574SEAT-J	Seat Jade, 1/bx
D574SEAT-R	Seat Rose, 1/bx
D574P-SEATLEG-BR	Seat with Legrest Blue Ridge, 1/ea
D574P-SEATLEG-J	Seat with Legrest Jade, 1/ea

PARTS	Description
D574P-SEATLEG-R	Seat with Legrest Rose, 1/ea
D574BACK-BR	Back Blue Ridge, 1/bx
D574BACK-J	Back Jade, 1/bx
D574BACK-R	Back Rose, 1/bx
D574P-1053	Seat Frame, 1/ea
D574P-1048	Headrest Adjuster, 1/ea
D574P-1004	Push Handle, 1/ea
D574P-1050	Back Lock Bar, 1/ea
D574P-1070BR	Armrest Blue Ridge, 1/ea
D574P-1070J	Armrest Jade, 1/ea
D574P-1070R	Armrest Rose, 1/ea
D574P-1020	Side Panels, 1/bx
D574P-1025B	Rear Caster, 1/bx
D574P-1068	Reclining Frame, 1/ea
D574P-1027	Front Caster, 1/bx
D574P-1038	Tray Table, 1/cv
D574P-1064	Table Slide Tube, 1/ea
D574P-1036	Table Fitting, 1/ea
D574P-BTN	Release Button, 1/ea
D574P-1037	Table Washer, 1/ea
D574P-1035	Reclining Mechanism, 1/bx
D574P-1060	Footrest Support Bar & Hdwr, 1/ea
D574P-1061	Hardware Only, 1/ea
D574P-1012	Protective Panels
D574-BBN	Back Release Button, 1/ea
D574P-1003	Hardware Set, 1/bx
D574EW-1003	Lower Push Handle, 1/ea
D574EW-1004	Hardware, 1/ea
D574EW-1010BR	Footrest Blue Ridge, 1/ea
D574EW-1010R	Footrest Rose, 1/ea
D574EW-1011	Support Bar & Hdwr, 1/ea
D574EW-1013	Seat Frame, 1/ea
D574EW-1017N	Frame New Style, 1/ea
D574EW-1035	Reclining Mechanism, 1/ea
D574EW-1050	Locking Bar, 1/ea
D574EW-1065	Tray, 1/ea
D574EW-1070BR	Armrest Old Style Blue Ridge, 1/ea
D574EW-1070J	Armrest Old Style Jade, 1/ea
D574EW-1070R	Armrest Old Style Rose, 1/ea
D574EWBACK-BR	Back Blue Ridge, 1/bx
D574EWBACK-R	Back Rose, 1/bx
D574EWSEAT-BR	Seat Blue Ridge, 1/bx
D574EWSEAT-J	Seat Jade, 1/ea
D574EWSEAT-R	Seat Rose, 1/bx
D574P-1048-R	HeadrestAdjBracket D574/D574EW Rght
D577P-1012	Protecting Panel, 1/ea
D574EW-1010J	Footrest Jade, 1/ea
D574EWBACK-J	Back Jade, 1/ea