



COVIDIEN

positive results for life™

Kangaroo ePump™ Enteral Feeding Pump

FEATURES

- Attitude independent
- "Smart" pump set technology
- "ATM-style" ease of operation
- Feeding history up to 72 hours
- DEHP-free pump sets
- Anti-free-flow pump sets
- Programmable flushing capabilities

REVOLUTIONIZING THE DELIVERY OF ENTERAL FEEDING BEGINS WITH e

When patients require nutrition through continuous feeding, intermittent feeding or feeding and flushing, the Kangaroo ePump™ feeding pump delivers in one compact, easy-to-use device.

WHAT DOES e MEAN TO YOU?

e means easy

Easy step-by-step prompts, intuitive user interface and simplified set loading.

e means everybody

Appropriate for virtually all patients, from infant through geriatric years.

e means everywhere

Easily transportable, with battery pack for up to 15 hours of backup power.

e means enteral

Continuous or bolus enteral feeding and flushing capabilities, accurately and reliably, all from one pump.



Easy-to-use pump set



Proximal Spike Set



1000 mL feed and flush pump set



1000 mL feed only pump set



Kangaroo™ safety screw spike

www.kangarooepump.com



FEATURES

- Volume to be Delivered –
1 – 3000 mL in 1 mL/hr increments
- Feed Rate –
1 – 400 mL/hr in 1 mL increments
- Flush Rate –
10 – 500 in 1 mL increments
- Flush Interval –
1 – 24 hrs in 1 hr increments
- 15 hour battery life at 125 mL/hr
- Weight –
2.34 lbs without pole clamp
2.69 lbs with pole clamp
- Multilingual – 15 languages
- Occlusion pressure 1 – 15 psi
- Alarms –8 options

ORDERING INFORMATION

Code Number	Description	Case Quantity
382400	Kangaroo ePump with Pole Clamp	1
382491	Power Adapter	1
382492	Pole Clamp	1
F010484WT	Battery	1

PUMP SETS

774655	Proximal Spike Set	30
772055	500 mL Pump Set	30
773656	1000 mL Pump Set	30
776255	600 mL Rigid Container Pump Set	30
776250	1200 mL Rigid Container Pump Set	30
774669	Proximal Spike Feed and 1000 mL Flush Set	30
773662	1000 mL Feed Set w/ 1000 mL Flush Set	30
672062	Spike Feed and 500 mL Flush Set	30
674668	500 mL Feed Set with 500 mL Flush Set	30
714650	Sterile Proximal Spike Set	30
712054	Sterile 500 mL Pump Set	30
713650	Sterile 1000 mL Pump Set	30
716154	Sterile 100 mL Burette Set	30

RECERTIFICATION SET

776150	100 mL Burette Recertification Set	5
--------	------------------------------------	---

KANGAROO FEEDING SETS FOR USE WITH NESTLÉ®* SPIKERIGHT™* SYSTEM

775655	Kangaroo Safety Screw Spike for Kangaroo ePump and Kangaroo Joey	30
775662	Kangaroo Safety Screw Spike w/1000 mL Flush Bag for Kangaroo ePump and Kangaroo Joey	30
775762	Safety Screw Spike Adapter Cap	30

©*Nestlé is a registered trademark of the Société des Produits Nestlé S.A., Vevey, Switzerland.
 ™*SpikeRight is a trademark of the Société des Produits Nestlé S.A., Vevey, Switzerland.
 COVIDIEN, COVIDIEN with logo and ™ marked brands are trademarks of Covidien AG or an affiliate.
 © 2009 Covidien. All rights reserved.

H6267-2.5M-0609