

How Safe is your Safety Needle?

Unlike other safety devices on the market, the Magellan™ Safety Needle features a robust shield that encases the entire cannula after activation, which protects downstream workers. Its proprietary integrated locking mechanism prevents deactivation of the device once the safety shield has been locked into place. Every Magellan™ Safety Needle may be activated using any of the following methods:



...thumb



...finger



...flat surface.

Needle Only

Code	Description	Packaging
8881850810	18 x 1"	50/Bx, 500/Cs
8881850815	18 x 1 1/2"	50/Bx, 500/Cs
8881850910	19 x 1"	50/Bx, 500/Cs
8881850915	19 x 1 1/2"	50/Bx, 500/Cs
8881850010	20 x 1"	50/Bx, 500/Cs
8881850015	20 x 1 1/2"	50/Bx, 500/Cs
8881850158	21 x 5/8"	50/Bx, 500/Cs
8881850110	21 x 1"	50/Bx, 500/Cs
8881850115	21 x 1 1/2"	50/Bx, 500/Cs
8881850210	22 x 1"	50/Bx, 500/Cs
8881850215	22 x 1 1/2"	50/Bx, 500/Cs
8881850358	23 x 5/8"	50/Bx, 500/Cs
8881850310	23 x 1"	50/Bx, 500/Cs
8881850558	25 x 5/8"	50/Bx, 500/Cs
8881850510	25 x 1"	50/Bx, 500/Cs

Pre-Attached Combinations

Code	Description	Packaging
1mL Combos		
8881811310	1mL 23 x 1"	50/Bx, 200/Cs
8881811558	1mL 25 x 5/8"	50/Bx, 200/Cs
8881811510	1mL 25 x 1"	50/Bx, 200/Cs
3mL Combos		
8881833810	3mL 18 x 1"	50/Bx, 400/Cs
8881833010	3mL 20 x 1"	50/Bx, 400/Cs
8881833015	3mL 20 x 1 1/2"	50/Bx, 400/Cs
8881833110	3mL 21 x 1"	50/Bx, 400/Cs
8881833115	3mL 21 x 1 1/2"	50/Bx, 400/Cs
8881833210	3mL 22 x 1"	50/Bx, 400/Cs
8881833215	3mL 22 x 1 1/2"	50/Bx, 400/Cs
8881833310	3mL 23 x 1"	50/Bx, 400/Cs
8881833558	3mL 25 x 5/8"	50/Bx, 400/Cs
8881833510	3mL 25 x 1"	50/Bx, 400/Cs

6mL Combos

8881866810	6mL 18 x 1"	50/Bx, 200/Cs
8881866015	6mL 20 x 1 1/2"	50/Bx, 200/Cs
8881866110	6mL 21 x 1"	50/Bx, 200/Cs
8881866115	6mL 21 x 1 1/2"	50/Bx, 200/Cs
8881866215	6mL 22 x 1 1/2"	50/Bx, 200/Cs

12mL Combos

8881822810	12mL 18 x 1"	50/Bx, 200/Cs
8881822015	12mL 20 x 1 1/2"	50/Bx, 200/Cs
8881822110	12mL 21 x 1"	50/Bx, 200/Cs
8881822115	12mL 21 x 1 1/2"	50/Bx, 200/Cs
8881822215	12mL 22 x 1 1/2"	50/Bx, 200/Cs



COVIDIEN

positive results for life™

COVIDIEN and COVIDIEN with Logo and TM marked brands are trademarks of Covidien AG or its affiliate.
© 2007 Covidien AG or its affiliate. All rights reserved. Patents Pending.

H-6069 10M 01/08

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

800-962-9888
508-261-8000

WWW.COVIDIEN.COM



Magellan™ Safety Needle



COVIDIEN

positive results for life™



100% inspection for
needle sharpness to
provide increased
patient comfort.



Magellan™ Safety Needle. The point at which technology & safety meet.

While the Needlestick Safety and Prevention Act has reduced accidental needle sticks, hundreds of thousands of unnecessary exposures and injuries still occur each year. To significantly lower this risk, the solution is clear — the Magellan™ Safety Needle, the point at which safety and technology meet.

- Each Magellan™ safety needle is individually inspected more than 20 times throughout the manufacturing process to ensure optimum quality.
- The unique Magellan™ safety needle design and ultra-durable shield ensure the needle remains protected even under the most extreme circumstances.
- Magellan™ safety needle strikes the perfect balance by providing a device that offers exceptional safety while still remaining affordable.

A Safety Needle Even OSHA would love.

- Easy to use, practical and intuitive one handed activation via multiple methods (finger, thumb or surface) enables clinical flexibility for safe delivery of medication.
- Robust safety shield covers entire needle, providing unparalleled protection after disposal.
- Integrated engineering controls incorporate the safety device as an integral part of the needle hub in conformance with OSHA standards, eliminating the possibility of removal.

*Technology & Engineering.
Safe & Sound.*