



CAREASSIST® BED **ADVANTA™ 2 BED** **VERSACARE® BED**

Feature	CAREASSIST® BED	ADVANTA™ 2 BED	VERSACARE® BED
IntelliDrive® Powered Transport: Allows a single caregiver to easily transport patient in bed	Optional	Optional	Optional
Advanced Three-Mode Bed Exit System (Patient Position Monitor): Bed exit system with three different sensitivities	Optional	Optional	Optional
Emergency CPR: Release lever lowers the head section with dampened descent	+	+	+
Brake Alarms: Alerts caregiver when the bed is plugged in and the brake is not set	+	+	+
Chair Position: Caregivers can place patient in a chair position with the touch of one button; bed remains at a low height	+	+	+
Safety Alerts at Foot End of Bed: Provides measures to monitor protocols for bed status (brakes not set, bed exit, bed not in lowest position, siderails down); lights flash an amber color alert when bed is not set in any of the programmed protocols	Optional	Optional	Optional
Integrated Scale	Optional	Optional	Optional
Battery Backup: Backup for all bed articulation (high/low position, head-of-bed, knee section)	+	+	+
Point-of-Care® Siderail Controls: All caregiver controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs	+	+	+
Line-of-Site® Feature: Head-of-bed and Trendelenburg indicators integrated into siderails	+	+	+
HandsFree® Controls: Foot pedals activate high/low and head up/down so your hands are available for your patient			+
Electric Vascular Foot: Eliminates the need for the nursing staff to manually lift the patient's legs during a cardiovascular event	+	+	
Drop Down Siderails: One-handed single action siderail operation	+	+	+
Boost® One-Button Function: Puts the bed in the ideal position for patient repositioning			Optional
Built-in Bed Extender: Easily adds length to the bed	+	+	+
FlexAfoot™ Mechanism: Allows caregiver to customize via siderail controls the overall length of the bed and mattress, in any bed position (including chair position) to fit the patient			+
Ability to send data to Nurse Call, EMR or other Clinical Applications	Optional	Optional	Optional
Fully Integrated Surface: Surface automatically adjusts to the patient's weight, size, movement and to the head-of-bed articulation, to improve patient comfort			Optional
Advanced Microclimate® Technology: Helps reduce heat and moisture at the patient sleep surface interface			Optional
Status Board Nurse Call Display via smart beds	Optional	Optional	Optional

[Click here to view detailed product descriptions and prices.](#) **Call 800-397-5899 to order today!**



Feature

TOTALCARE® BED

**TOTALCARE®
BARIATRIC PLUS BED**

IntelliDrive® Powered Transport: Allows a single caregiver to easily transport patient in bed	Optional	Optional
Advanced Three Mode Bed Exit System (Patient Position Monitor): Bed exit system with three different sensitivities	Optional	Optional
Emergency CPR: One-step CPR flattens the frame and max inflates the air surface	+	+
Brake Alarms: Alerts caregiver when the bed is plugged in and the brake is not set	+	+
Seat Deflate: Seat and foot deflate offers low bed height for easier exit	Optional	Optional
Turn Assist: Uses the bed's pressure redistribution surface to help turn the patient, making it easier to perform linen changes and wound checks	Optional	Optional
CLRT and Percussion & Vibration: Integrated continuous lateral rotation therapy and percussion and vibration modules	Optional	Optional
Preliminary Tilt Table: Extends Reverse Trendelenburg to achieve 20° to gradually acclimate your patient to the upright position	+	+
FullChair® Position: One-button FullChair® position means patients are up in upright FullChair® position with dependent leg positioning for "up in chair" orders	+	+
Chair Egress with Stand Assist: Use the bed to help stand heavier or less sturdy patients	+	+
Integrated Scale	Optional	Optional
Battery Backup: Backup for all bed articulation (high/low position, head-of-bed, knee section)	+	+
	<i>(Dual Electric and Manual Foot Pedal)</i>	<i>(Dual Electric and Manual Foot Pedal)</i>
Point-of-Care® Siderail Controls: All caregiver controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs	+	+
Head-of-Bed Compliance: Multiple digital head-of-bed angle indicators and a head-of-bed alarm make it easier to maintain head-of-bed compliance	+	+
	<i>(Able to store up to 8 days of statistics)</i>	<i>(Able to store up to 8 days of statistics)</i>
Drop Down Siderails: One handed single action siderail operation	+	+
Boost® One-Button Function: Puts the bed in the ideal position for patient repositioning	+	+
FlexAfoot™ Mechanism: Allows caregiver to customize via siderail controls the overall length of the bed and mattress, in any bed position (including chair position) to fit the patient	+	+
Ability to send data to Nurse Call, EMR or other Clinical Applications	Optional	Optional
Fully Integrated Surface: Surface automatically adjusts to the patient's weight, size, movement and to the head-of-bed articulation, to improve patient comfort	Optional	Optional <i>(40" Standard Width)</i>
Advanced Microclimate® technology: Helps reduce heat and moisture at the patient sleep surface interface	Optional	Optional
Status Board Nurse Call Display via smart beds	Optional	Optional

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER 173672 rev 2 04-MAR-2014 ENG – US

Zinc Number US-AC-1213-0662

www.hill-rom.com

USA 800-445-3730
CANADA 800-267-2337

Enhancing outcomes for
patients and their caregivers:

Hill-Rom