

Nursing Products



DYND70104

Enema Bucket Sets

The complete kit contains an enema bucket, 54-inch (137 cm) tube with plastic shut-off clamp, tube cap, castile soap packet and moisture-proof underpad.

- » 1,500 cc bucket is graduated in both ounces and milliliters
- » 54-inch (137 cm) soft vinyl tube has a pre-lubricated, smooth, non-traumatic tip with rounded eyes for safety and easy insertion
- » Clamp on tube is designed for one-handed operation
- » Individually wrapped in a polybag
- » Non-sterile

Item No.	Description	Pkg
DYND70104	Enema Bucket Set with Clamp	50/cs



DYND70125

Perineal Cleansing Bottle

Regulate water flow by lifting cap and turning.

- » 8-oz. (236.5 mL) perineal cleansing bottle with screw top
- » Graduated
- » Not made with natural rubber latex
- » Non-sterile

Item No.	Description	Pkg
DYND70125	Perineal Cleansing Bottle	50/cs