

Dimensional & Technical Data

Product Dimensions			
Component	Size	Dimension mm (Inches)	
Sleepform Mattress	One size	950 X 950mm (37.25" X 37.25")	
Airflow Mattress	1 & 2	960 X 1180 X 15mm(37.75" X 46.5" X 0.6")	
	3 & 4	1000 X 2000 X 30mm (39.4" X 78.75" X 1.2")	
Airflow Knee Pillow	1 & 2	230 X 330 X 60mm (9" X 13" X 2.4")	
	3 & 4	450 X 650 X 60mm (17.7" X 25.6" X 2.4")	
Knee Roll	1 & 2 Small	140 X 50 X 500mm (5.5" X 2" X 19.7")	
	1 & 2 Large	140 X 100 X 500mm (5.5" X 2" X 19.7")	
	3 & 4 Small	195 X 75 X 900mm (7.7" X 3" X 35.4")	
	3 & 4 Large	195 X 125 X 900mm (7.7" X 5" X 35.4")	
Cushioned Chest & Hip Guides		min – max (mm)	min – max (inches)
1	Width	115 - 160	4.6 – 6.6
	Circumference	345 - 485	13.8 – 19.4
2	Width	160 - 196	6.4 – 7.8
	Circumference	478 - 601	19.1 - 24
3	Width	160 - 196	6.4 – 7.8
	Circumference	478 - 601	19.1 - 24
4	Width	175 - 328	7 – 13.1
	Circumference	550 - 1007	22 – 40.2
Cushioned Leg & Knee guides		min – max (mm)	min – max (inches)
1	Circumference	137 - 271	5.5 – 10.8
2	Circumference	181 - 322	7.2 – 12.8
3	Circumference	181 - 322	7.2 – 12.8
4	Circumference	230 - 553	9.2 – 22.1
Under mattress Straps	1 & 2	630mm (24.8") Strap length	
	3 & 4	1020mm (40.2") Strap length	
Temperature Control Sheet	1 & 2	810 X 1220mm (32.4" X 48.8")	
	3 & 4	1000 X 2000mm (40" X 80")	
Waterproof Mattress Protector	One size	860 X 760mm (34" X 30")	

Technical Considerations

Medical Devices Regulations 2002	
Component	Complies
Sleepform Mattress	BS-EN 1021-1 and 2; 1994 BS EN 597-1:1995 BS EN 597-2:1995 BS 7177:1996 Complies with NHS Firecode Document HTM 87
Knee Rolls	BS EN ISO 12952-2 :1999 clause 3.3.4 (Pillows/Bolsters)
Air Flow Mattress Air Flow Knee Pillow	BS EN 597-1 :1995
Inner Cover	BS7175:1989:Sect.3, Ignition sources 0,1,5 100% Vapour Permeable Does not contain Latex
Outer Cover (100% Brushed Cotton)	100% Vapour Permeable Does not contain Latex
Foam Knee Roll	All foam used is CMHR (combustion modified high resilience) and complies with BS 5852, part 2, ignition source 5 and CAA Spec 8 Issue 2 Part 2.2b Foam is class V (BS3379).
Fibre Air Flow Mattress, Air Flow Knee Pillow, Leg & Knee Guides Chest & Hip Guides	Inherently low flame response – Schedule 2 Part 1 The Furniture and Furnishings (Fire) (Safety) Regulations 1988 (based on BS 5852:1982 Part 2) Performance Requirements for Air Flow - BS 6595 : 1985 : Clause 3.3 BS 4578 : 1970 (91)
Maximum User Weight Size 1&2 Air Flow Mattress Size 3&4 Air Flow Mattress	19kg (42 lbs) 76kg (167.2 lbs)
Cleaning Cycles	There are no restrictions on the number of cleaning cycles for this product
Product Warranty	2 years
Life Expectancy	Should exceed 5 years
Labelling	> CE marked for Class 1 Medical Devices (MDD 93/42 EEC) > Year of manufacture > Name of manufacturer > Material composition > Care Instructions

LECKEY[°]
sleepform

Leckey
Kilwee Business Park
Dunmurry
BT17 0HD
Northern Ireland
United Kingdom

T: 028 9060 2277
F: 028 9060 0795
E: info@leckey.com
W: www.leckey.com