

# InnoSpire Elegance



## InnoSpire Elegance and SideStream technology

Elegance and SideStream disposable and reusable nebulizers combine to provide fast and efficient drug delivery. SideStream nebulizers feature an active venturi system that provides air flow in addition to that supplied by the compressor, resulting in faster drug delivery and shorter treatment times. They are easy to use with fewer parts to assemble or clean than other conventional nebulizers.



**Venturi system** provides a faster delivery rate resulting in shorter treatment times than other conventional nebulizers

**Innovative design** consistently produces high-quality aerosol

The breath-enhanced SideStream Plus nebulizer features an easy-action inspiratory valve that opens on inspiration to boost medication delivery and closes on exhalation, preserving medication and reducing waste.<sup>1</sup> Less pressure is required to open the inspiratory valve compared to other breath-enhanced nebulizers.\* And there's no need to remove the valve when using a mask, making treatments easy for pediatric and elderly patients.



SideStream Plus

<sup>1</sup> Performance characterization of the new SideStream Plus breath enhanced jet nebulizer  
\* Data on file

Product specifications	
Weight	3.3 lbs.
Size	4.2" x 6.5" x 6.5"
Electrical rating/power	120V/60Hz, 1.7 Amp
Nebulization capacity	8 ml
MMAD	(80%) <5 microns
Treatment time	5-7 minutes
Maximum pressure	44 psi
Maximum flow rate	9.4 LPM
Operating flow rate	7 LPM @ 10 psi
Mode of operation	Continuous
Warranty	5 years

Ordering information		
Product description	Reorder number	NDC number
InnoSpire Elegance with SideStream disposable nebulizer	1100134	08373-0134-00
InnoSpire Elegance with SideStream disposable and reusable nebulizers	1099969	08373-9969-00
InnoSpire Elegance with SideStream Plus	1102828	08373-2828-00
Replacement filter (4-pack)	1102088	08373-2088-00
Standard carrying case	1101452	08373-0452-00
Premium carrying case	1026568	08373-0568-00
InnoSpire Elegance compressor with SideStream reusable nebulizer DVD	1101988	N/A
InnoSpire Elegance compressor with SideStream Plus nebulizer DVD	1101989	N/A

## Economical

Product description	HCPCS
InnoSpire Elegance compressor nebulizer system	E0570
SideStream disposable nebulizer	A7003
SideStream reusable nebulizer	A7005
SideStream Plus reusable nebulizer	A7005

Please visit [www.philips.us/elegance](http://www.philips.us/elegance) for more information



Respironics, InnoSpire, Elegance, and SideStream are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts this device to sale by or on the order of a physician.

DD MCI 4104969 PN 1101728

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International  
Headquarters  
+33 1 47 28 30 82

Philips Respironics Asia Pacific  
+65 6882 5282

Respironics New Jersey Inc.  
5 Wood Hollow Road  
Parsippany, NJ 07054  
+1 800 345 6443  
+1 724 387 4000  
[www.philips.com/respironics](http://www.philips.com/respironics)

[Click here to view detailed product information and pricing.](#) **Call 800-397-5899 to order today!**