



PRM256000

Drain Sponges

- » Rayon and polyester blend
- » Precisely fenestrated for easy and consistent application
- » Unique fenestration allows sponge to lie flat when wrapped around a drainage tube
- » Pre-cut sponge eliminates the need for sterile scissors and conserves nursing time
- » 6-ply, highly absorbent, low-linting, fast wicking
- » Not made with natural rubber latex

Item No.	Description	Pkg.
PRM256001	2 x 2" (5.1 x 5.1 cm)	35/bx, 20 bx/cs, 700 pk/cs
PRM256000	4 x 4" (10.2 x 10.2 cm)	25/bx, 12 bx/cs, 300 pk/cs
PRM256000H	4 x 4" (10.2 x 10.2 cm)	25/bx



PRM25710

Non-Adherent Pads

- » Sterile, cotton and polyester blend
- » Protects granulating tissue by not sticking to the lightly draining wound
- » Soft, perforated mylar film is porous, allowing air to circulate
- » Absorbent pad allows fluid to collect
- » Not made with natural rubber latex

Item No.	Description	Pkg.
PRM25700	2 x 3" (5.1 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25710	3 x 4" (7.6 x 10.2 cm)	100/bx, 12 bx/cs, 1200/cs



PRM21450

Abdominal (ABD) Pads

- » Extra absorbency and wicking for heavily draining wounds
- » Soft non-woven outer layer quickly wicks fluid to a cellulose center
- » Cellulose quickly absorbs and disperses fluids laterally to prevent pooling
- » Hydrophobic back impedes strike-through
- » Sealed edges prevent linting
- » Not made with natural rubber latex

Item No.	Description	Pkg.
----------	-------------	------

Sterile 1's

PRM21450	5 x 9" (12.7 x 22.9 cm)	25/bx, 16 bx/cs, 400/cs
PRM21453	8 x 7.5" (20.3 x 19.1 cm)	20/bx, 12 bx/cs, 240/cs
PRM21454	8 x 10" (20.3 x 25.4 cm)	16/bx, 20 bx/cs, 320/cs

Non-Sterile Bulk

PRM21451	5 x 9" (12.7 x 22.9 cm)	576/cs
PRM21456	8 x 10" (20.3 x 25.4 cm)	432/cs



NON255000

100% Cotton Drain Sponges

- » Natural breathability
- » Tightly woven into fenestrated shape for your convenience
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON255000	4 x 4" (10.2 x 10.2 cm), 16-ply, Sterile	1/pk, 25/bx, 12 bx/cs