

POC COMPARISON CHART

MANUFACTURER

3B Medical



GCE Zen O



Inogen



Chart



Respironics



POC

Aer X

Life

G3

Freestyle 5

Simply Go Mini

Height

8.38"

9.25"

8.25"

10.7"

9.4"

Depth

2.64"

3.8"

3"

4.4"

3.6"

Length

7.22"

9.8"

8.75"

6.6"

8.3"

Weight (lbs)
(including 1 Battery)

4.25

5.5

4.8

6.7

5

Hours
of Use
(1 battery)

Up to 4 hours
on pulse
setting of 2pm

Up to 4 hours on
pulse setting of
2pm

Up to 4.7 hours
on setting of
2pm

Up to 2.5 hours
on a setting of
2pm

Up to 4.5 hours
on setting of
2pm

DbA

40dbA
at 21pm

37dbA
at 21pm

39DbA
at 21pm

48DbA
at 21pm

43DbA
at 21pm

Patient
Changeable Sieve

Yes

Yes

Yes

No

No

Pulse Settings

1-5

1-5

1-5

1-5

1-5

Bolus

1100ml

1050ml

1050ml

1000ml

1000ml