

Adhesive Tapes

Re-engineered for just the right adhesion

CURAD® Paper Tape **Compare to Micropore®***

A gentle and breathable tape for sensitive skin, recommended especially for pediatric and geriatric patients. Flexible and conformable. Adheres well to damp skin. Single use rolls are great for reducing waste and cross contamination. **Latex-Free.**

Item#	Description	Quantity
CURAD® Paper Tape, Rolls		
NON270012	½" x 10 yds (1.3 cm x 9.14 m)	24 rolls/bx, 10 bx/cs, 240 rolls/cs
NON270001	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 10 bx/cs, 120 rolls/cs
NON270002	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 10 bx/cs, 60 rolls/cs
NON270003	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 10 bx/cs, 40 rolls/cs

Item#	Description	Quantity
CURAD® Paper Tape, Single Use Rolls		
NON270001S	1" x 1.5 yds (2.54 cm x 1.37 m)	100 rolls/bx, 5 bx/cs, 500 rolls/cs
NON270002S	2" x 1.5 yds (5.08 cm x 1.37 m)	50 rolls/bx, 5 bx/cs, 250 rolls/cs

CURAD® Cloth "Silk" Tape **Compare to Durapore®***

With its strong woven construction, this cloth "silk-like" tape is perfect for affixing dressings. It's water repellent and easy to tear. **Latex-Free.**

Item#	Description	Quantity
CURAD® Cloth "Silk" Tape, Rolls		
NON270112	½" x 10 yds (1.3 cm x 9.14 m)	24 rolls/bx, 10 bx/cs, 240 rolls/cs
NON270101	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 10 bx/cs, 120 rolls/cs
NON270102	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 10 bx/cs, 60 rolls/cs
NON270103	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 10 bx/cs, 40 rolls/cs


Item#	Description	Quantity
CURAD® Cloth "Silk" Tape, Single Use Rolls		
NON260101S	1" x 1.5 yds (2.54 cm x 1.37 m)	100 rolls/bx, 5 bx/cs, 500 rolls/cs



NON270001
1" x 10 yds Paper Tape

NON270002S 2" x 1.5 yds
Paper Tape, Single Use

Unique adhesive coating is gentle on sensitive skin.



NON260001E
1" x 10 yds Sensitive Skin Paper Tape

CURAD® Sensitive Skin Paper Tape **Revolutionary Gentle Adhesion**

Adhesive is present only on every alternate wave. This unique wave pattern provides plenty of surface area for the tape to stick, yet the smaller amount of adhesive is gentler to the skin when the tape comes off. The gaps in the adhesive also allow for breathability. **Latex-Free.**

Item#	Description	Quantity
CURAD® Sensitive Skin Paper Tape, Rolls		
NON260001E	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
NON260002E	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
NON260003E	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs



NON270101 1" x 10 yds
Cloth "Silk" Tape

NON260101S 1" x 1.5 yds
Cloth "Silk" Tape, Single Use

* This product is not manufactured or distributed by 3M Healthcare, makers of Micropore® and Durapore®.