



ESSENTIAL[®]
MEDICAL SUPPLY, INC.

PRODUCT CATALOG

2019

In this Catalog

We set ***THE BAR*** for Product Quality...
...***Because Quality Pays - It Doesn't Cost!***

AMBULATORY PAGES 3-19



MOBILITY ACCESSORIES PAGE 15-17



Everyday **ESSENTIALS** ADL'S PAGES 20-26



ESSENTIAL BATH SAFETY PAGES 27-33



BEDDING PAGES 34-37



CUSHIONS PAGES 38-46



SHEEPETTE® AND DECUBITUS PAGES 47-48



ENDURANCE HD **HEAVY DUTY** PAGES 49-51



Quik-Sorb™ INCONTINENT PRODUCTS PAGES 52-54



PT/OT PAGES 55-57



PATIENT CARE PAGES 58-63



PLANOGRAMS PAGE 64-65



ESSENTIAL things to know about **Essential Medical Supply, Inc.®**



ORLANDO HEADQUARTERS

2019 is our 35th year in business and we have **YOU**, our customers, to thank for it!


We ship from **TWO** Essential Medical Supply® Warehouse locations:
Orlando, FL | Las Vegas, NV



LAS VEGAS WAREHOUSE

We have a top notch Customer Service Department that truly cares.

We are Family Owned and Operated.





Contact us for information on our Marketing Program to increase your sales and rentals.

FREE Spirit[®]

KNEE and LEG WALKER

Superbly Engineered and Patented

It Is Truly State of the Art!

With its unique features and outstanding quality, the Free Spirit[®] will maximize your rental returns and minimize your problems.

- Large, soft and comfortable knee pad
- Easy to use tiller with one hand locking mechanism
- Folding strap to secure the unit during transport
- Includes carrying pouch with zippered pocket
- Extra wide tires and rack and pinion steering design prevents overturning and increases stability

P4000 Free Spirit[®] Knee and Leg Walker

Adjusts for users 5'3" to 6'8"

400 lb. User Weight Limit!


Patent #7,780,180,132

SPECIFICATIONS

Overall Length:	31¼"
Width:	17¼"
Range of Height Adjustments:	
Knee Rest:	17½" - 23½"
Handle Bars:	34" - 42"
Weight Limit:	400 lbs.
Weight:	22.5 lbs.
Case Pack:	1



Folds for Transportation



Unique Steering



Height Adjustable

SHEEPETTE[®] KNEE WALKER PADS

D5011 Sheepette[®] Knee Walker Pad
11" (original) Free Spirit[®] Knee Walker Pad

D5012 Sheepette[®] Knee Walker Pad
14" (current) Free Spirit[®] Knee Walker Pad

- Thick, plush Sheepette[®] material provides superior user comfort
- Specifically designed for the Free Spirit[®] Knee Walker.
- Elastic material ensures snug fit while protecting the vinyl knee pad underneath
- Sold one (1) per retail package



Increase Profit On Every Sale and Rental!

FASHION CANES



Couture™ Offset Canes with Matching Super Bigfoot™ Cane Tip

- First and ONLY cane with a matching standing cane tip
- Adjustable, 7/8" aluminum shaft
- Thick, comfortable foam handle with wrist strap is great for arthritic hands



W1343C
W1343F
W1343H

Celebration
Sunflower
Houndstooth

W1343P
W1343S
W1343MC

Paradise
Swirl
Mixed Carton

W1343 SPECIFICATIONS
Height Adjustment: 31" - 40"
Weight Limit: 250 lbs.
Sold in cartons of 2

*One of each design
- 5 pc per carton*




You Can Feel The Difference!

- Ergonomic hand grip with finger notches
- Comfortable handle is latex free
- Great choice of vibrant designs or a matte color finish

Use **T10034** replacement tips

- W1346DB Danube Blue**
- W1346G Graphite**
- W1346PA Pop Art**
- W1346S Stone**
- W1346T Tropical**
- W1346ABK Black**
- W1346AB Blue**
- W1346APK Pink**
- W1346APW Pewter**



W1346 and W1347 SPECIFICATIONS
W1346 Height Adjustment: 30 1/4" - 39 1/4"
W1347 Height Adjustment: 34" - 38"
Weight Limit: 250 lbs.
Sold in Eaches



- All the features of Gentle Touch™ Offset Canes, but with a folding shaft
 - Folds down to 13"
 - Includes plastic carrying case
- Use **T10034** replacement tips

W1347DB Danube Blue



Laser Cuit™
fashion canes



Detail of Shaft

Truly Dazzling Canes... at a Great Price!

- Adjustable, 7/8" aluminum shaft
- Thick, comfortable foam handle with wrist strap is great for arthritic hands

Use **T10034** Replacement Tips

- | | | | |
|----------------|--------------|---------------|---------------|
| W1580B | Blue | W1580P | Purple |
| W1580BK | Black | W1580R | Red |
| W1580G | Green | | |

W1580 SPECIFICATIONS

Height Adjustment: 30" - 39"
Weight Limit: 250 lbs.
Sold in Eaches

Designer Offset Canes With Ribbed Foam Handle

- Adjustable, 7/8" aluminum shaft
- Ribbed foam grip for comfort
- Unique tube within a tube handle for added strength
- **300 lb. weight capacity**

Use **T10034** Replacement Tips

- | | | |
|---------------|-----------------|---|
| W1344A | Amethyst |  |
| W1344B | Bouquet |  |
| W1344H | Halo |  |
| W1344O | Ocean |  |
| W1344P | Paradise |  |
| W1344R | Retro |  |



W1344 SPECIFICATIONS

Height Adjustment: 29" - 38"
Weight Limit: 300 lbs.
Sold in Eaches


Designer Offset Canes

Great Designs... Guaranteed to FLY Out of Your Store!

- Adjustable, 7/8" aluminum shaft
- Unique tube within a tube handle for added strength
- Thick, comfortable foam handle with wrist strap is great for arthritic hands

Use **T10034** Replacement Tips

- | | | | |
|----------------|--------------------|----------------|----------------------|
| W1345B | Butterfly | W1345M | Morning Glory |
| W1345CL | Celestial | W1345P* | Paisley |
| W1345F | Floral | W1345SP | Scotch Plaid |
| W1345H | Houndstooth | W1345Y | Burgundy |




W1345 SPECIFICATIONS

Height Adjustment: 29 1/4" - 38 1/4"
Weight Limit: 300 lbs.
(excluding W1345P)
*W1345P 250 lbs.
Sold in Eaches

SPRING GARDEN COLLECTION®

LADY'S ADJUSTABLE FASHION CANES




Adjustable Straight Canes

- Lighter 5/8" shaft
 - Genuine wood derby handle sized for smaller hands
- Use **T10058** Replacement Tips

- W1550L Lilac**
- W1550R Rainbow**
- W1550T Tiger**
- W1550W Wildflower**

W1550 SPECIFICATIONS

Height Adjustment: 30¼" - 39¼"
 Weight Limit: 250 lbs.
 Sold in Eaches



Adjustable Folding Canes

- Lighter 5/8" shaft
 - Genuine wood derby handle sized for smaller hands
 - Includes plastic carrying case and folding lock
- Use **T10058** Replacement Tips

- W1400B Butterfly**
- W1400L Lilac**
- W1400W Wildflower**

W1400 SPECIFICATIONS

Height Adjustment: 33" - 37"
 Weight Limit: 250 lbs.
 Sold in Eaches



Demi Folding Canes

- The first cane designed for users 4'8" to 5'4" tall
 - Genuine wood derby handle sized for smaller hands
 - Includes plastic carrying case and folding lock
 - Compact design easily fits into a purse
- Use **T10058** Replacement Tips

- W1351L Lilac**
- W1351W Wildflower**



W1351 SPECIFICATIONS

Height Adjustment: 27½" - 32½"
 Weight Limit: 250 lbs.
 Sold in Eaches

STEPPIN' OUT[®]

FASHION CANES


Adjustable Straight Canes

- 2 fashion finishes to choose from:
Textured or Carbon Finish
- Comfortable wooden T-handle
- Lightweight 7/8" aluminum shafts
- Great selection of "fashion colors"

Use **T10034** Replacement Tips

- W1430BL** Textured - Black
- W1430DB** Textured - Blue
- W1430R** Textured - Red
- W1442L** Carbon Finish - Charcoal Gray

*Great Patterns That
Your Customers Will Love!*



W1430 and W1442 SPECIFICATIONS

Height Adjustment: 30" - 39"
Weight Limit: 250 lbs.
Sold in Eaches

Adjustable Folding Canes

- Includes plastic carrying case
- Comfortable wooden derby handle
- Lightweight 7/8" aluminum shaft

Use **T10034** Replacement Tips

- W1350BL** Textured Black 
- W1350DB** Textured Blue 



W1350 SPECIFICATIONS

Height Adjustment: 33" - 37"
Weight Limit: 250 lbs.
Sold in Eaches



W1540 Adjustable Offset Cane

- Adjustable, 7/8" aluminum shaft
- Unique tube within a tube handle for added strength
- Thick, comfortable foam handle with wrist strap is great for arthritic hands

Use **T10034** Replacement Tips



W1540 SPECIFICATIONS

Height Adjustment: 30" - 39"
Weight Limit: 250 lbs.
Sold in Eaches

QUAD CANES

For Endurance® Heavy Duty Quad Canes, see page 50

Designer Quad Canes



- Comfortable foam handle is great for arthritic hands
 - Use with either hand by rotating shaft
- Use **T50012** Replacement Tips

Small Base - 6" x 8"

- W1303B** **Burl**
 - W1303D** **Desert**
 - W1303F** **Flower**
 - W1303I** **Indigo**
 - W1303M** **Mixed Case**
- Contains: One (1) of each design



BURL **DESERT** **FLOWER** **INDIGO**

W1303 SPECIFICATIONS

Height Adjustment:	28¾" - 37¾"
Weight Limit:	250 lbs.
Weight:	1.7 lbs.
Sold in Case Pack of:	Same Design - 2 Mixed Carton - 4

Endurance® Quad Canes

- Comfortable foam handle is great for arthritic hands
 - Anodized aluminum shaft with welded steel base
 - Use with either hand by rotating shaft
- Use **T50012** Replacement Tips for Small Base
Use **T50058** Replacement Tips for Large Base

Small Base - 6" x 8"

- W1301S** **Silver**
- W1301BL** **Black**
- W1301BR** **Bronze**
- W1301M** **Mixed**

Contains: 2 Silver, 1 Black, 1 Bronze

Large Base - 8" x 12"

- W1302S** **Silver**
- W1302BL** **Black**
- W1302BR** **Bronze**

SPECIFICATIONS

Height Adjustment:	28¾" - 37¾" in 1" increments
Weight: Small Base	1.7 lbs.
Large Base	2.1 lbs.
Weight Limit:	300 lbs.
Sold in Case Pack of:	4

GENUINE WOOD CANES

GENUINE WOOD CANES



Black

Mahogany

Natural

Endurance® Wood Canes

- 1" diameter shaft for strength and durability
- Decorative brass ring
- Large bend for hand comfort

CURVED HANDLE CANES

- W1539BK** **Black**
- W1539M** **Mahogany**
- W1539N** **Natural**

W1539 SPECIFICATIONS

Length:	36"
Weight Limit:	250 lbs.
Sold in Eaches	

STANDARD CANES

For Endurance® Heavy Duty Offset Canes, see page 50

Offset Handle Canes

- Height adjustable in 1" increments
- Lightweight 7/8" anodized aluminum shaft
- Comfortable foam handle is great for arthritic hands

Use **T10034** Replacement Tips

- W1340BL Black**
- W1340BR Bronze**
- W1340S Silver**
- W1340M Mixed Case**

Contains: 2-Silver, 2-Bronze, 2-Black



W1340 SPECIFICATIONS

Height Adjustment: 29" - 38"
Weight Limit: 300 lbs.
Sold in Case Pack of: 6

Curved Handle Canes

- Height adjustable in 1" increments
- Lightweight 7/8" anodized aluminum shaft
- Comfortable foam handle is great for arthritic hands

Use **T10034** Replacement Tips

- W1510BL Black**
- W1510BR Bronze**
- W1510S Silver**



W1510 SPECIFICATIONS

Height Adjustment: 31" - 40"
Weight Limit: 250 lbs.
Sold in Case Pack of: 6

Endurance® Vinyl Handle Canes

- Choice of 3/4" or 7/8" diameter shaft
- Contour vinyl handle for firm grip
- Locking ring for added safety
- Height adjustable in 1" increments

W1505 Lighter 3/4" anodized aluminum shaft (for women)

Use **T10058** Replacement Tips for W1505

W1506 Anodized 7/8" aluminum shaft (for men and women)

Use **T10034** Replacement Tips for W1506



W1505 & W1506 SPECIFICATIONS

Height Adjustment: 31" - 40"
Weight Limit: 250 lbs.
Sold in Case Pack of: 6

W1420BL Black Deluxe Folding Cane

- Compact cane folds into 15" plastic carrying pouch
- Lightweight 7/8" anodized aluminum shaft
- Unique plastic handle has ergonomic design

Use **T10034** Replacement Tips



W1420 SPECIFICATIONS

Height Adjustment: 30" - 39"
Weight Limit: 250 lbs.
Sold in Eaches

FOLDING SEAT CANES

THREE LEGS

W1450 Folding Seat Cane

- Lightweight, anodized aluminum frame
- High impact plastic seat and handle

Weights 2 lbs.



W1450 SPECIFICATIONS

Seat Size: approx. 7½" x 7½"
 Seat Height: 21½"
 Cane Length: 34¼"
 Weight Limit: 225 lbs.
 Sold in Eaches

FOUR LEGS

W1451 Folding Seat Cane

- Four legs for added stability
- Large, comfortable canvas seat
- Lightweight, anodized aluminum frame

Weights 2½ lbs.



W1451 SPECIFICATIONS

Seat Size: 14" L x 6" D"
 Seat Height: 21"
 Cane Length: 35"
 Weight Limit: 225 lbs.
 Sold in Eaches



ESSENTIAL BLIND CANE



W1311 Folding Aluminum Blind Cane

- 52" long
- Red stripe on bottom
- Wrist strap on handle
- Sold one (1) per retail package

CANE ACCESSORIES



W1340-FHS



W1340-FHBR



W1340-FHBL



W1340-SBL

Offset Cane Replacement Handles and Wrist Straps

- Will fit most brands of offset canes
- Latex free
- Sold one (1) per retail box

Foam Handles

W1340-FHBL

Black

W1340-FHBR

Bronze

W1340-FHS

Silver

W1344-FHBL

Ribbed



W1344-FHBL

Wrist Strap

W1340-SBL

Black

CANE ACCESSORIES

W1319 Cane Holder

- Easily attaches to any cane
- Allows user to suspend cane from flat surfaces
- Sold one (1) per retail package



W1340-USBL Black Universal Cane Strap

- Fits all canes; simply slip over shaft
- Sold one (1) per retail package



T20009 BigFoot™ Standing Cane Tip

- Enables canes to stand on their own
- Fits 3/4" shaft canes
- 250 lbs. weight limit
- Sold one (1) per retail package



T20005 Ice Attachment

- Easily attaches to most canes
- Five (5) stainless steel prongs for extra stability
- Sold in eaches
- Sold one (1) per retail package



Super Bigfoot™ Standing Cane Tip

- Six feet provide superb stability
- Fits 3/4" shaft canes
- Enables canes to stand on their own
- 250 lbs. weight limit
- Sold one (1) per retail package



T20010B Black
T20010P Pink

REPLACEMENT TIPS

- Choice of tips in retail boxes for easy display
- All tips include metal washers for long wear



Cane Replacement Tips

Packed two (2) per retail box

- T10034 Fits 3/4" Shaft
- T10058 Fits 5/8" Shaft
- T10078 Fits 7/8" Shaft



Walker and Commode Tips

Packed four (4) per retail box

- T40100 Fits 1" Shaft
- T40118 Fits 1 1/8" Shaft



SPECIFICATIONS

Retail Box

Cane Tips:	2 Packed
Quad Tips:	4 Packed
Walker, Commode, Shower Bench Tips:	4 Packed



Quad Cane Tips

Packed four (4) per retail box

- T50012 Fits 1/2" Shaft
- T50058 Fits 5/8" Shaft



Shower Bench Tips

Packed four (4) per retail box

- T60118G Fits 1 1/8" Shaft - Gray Only

Tips Are Available In Gray Or Black

For Replacement Tips
Add: G - for Gray
BL - for Black to code Number
(e.g. T10034G)

CANE DISPLAYS



Couture[™]
Offset Canes
Display

SW1343-Display Couture™ Offset Canes Display

Includes:

- Dynamic **FREE** display that shows your Couture™ Canes at their very best!
- Graphics that highlight the features and benefits of these unique canes
- The **FIRST** and **ONLY** canes that feature a matching pattern on Super BigFoot Tips™
- A mix of 12 Couture™ Canes



SW726

Fashion Cane Displays

Maximize your cane sales by letting your customers see what you have!

- Open display shows your canes at their best!
- Compact, takes up a minimum of floor space
- Includes **FREE** deluxe cane stand

SW711 18 pc. Laser Cut™, Designer Offset and Spring Garden Collection® Canes with **FREE** W1318 Cane Stand

SW718 12 pc. Designer Offset Canes with **FREE** W1316 Metal Cane Stand

W1318 Deluxe Wood Cane Stand, holds 18 canes, easy to assemble, includes **FREE** sign

SW726 Rolling Mobility Center

- Four (4) standard offset canes
- Fourteen (14) offset fashion canes
- Twelve (12) folding canes
- Six (6) styles of replacement tips
- **FREE** rolling display rack and signage



SW711

Seasonal Cane Displays

Maximize your cane sales by using our colorful displays designed to **highlight the season**. All displays contain twelve (12) canes and **FREE** display stand.



SW730-Father



SW730-Holiday



SW730-Summer

SW730-Father

- Twelve masculine canes that make perfect Father's Day Gifts

SW730-Holiday

- Colorful canes designed for the holiday season

SW730-Summer

- Canes that make great Mother's Day gifts and a super way to greet spring and summer

CRUTCHES AND CRUTCH ACCESSORIES

ENDURANCE® Crutches

- Reinforced center tube for superior weight bearing strength
- Silencing collar improves stress tolerance at all angles
- Easy to read calibrated height scale for adjustment to a wide range of heights
- Easy one button height adjustment
- 250 lb. weight limit
- Sold one (1) pair per retail package



		User Height
W4000	Child	4' 0" - 4' 8"
W4001	Youth	4' 6" - 5' 2"
W4002	Medium Adult	5' 2" - 5' 10"
W4003	Tall Adult	5' 10" - 6' 6"

All Endurance® Crutches come completely accessorized with attractive gray pads - no assembly required

D5009 Sheepette® Crutch Cover Set

- Thick, plush Sheepette® material provides superior user comfort
- Set includes pair of crutch arm and grip covers
- Each set outfits one (1) pair of crutches
- Sold one (1) set per retail package

INCREASE YOUR PROFIT ON EVERY CRUTCH SALE OR RENTAL!



T70040 Crutch Accessory Kit

- Everything you need to rehab your old crutches
- Packed in retail box for easy display
- Sold one (1) kit per retail box

Each kit contains 2 each: T70034 Crutch Tips, T70010 Crutch Arm Pads, T70020 Crutch Hand Grips



Below sold two (2) per retail box, for easy display

T70010 Crutch Arm Pads **NEW PACKAGING**

- Fits a standard crutch
- Soft and cushioned to ensure patient comfort
- Reinforced for longer wear



T70020 Crutch Hand Grips **NEW PACKAGING**

- Fits a standard crutch
- Soft and cushioned to ensure patient comfort
- Reinforced for longer wear



T70034 Crutch Tips **NEW PACKAGING**

- Fits either ¾" or 1" shafts
- Reinforced with metal washers



NEW PACKAGING

T70050 Forearm Crutch Tips

- Reinforced with metal washers
- 7/8" diameter



For Endurance® Heavy Duty Walkers, see page 51



ENDURANCE®
— BUILT TO LAST! —

W1260 Trigger Release Walker
W1260-5F Trigger Release Walker with 5" Fixed Wheels

- Easy to open and fold . . . fingers, palm or hand can be used
- Available with and without 5" fixed wheels

W1260 SPECIFICATIONS	
Width:	23"
With Wheels:	24½"
Inside Width:	17"
Weight:	5 lbs 1 oz
Tubing Size:	1" aluminum
Height Adjustment:	32½" - 39½"
User Weight Limit:	300 lbs.
Case Pack:	4
HCPCS:	W1260 - E0135 W1260-5F - E0143

All Endurance® Trigger Release Walkers include a free pair of Glide Caps and Carpet Glides...
One for use outdoors; one for use on carpets



ENDURANCE®
— BUILT TO LAST! —

W1258-5F Trigger Release Junior Walker with 5" Fixed Wheels

- A true junior walker designed to meet the needs of users under 5'4" tall
- Available in silver with 5" fixed wheels

W1258 SPECIFICATIONS	
Width:	22½"
With Wheels:	24"
Inside Width:	16¾"
Weight:	4 lbs 14 oz
Tubing Size:	1" aluminum
Height Adjustment:	25¼" - 32¼"
User Weight Limit:	300 lbs.
Case Pack:	4
HCPCS:	W1258 - E0135 W1258-5F - E0143

W1262F Silver 5" Fixed Wheels for Endurance® Trigger Release Walkers

- Sold one (1) pair per package



SPECIFICATIONS	
Width:	17"
Depth:	13"
Height Adjustment:	28" - 33"
User Weight Limit:	250 lbs.
Case Pack:	2
HCPCS:	E0135

W1300 Pyramid Cane

- Sturdy, ¾" aluminum tubing for lightweight stability
- The advantages of a standard cane with the support of a walker
- Contoured foam hand grips for patient comfort

Use **T10078** Replacement Tips

Excellent Product For Stroke Victims

FOLDING WALKER ACCESSORIES

W601 Walker Skis

- Will fit most walkers with 1" tubing
- Easy to install
- Constructed of high density plastic
- Sold one (1) per retail box



W604 Universal Walker Skis

- Adjustable cylinder with tightening screw fits both 7/8" and 1" tubing
- Easy to install
- Sold one (1) per retail box



W605 Pop 'N Glide™ Walker Skis

- Will fit most walkers with 1" tubing
- Constructed of high density plastic with rubber tip
- Rear legs just "pop" into the skis
- Sold one (1) per retail package



W602 Glide Caps

- Rear leg glide caps for 1" walkers
- Allows smooth, effortless walking when used with front wheels
- Sold one (1) per retail package



W603 Carpet Glides

- Great for use indoors
- Slip over rear walker tips
- Sold one (1) per retail package



W1245 Universal Walker Wheels

- 5" fixed wheels will fit most brands of 1" walkers
- Multiple hole sizes fit various pins
- Sold one (1) per retail box



Universal wheels also include one (1) pair of W602 Glide Caps and one (1) pair of W603 Carpet Glides

W1617 Deluxe Walker Tray

- Fits MOST folding walkers
- Raised edges help prevent spillage
- Two beverage holders
- Great consumer packaging that takes up less shelf space
- Sold one (1) per retail package



W1216 Deluxe 3 Pocket Walker Pouch



- Will fit all walkers with two front cross braces
- Hook and loop closure on pockets
- Black color
- 11" L x 13" W capacity
- Sold one (1) per retail package

W1616 Easy Basket™

- Will fit all brands of walkers with 2 front cross braces
- Easy to slip on and lift off
- No tools or difficult straps required
- White vinyl coated metal for easy cleaning
- 14½" L x 9" D x 8" W capacity
- Sold one (1) per retail package



W1614 Universal Walker Basket

- Will fit most brands of walkers via wire arms
- White vinyl coated metal for easy cleaning
- 16" L x 5½" D x 7" W capacity
- Sold one (1) per retail package



W1619 Walker Basket

- Will fit all brands of walkers with 2 front cross braces
- Attaches via hook and loop straps
- White vinyl coated metal for easy cleaning
- 16" L x 7" D x 5½" W capacity
- Includes plastic tray insert and cup holder
- Sold one (1) per retail package



D5008 Sheepette® Walker Grip Covers

- Thick, plush Sheepette® material provides superior user comfort
- Great for users with arthritic hands
- Sold one (1) pair per retail package

D5010 Sheepette® Platform Walker Pad Cover



- Thick, plush Sheepette® material provides superior user comfort
- Will fit almost all platform walker brands
- Sold one (1) per retail package

W1219 Platform Walker Attachment



- Designed for users who cannot use one hand
- Will fit most brands of walkers with two side cross bars
- 250 lb. user weight capacity
- Sold in eaches

WHEELCHAIR ACCESSORIES

H1301 Wheelchair Bag

- Large, waterproof nylon bag with zippered closure
- Has loop straps to fit over wheelchair handles
- Reflector strip on the back for added user safety!
- 13" H x 13" W x 6" D capacity
- Sold one (1) per retail package



H1201 Deluxe Rosewood Wheelchair Tray

- Extra wide, adjustable straps
- Raised lip to keep items on the tray
- Lighter weight than most other wheelchair trays
- 24" W x 20" D surface area
- Sold one (1) per retail package



H1303 Universal Cup Holder

- Easy to attach to wheelchair or transport chair arms
- Can also be used with four wheel and folding walkers
- Sold one (1) per retail package




Deluxe Wheelchair Pouches

- Can be used on transport chairs, wheelchairs and folding walkers
- Soft, quilted material
- Attaches via loops or hook and loop flap
- Pouch closes via hook and loop
- Sold one (1) per retail package

W4550 Confetti

W4551 Pinpoint



Detail of Confetti



Detail of Pinpoint



WHEELCHAIR TRANSFER BOARDS

P2300 Hardwood Transfer Board

- Highly polished and sealed for easy transfer
- 8" W x 30" L; two hand cut out
- 250 lbs. weight capacity
- Sold one (1) per retail box



P2301 Plastic Transfer Board

- Smooth surface for ease of transfer
- 8" W x 27½" L
- 250 lbs. weight capacity
- Sold one (1) per retail box





Includes both a pouch and a basket



Unique folding strap secures the unit during transportation and allows it to stand on its own



Curved flip up back and loop brakes





Large adjustment knobs for arthritic hands



Extra thick soft padded seat

Featherlight™

FOUR WHEEL WALKER

- W1650B-1  Blue
- W1650FG-1  Forest Green
- W1650R-1  Red
- W1650S-1  Silver
- W1650T-1  Teal



demi Featherlight™

FOUR WHEEL WALKER

- Designed for users 4'8" to 5'4" tall
- All of the great features of our W1650-1 Standard Featherlight™ plus:
 - Lower seat that enables users under 5'4" to rest their feet on the ground
 - Lower handle adjustment for user comfort and safety

- W1660B-1  Blue
- W1660R-1  Red
- W1660T-1  Teal



W1650 Standard

W1660 Demi

The Blazer®

ALUMINUM FOUR WHEEL WALKER

- W1630B-12**  **Blue**
- W1630BK-12**  **Black**
- W1630R-12**  **Red**

Large adjustment knobs for arthritic hands

Curved, flip up back

Securing strap for storage and transportation

Large 8" x 1½" wheels, great for use in and outdoors



Loop hand brakes for safety and security

Large, padded seat

Includes both basket and carrying pouch

Essential Rollator

LIGHTWEIGHT STEEL FOUR WHEEL WALKER

- W1700B-14**  **Blue**
- W1700BK-14**  **Black**
- W1700R-14**  **Red**

Loop hand brakes

Curved, padded back for comfort

Padded seat for added comfort



Includes carrying pouch

6" spoke wheels with sealed bearings

	Featherlight™ Walker (W1650-1)	Demi Featherlight™ Walker (W1660-1)	Blazer® Walker (W1630-12)	Essential Rollator (W1700-14)	Endurance HD® Walker (on page 51) (W1802B-14)
Wheel Size	6"	6"	8"	6"	7½" x 1½"
Weight Capacity	250 lbs.	250 lbs.	250 lbs.	300 lbs.	500 lbs.
Height Adjustments	32¼" - 37½"	28¼" - 32½"	31" - 37"	30½" - 35"	31¾" - 36"
Overall Width	24¾"	24¾"	23½"	23¾"	26¾"
Width Between Handles	17½"	17½"	17½"	17¾"	22"
Weight	15½ lbs.	15½ lbs.	18 lbs.	16 lbs.	31 lbs.
Seat Size	14" W x 13" D	14" W x 13" D	14" W x 13" D	12" W x 12" D	18½" W x 14¼" D
Seat Height from the Floor	22"	18"	22"	20½"	22"
Case Pack	1	1	1	1	1

All Everyday Essentials[®] sold one (1) per retail package



Dressing / Bathing



L3008 Standard Sock Aid

- Flexible to help put on socks or stockings
- Adjust straps by cutting or knotting to desired length



L3020 Wooden Dressing Stick

- Facilitates positioning and putting on clothing
- Features plastic coated "C" hook at one end; specially shaped push/pull hook at the other end
- 27" Long



L3009 Deluxe Sock Aid

- Molded design makes it easier to put on socks or stockings
- 48" strap has handle for ease of application
- Includes garter clips



L3023 Plastic Shoehorn

- 18" Long
- Has hook at one end to position shoes
- Hanging loop for storage



L3010 Terry Cloth Covered Sock Aid with Foam Handles

- 24" straps with two foam grips for easy sock application
- Terry cloth helps secure socks during application
- Will fit any size sock



L3018 Shoehorn/Dressing Stick

- Extra thick handle is easy to hold
- Features specially shaped push/pull hook at one end and shoehorn on the other.



L3012 Deluxe molded Sock Aid with Foam Handles

- 31" straps with two foam grips for easy sock application
- Molded design makes it easier to put on socks or stockings



L3024 Deluxe Metal Shoehorn

- Chrome plated, 23" long
- Head bends for easier application
- Large, easy to grip handle with hanging loop



L3019 Adjustable Shoehorn

- Adjustable between 19" - 28"
- Easy to lengthen. Just twist to lengthen and then twist again to tighten
- Easy to grip handle
- Great for travel. Fits into the smallest suitcase



L3027 Jumbo Metal Shoehorn

- Solid steel construction
- 23½" Long
- Hanging loop for storage



Dressing / Bathing

Elastic Shoelaces

- Tie once, and sew thru knot
- Enables users to slip shoes on/off without untying laces
- Each box has three (3) pairs

L3025BK Black 24"
L3025BR Brown 24"
L3025WH White 24"

L3026BK Black 32"
L3026BR Brown 32"
L3026WH White 32"



L3022 Plastic Handle Button/Zipper Pull

- Stainless steel hook at one end for buttons and a coated "C" hook at the other end for zippers



L3046 Cleansing Body Puff

- Extra long handle
- Gently exfoliates the skin
- Great retail packaging



L3040 Bathing Sponge

- Round, firm sponge
- Long handle for ease of use
- Hanging loop on handle
- 20" long
- Sponge is 4³/₄" diameter



L3041 Deluxe Bathing Sponge

- Large, hourglass shape sponge
- Long, 20" angled handle to reach the back
- Hanging loop on handle



L3042 Lotion EZE™ Applicator

- Soft, washable pad applies lotion to hard to reach areas
- Comfortable, ergonomic handle
- 14" Long
- Latex free



L3043 Lotion EZE™ L3042 Replacement Pads

- Package of two (2) foam replacement pads
- Latex free



L3045 Deluxe Shampoo Basin

- Inflatable heavy duty vinyl with two air chambers
- Long, attached drainage hose
- Extra large size: 24" L x 28" W x 6" H



L3050 Deluxe Handy Nail Brush

Absolutely the BEST Nail Brush you'll ever use!

- Easy to grasp
- Handy suction cups allow you to attach it to the wall



L3060 Personal Hygiene Wipe

- Easy to use for personal hygiene
- Extends the users' reach more than 12"
- Works with toilet paper or pre-moistened wipes
- Features an easy release button to discard paper after use





Household



L2002 Tablet Splitter

- Cuts tablets in half for smaller dosage, easier swallowing
- Front tray stores extra tablets



L2003 Tablet Crusher

- Screw top crushes tablets into a fine powder
- Compartment in front stores spare tablets



L2006 Key Turner

- Provides large grip and improved leverage
- Holds three (3) metal keys
- Keys fold in when not in use



L2008 Multi Opener

- Makes it a **snap** to open:
 - Plastic Bottles
 - Can Tabs
 - Bottle Caps
- Can also be used to break the vacuum seal on glass jars



L3007 Leg Lifter

- Enables user to lift their leg on/off beds, footrests, etc.
- Reinforced loop fits around foot
- Total length is 41"



Reading / Writing



L4001 Lighted Magnifier

- Rimless design allows for MAXIMUM viewing area
- Provides 3X magnification in main viewing area
- Comes with protective cloth pouch and batteries



L4002 Folding Magnifier

- Folding handle provides adjustable angle
- 2X and 4X times magnification lenses
- Large size: 4" L x 2" W
- Soft plastic carrying case

Contact us for
information on our
Marketing Program to
increase your referrals
and sales



Created by
Essential Medical to make
YOUR EVERYDAY
easier and more enjoyable!





Eating / Drinking



L5001 Bendable Spoon

- Soft cushion grip is easy to hold
- Flexible section allows it to be bent to desired angle
- Dishwasher safe
- Latex free



L5002 Bendable Fork

- Soft cushion grip is easy to hold
- Flexible section allows it to be bent to desired angle
- Dishwasher safe
- Latex free



L5003 Rocker Knife

- Soft cushion grip is easy to hold
- Rounded knife blade utilizes a rocking motion for easier cutting
- Dishwasher safe
- Latex free



L5004 Utensil Strap

- Ideal for individuals with limited grip
- Attach strap to utensil and slip hand thru loop



L5005 Utensil Set

- Complete set that includes bendable spoon, bendable fork and rocker knife
- Dishwasher safe
- Latex free

At Essential, we can help you create planograms to maximize your store space and sales!



AS60 Everyday Essentials[®] Planogram

Power of RED[™]

See the
Difference!



Studies have shown that red plates and utensils increase eating for Alzheimer's patients.

L5031 Scoop Bowl

L5032 Scoop Dish

- Rim and reverse curve aid in scooping food onto utensils
- Includes suction base to prevent sliding
- Bowl and dish (not suction bases) are dishwasher safe
- Sold one (1) per retail box



L5033 Nose Cutout Cup

- Shaped to provide room for nose
- Ideal for people who cannot tilt their head back due to neck injuries or choking
- 8 oz. fluid capacity
- Dishwasher safe
- Sold one (1) per retail box



Individual Utensils

- All utensils feature large, easy to grip handles
- Fork and spoon are bendable to desired position
- Dishwasher safe
- Latex Free
- Sold one (1) per retail package

L5041 Bendable Spoon

L5042 Bendable Fork

L5043 Rocker Knife



L5041

L5042

L5043

L5045 Utensil Set

- Utensil set includes bendable fork, spoon and rocker knife
- Dishwasher safe
- Latex free
- Sold one (1) set per retail package



L5046 Complete Table Setting

Complete Table Setting

- Set includes bendable fork, spoon, rocker knife, cup, bowl and dish
- Sold one (1) set per retail box



Please see Page 65 for the Power of Red[™] Planogram

Please contact us for our Power of Red[™] Marketing Program that includes:

Planogram | Signage | Brochures | Marketing Tips

HIP AND KNEE KITS

- All Kits sold one (1) per retail package
- Everything needed to make life easier after surgery or if there is a limited range of motion due to other factors.

Custom Hip Kits Available
Please see your Sales Rep or call
Customer Service for details

STANDARD KIT

L3030 Standard Hip / Knee Kit

Each kit contains one each:

- P2226 26" aluminum reacher with magnetic tip
- L3008 Standard sock aid
- L3040 Long handle round bath sponge
- L3020 Wooden dressing stick
- L3023 18" long plastic shoehorn



Sock Aid



Round Sponge



Shoehorn



26" Reacher



Dressing Stick

DELUXE KITS


L3031 Deluxe Hip / Knee Kit

Each kit contains one each:

- P2232 32½" aluminum reacher with magnetic tip
- L3041 Deluxe hourglass bath sponge
- L3020 Wooden dressing stick
- L3009 Deluxe molded sock aid
- L3024 Deluxe 23" metal shoehorn



Sock Aid



Hourglass Sponge



Shoehorn



32½" Reacher



Dressing Stick

L3032 Deluxe Hip / Knee Kit

Each kit contains one each:

- P2232 32½" aluminum reacher with magnetic tip
- L3041 Deluxe hourglass bath sponge
- L3020 Wooden dressing stick
- L3010 Terry lined sock aid
- L3023 18" long plastic shoehorn



Sock Aid



Hourglass Sponge



Shoehorn



32½" Reacher



Dressing Stick

L3033 Deluxe Hip / Knee Kit

Each kit contains one each:

- P2232 32½" aluminum reacher with magnetic tip
- L3041 Deluxe hourglass bath sponge
- L3020 Wooden dressing stick
- L3009 Deluxe molded sock aid
- L3026BK 1 pair 32" black shoelaces
- L3026WH 1 pair 32" white shoelaces
- L3027 Jumbo shoehorn



Sock Aid



Hourglass Sponge



Shoehorn



32½" Reacher



Shoelaces



Dressing Stick

ESSENTIAL REACHERS

Looking for a reacher with incredible features?

Then the **EVERYDAY PICK™** should be your **First Choice!**



Easy!
Unique!



360° Rotating Head
that enables you to get to items
in hard to reach areas.


Pressure Sensitive Device
that prevents crushing even the most
delicate items.

Everyday Pick™ Reachers in Blue

- Sold one (1) per retail package

P2227B 26" Length

P2233B 32" Length



With Rotating Head

Deluxe Reacher / Grabber

- Rotating head to reach into narrow areas
- Head has both a magnet and built-in hanging device for storage
- Anodized aluminium
- Sold one (1) per retail package

P2224 26" Length

P2225 32½" Length



Raised magnetic post
to pick up/secure items



Aluminum Reacher

- Anodized aluminum
- Trigger activates strong plastic jaw with heavy duty nylon cord
- Raised post magnetic tip on end to pick up lightweight metallic objects and clothing
- Sold one (1) per retail package

P2226 26" Aluminum Reacher

P2232 32½" Aluminum Reacher



P2222 Folding Reacher

- 22" long
- Anodized aluminum
- Magnetic tip on end
- Trigger activates strong plastic jaw with heavy duty nylon cord
- Sold one (1) per retail package



Bath Safety Brochures help you show your customers all of their options with Essential Products!

For Endurance HD® Shower Benches, see page 49

SHOWER BENCHES/ STOOLS

Deluxe Tool Free Shower Benches

- With and without Back
- Full color retail box
- Tool free assembly
- Large, textured non-slip seat
- Height adjustable legs in 1" increments with drain holes
- Rust resistant aluminum frame

B3002-S Shower Bench

B3003-S Shower Bench with Back



SPECIFICATIONS:	B3002-S	B3003-S
Height Adjustments:	14¼" - 21¼"	14¼" - 21¼"
Seat Dimensions:	12" D x 19½" W	12" D x 19½" W
Overall Width:	19½"	19½"
Width at Bottom of Tips:	17½"	17½"
User Weight Limit:	300 lbs.	300 lbs.
Sold in Case Pack of:	2	2

Deluxe Shower Benches with Arms

- With and without Back
- Padded, removable arms
- Full color retail box or pre-assembled
- Large, textured non-slip seat
- Height adjustable legs in 1" increments with drain holes
- Rust resistant aluminum frame



SPECIFICATIONS: ALL MODELS

Height Adjustment:	16¼" - 21½"
Seat Dimensions:	13" D x 15¾" W
Overall Width:	20¾"
Width at Bottom of Tips:	17½"
Back Dimensions:	
B3011, B3013	16½" W x 6¾" H
User Weight Limit:	300 lbs.
Sold in Case Pack of:	
B3010, B3011	4
B3012, B3013	2

B3010 Shower Bench with Arms
without Back, Retail Box, Tool Free

B3011 Shower Bench with Arms
with Back, Retail Box, Tool Free

B3012 Shower Bench with Arms
without Back, Pre-Assembled

B3013 Shower Bench with Arms
with Back, Pre-Assembled



B3100 Round Bath Stool

- Tool free assembly
- Compact size will fit most tubs and showers, especially narrow ones found in RVs or manufactured homes
- Height adjustable legs in 1" increments with drain holes
- Round, textured non-slip seat
- Full color retail box
- Sold one (1) per retail box

B3100 SPECIFICATIONS

Height Adjustment: 17¼" - 21¼"
 Seat Dimensions: 12½" Diameter
 Overall Width: 13"
 Width at Bottom of Tips: 13"
 User Weight Limit: 300 lbs.
 Sold in Eaches

TUB TRANSFER BENCHES



B3005 Three Panel Bath Tub Transfer Bench

- White color
- Three panel seat w/drain holes
- Suction cups to secure unit to bathtub surface
- Tool free, snap together assembly

B3005 SPECIFICATIONS

Seat Height: 20" - 24"
 Seat Dimensions: 26¾" L x 16" W
 User Weight Limit: 250 lbs.
 Sold in Case Pack of: 2

HAND HELD SHOWER SPRAY

The **PERFECT** addition to Shower and Transfer Benches




B3501 Deluxe Hand Held Shower Standard Sprayer

- Easy to use on/off button for added control
- 84" reinforced hose
- Diverter valve
- Wall mount with screws
- Sold in eaches
- Full color retail box
- Sold one (1) per retail box

SHOWER AND TUB SAFETY ACCESSORIES

Deluxe Shower Mats

- 21" x 21"
- Ideal for walk-in showers
- Non skid surface
- Secure grip suction cups
- Sold one (1) per retail package



Unique Packaging
To Increase Sales!

- B3417B Blue Deluxe Shower Mat**
B3417C Cream Deluxe Shower Mat

AS34 Shower Mat Display

Everyone worries about bathroom **SLIP AND FALL**.

Presenting the first effective way to ensure that Shower Mats are highly visible to your customers.

- Display box is a compact 8" x 8"
- Contains five each B3417B Blue Shower Mats and four each B3417C Cream Shower Mats, 21" x 21" with center drain
- High visibility - good for retail shelf or counter



B3416C Deluxe Bath Safety Mat - Cream

- **Super thick - not the FLIMSY mats found in discount stores**
- Unique suction cups provide a secure grip
- Extra large size for safety - 41 1/4" L x 16" W
- Full color consumer packaging
- Latex free
- Sold one (1) per retail package



B3200 Essential Bathtub Safety Bar

- Adjustable clamp fits tub walls 3" to 7" wide
- Clamp has rubber pads to protect tub walls, ensure a firm grip
- Bar is made from 1 3/8" enamel coated steel
- Full color retail box
- Sold one (1) per retail box

SPECIFICATIONS

Height: 14"
Width: 6"
Color: White
Sold in Eaches



NOTE: B3200 should not be used on fiberglass or plastic tubs



Essential Chrome Plated Grab Bars
Professional Installation Required

- Chrome plated “knurled” finish for a secure grip
- Stainless steel hardware is included
- Diameter is 1" and mount is 1 1/2" from the wall
- Sold one (1) per retail box

B3212C	12" Long Bar
B3216C	16" Long Bar
B3218C	18" Long Bar
B3224C	24" Long Bar



Essential White Enamel Grab Bars
Professional Installation Required

- Rugged, enamel coated steel construction
- Stainless steel hardware is included
- Diameter is 1" and mount is 2 3/4" from the wall
- Sold one (1) per retail box

B3212	12" Long Bar
B3216	16" Long Bar
B3218	18" Long Bar
B3224	24" Long Bar

TOILET SAFETY RAILS



B5040 Essential Adjustable Toilet Safety Rails

- Anodized aluminum frame
- Full length waterfall arms
- Locking bracket attaches under toilet seat
- Textured plastic armrest
- Height and width adjustable
- Full color retail box

SPECIFICATIONS

Width between arms:	16 3/4" to 19 3/4"
Height Adjustment:	23 3/4" to 31 3/4" in 1" increments
User Weight Limit:	250 lbs.
Sold in Case Pack of:	4 pairs



B5041 Standing Adjustable Toilet Safety Rails

- Lightweight steel frame
- **Height adjustable**
- Padded arms for comfort
- Magazine rack included
- Easily assembled using wing nuts
- Pads on the bottom to help prevent sliding
- Full color retail box

SPECIFICATIONS

Total Width:	24 1/2"
Width Between Handle Bars:	20 3/4"
Height Adjustment:	24" to 29"
User Weight Limit:	300 lbs.
Sold in Case Pack of:	2 ea



TOILET SEAT RISERS AND RAISED TOILET SEATS

Toilet Seat Risers

- Enables users to use their existing toilet seat
- Includes 6" bolts and wing-nuts for easy installation
- Molded plastic design raises toilet seat 3½"
- Can be combined with Adjustable Toilet Safety Rails (B5040 & B5041)
- Full color retail box



B5080 Standard Size Bowl
B5081 Elongated Size Bowl

SPECIFICATIONS for:	B5080	B5081
Height:	3½"	3½"
Overall Width:	13½"	14"
Depth:	17½"	19½"
User Weight Limit:	300 lbs.	300 lbs.
Sold in Case Pack of:	4	4

Toilet Seat Risers with Removable Arms

- Padded, removable arms provide support while sitting down or getting up
- Enables users to use existing toilet seats
- Includes 6" bolts and wing-nuts for easy installation
- Molded plastic design raises toilet seat 3½"
- Full color retail box



B5082 Standard Size Bowl
B5083 Elongated Size Bowl

SPECIFICATIONS for:	B5082	B5083
Height:	3½"	3½"
Overall Width:	23"	22¾"
Width Between Arms:	19½"	19"
Depth:	17½"	19"
Arm Height:	9¾"	10"
User Weight Limit:	300 lbs.	300 lbs.
Sold in Case Pack of:	2	2

Hinged Toilet Seat Risers

- Unique hinged design makes cleaning under the riser easy
- Enables users to use existing toilet seats
- Includes 6" bolts and wing-nuts for easy installation
- Molded plastic design raises toilet seat 3½"
- Can be combined with Adjustable Toilet Safety Rails (B5040 & B5041)
- Full color retail box



B5084 Standard Size Bowl
B5085 Elongated Size Bowl

SPECIFICATIONS for:	B5084	B5085
Toilet Set Elevation:	3½"	3½"
Seat Dimensions:	17¼" x 13¾"	19¼" x 14"
Weight Capacity:	300 lbs.	300 lbs.
Sold in Case Pack of:	3	3



B5050 without Arms

Locking Raised Toilet Seats

- Raises toilet seat 5" for easier use
- Choice of styles with or without arms
- Fits most elongated and standard size toilets
- Adjustment knob allows it to lock securely on toilet bowl
- Full color retail box



B5051 with Arms

B5050 Without Arms
B5051 With Removable Arms

SPECIFICATIONS for:	B5050	B5051
Height:	5"	5"
Width:	17"	20"
Depth:	17½"	17½"
Arm Height:	—	5½"
Width Between Arms:	—	17½"
User Weight Limit:	300 lbs.	300 lbs.
Sold in Case Pack of:	2	2



B5060 Raised Toilet Seat for Standard Bowl

- Raises toilet seat 5" for easier use
- Rubber tabs on its bottom for added stability
- Tool free installation
- User weight capacity 300 lbs.
- Full color retail box

SPECIFICATIONS:	
Height:	5"
Width:	15"
Depth:	15½"
User Weight Limit:	300 lbs.
Sold in Case Pack of:	3




2"

Padded Toilet Seat Cushions

- Soft, padded seat attaches to toilet seat via hook and loop straps
- Covered with easy to clean vinyl
- Sold one (1) per retail package

B5070 2" Thick
B5071 4" Thick



4"

For Endurance® Heavy Duty Commodes, see page 49

ADJUSTABLE COMMODES

B4100 ESSENTIAL All-in-One Steel Commode

*Inexpensive —
But Not Cheap!*

- Removable back
- Covered armrests
- Consistent, gray powder coated finish
- Includes bucket with lid and splash guard

Use **T40100** Replacement Tips

SPECIFICATIONS

Height Adjustment:	16¼" - 22¼"
Overall Width:	23"
Width Between Arms:	18½"
Overall Depth:	18"
User Weight Limit:	300 lbs.
Sold in Case Pack of:	4



COMMUNE BUCKETS / SPLASH GUARDS

B5002 Buckets with Lid

B5003 Splash Guards for 3-in-1's

- Will fit most commodes / 3-in-1's
- Sold in Carton Lots of Twelve (12) ONLY



B5002

B5650 Replacement Buckets for B5600 H.D. Commode

- Sold in eaches



B5003

At Essential, we can help you create planograms to maximize your store space and sales!



AS67 Bath Safety 8 ft Planogram

ENDURANCE® HEIGHT ADJUSTABLE HAND BED RAILS AND BED ACCESSORIES

Give your customers the **Helping Hand** they need to get in and out of bed safely



P1410 Hand Bed Rail
P1410-P Hand Bed Rail with Pouch




- 20" wide, non-slip hand grip
- Tool free assembly
- Slides up to 36" under mattress with easy adjustments and has non-slip sleeves to keep it in place with an adjustable strap to secure to frame
- Fits most beds
- Made of chrome plated steel
- Sold one (1) per retail box

Both P1410 and P1411 rails include securing strap to attach rail to bed frame.

SPECIFICATIONS

Height Adjustment: 17" - 21"
 Sold in Eaches

Both P1410 and P1411 are for the exclusive use of adults as an assist device and not intended for use by or for children, infants or toddlers. Moreover, it is not a restraining device meant to keep a person in bed.



P1411 Hand Bed Rail with Floor Support
P1411-P Hand Bed Rail with Floor Support with Pouch



- 20¼" wide handle has ¼" thick foam tubing to provide a secure grip
- Slides up to 43" under mattress with easy adjustments and has non-slip sleeves to keep it in place with an adjustable strap to secure to frame
- Fits most beds
- Made of extra strong steel
- Sold one (1) per retail box

SPECIFICATIONS

Height Adjustment: 19¼" - 21¼"
 Leg Adjustment: 13¾" - 21"
 Sold in Eaches



P1409 Hand Bed Rail Pouch

- Three pockets will hold remotes, cell phones, reading materials, etc.
- Fits almost any bed rail
- Attaches via adjustable straps with snap buttons
- Sold one (1) per retail box



P1460 Home Bed Rails

- Will fit twin, full, queen and California king size beds via extension tubes included at **no extra charge**
- Constructed of 7/8" diameter steel tubing for strength
- Rails are fully adjustable with choice of two height adjustments (15¾" or 22¾") and three length adjustments (37¾", 43" or 56¾")
- Sold one (1) pair per retail box

P2600 Automatic Overbed Table

- Takes only fingertip pressure to height adjust
- Easy to assemble
- Melamine laminated top
- Four 2" swivel casters

SPECIFICATIONS

Table Top: 15¼" W x 30¼" L
 Height Adj: 27½" - 39¼"
 Sold in Eaches



P2601 Tilt Top Overbed Table

- Easy to assemble and adjust
- Laminated top
- Four 2" swivel casters
- Tilts for reading in bed

SPECIFICATIONS

Table Top: 16" W x 23¾" L
 Height Adj: 28" - 47"
 Sold in Eaches



P2400 Aluminum Blanket Support

- Height adjustable to 20"
- Fits under end of mattress
- Top is 20" wide, bottom is 15" wide
- Keeps blanket off feet and legs
- Sold one (1) per retail box



VINYL MATTRESS PROTECTORS

Contour Style



- High quality vinyl material
- Elasticized edges keep protector in place

For Mattresses: 36½" x 80½" x 6½"

C4001 12 Bagged Protectors/Case
C4100 Retail Package

Home Bed Sizes - Individually retail packed

C4300T	Twin Size	39" x 75" x 7½"
C4300F	Full Size	54" x 75" x 7½"
C4300Q	Queen Size	60" x 80" x 7½"
C4300K	King Size	78" x 80" x 7½"

Zippered Style



- High quality vinyl material
- Rust proof zipper
- Gives full mattress protection

For Mattresses: 36½" x 80½" x 6½"

C4002 12 Bagged Protectors/Case
C4200 Retail Package

Home Bed Sizes - Individually retail packed

C4400T	Twin Size	39" x 75" x 9½"
C4400F	Full Size	54" x 75" x 9½"
C4400Q	Queen Size	60" x 80" x 9½"
C4400K	King Size	78" x 80" x 9½"

Vinyl Pillow Protectors

- Zippered style for full protection
- Sold one (1) per retail package

C4800S	Standard Size	20" x 27"
C4800K	King Size	20" x 37"

HOSPITAL BED SHEETS



Cotton/Poly Pillow Case - White Muslin

- C3050** Retail Pack
- C3050B** Bulk Pack

PILLOW CASE SPECIFICATIONS

Material: 50/50 Cotton/Poly
 Size: 42" x 34"
 Retail Package: 1
 Bulk Pack: 12



Cotton/Poly Blend Sheets

- 55% Cotton/45% Poly

Contour Fitted - White

- C3051** Retail Pack
- C3051B** Bulk Pack

Top Sheets - White

- C3053** Retail Pack
- C3053B** Bulk Pack

COTTON/POLY SPECIFICATIONS

Contour Fitted:
 Size: 36" x 80" x 6"
 Retail Package: 1
 Bulk Pack: 12

Top Sheets:
 Size: 66" x 104"
 Retail Package: 1
 Bulk Pack: 12



Jersey Knit Sheets - White

- 55% cotton/45% poly blend
- Contour fitted sheet will fit all hospital mattresses
- Unique knit allows maximum air circulation for patient comfort
- Stretches to allow use of convoluted bed pad under sheet

- C3052** Retail Pack
- C3052B** Bulk Pack

KNIT SPECIFICATIONS

Size: 36" x 80" x 6"
 Retail Package: 1
 Bulk Pack: 12

HOSPITAL BEDDING SETS

Great delivery van item when doing bed set ups!



Sold one (1) per retail package (zippered blanket bag)

C3054 Pillow Protector Bed Set

- Includes flat and fitted cotton/poly hospital bed sheets, zippered vinyl pillow protector

C3055 Standard Hospital Bed Set

- Includes flat and fitted cotton/poly hospital bed sheets, pillow case

C3055K Hospital Bed Set, w/Jersey Knit Fitted Sheet

- Includes flat cotton/poly and fitted jersey knit sheets, pillow case



C3056 Deluxe Hospital Bed Set

- Includes flat and fitted cotton/poly hospital bed sheets, pillow case, vinyl zippered mattress protector, zippered vinyl pillow protector

COLORFUL HOSPITAL BED SHEETS AND BEDDING SETS



Sold one (1) per retail package

Fitted Hospital Bed Sheets

- C3057L** Light Blue
- C3057Y** Burgundy

Hospital Bed Bedding Sets

(Includes Fitted Sheet, Flat Sheet and Pillowcase)

- C3058L** Light Blue

ESSENTIAL BED WEDGES

Outstanding quality at a very attractive price!

All Essential Medical Bed Wedges are Fire Retardant!

Bed Wedges

- Use for head, foot or leg elevation
- Blue cotton/poly blend zippered cover
- 24" x 24"

SPECIFICATIONS

Sold in Case Pack of: 2

F1575-2	7½" Deep, with 17° incline
F1510	10" Deep, with 22° incline
F1512	12" Deep, with 25° incline



Zippered Replacement Covers

F1575-24-Cover	For 7½" Deep Wedges
F1510-24-Cover	For 10" Deep Wedges
F1512-24-Cover	For 12" Deep Wedges

Folding Bed Wedges

- Gives the same great support as our standard bed wedges
- Folding feature is ideal for users who travel
- Takes up less room on retail shelves
- Smaller carton means lower shipping cost
- 24" x 24"
- Packed one (1) per carton

F1575F	7½" Deep, with 17° incline
F1510F	10" Deep, with 22° incline
F1512F	12" Deep, with 25° incline



Zippered Replacement Covers

F1575F-Cover	Cover for 7½" Folding Wedge
F1510F-Cover	Cover for 10" Folding Wedge
F1512F-Cover	Cover for 12" Folding Wedge



CONVOLUTED FOAM BED PADS



All Essential foam bed pads are **fire retardant**, **latex free** and meet the requirements of **California Technical Bulletin #117**.

Hospital Bed Size

- Compact rolls sealed in plastic w/sealed ends to preserve shelf life
- 33" W x 76" L
- Sold in case lots ONLY

F7220	2" Thick	12 per case
F7201	3" Thick	8 per case
F7243	4" Thick	6 per case



HOME BED SIZE

Home Bed Size

- Packed one (1) per carton in **ZIPPERED BLANKET BAGS**

	Twin Size	Full Size	Queen Size	King Size
2"	F9010E	F9001E	F9002E	F9003E
3"	F9011E	F9004E	F9005E	F9006E
4"	F9012E	F9007E	F9008E	F9009E

Home Bed Size Dimensions

Twin:	36" W x 76" L
Full:	50" W x 76" L
Queen:	56" W x 76" L
King:	72" W x 76" L

Memory P.f.[®]

High quality memory foam designed to give your customers a

*Personal Fit (P.f.)
of Comfort!*



Heavyweight 6 lb. Density — not the lighter version sold by many others

All Memory P.f.[®] products are latex free!



N3001 Memory Foam PF[®] Knee Separator

- Luxurious zippered blue velveteen cover
- Includes strap w/hook and loop closure
- Dimensions: 10" L x 8" W
- Sold one (1) per retail box



N3002 Memory Foam PF[®] Round Travel Pillow

- Crescent design to "cradle" the neck
- Luxurious zippered blue velveteen cover
- Sold one (1) per retail box



N3003 Memory Foam PF[®] Full Size Contour Sleep Pillow

- Luxurious zippered velveteen cover
- Dimensions: 20" L x 12" W
- Sold one (1) per retail box




N3004 Memory Foam PF[®] Lumbar Support Cushion

- Attached strap holds cushion in place
- Dimensions: 13¼" L x 16" W
- Sold one (1) per retail box

Memory P.f. **COMFORT CUSHIONS**

Comfort cushions are designed to relieve pressure, provide a comfortable seat for individuals sitting in chairs, riding in cars, etc. **Great for use in Scooters!**

• All *Memory P.f.* Comfort Cushions are sold one per retail box. •



N3008 Molded Memory Foam P.F.® Wedge "Comfort" Cushion

- Slope provides added user comfort, especially while riding in cars
- Luxurious zippered blue velveteen cover
- Dimensions: 18" x 16" x 3" x 1½"



N3009 Memory Foam PF® Sculptured "Comfort" Cushion

- Sculpted design is anatomically correct and provides maximum comfort
- Luxurious zippered blue velveteen cover
- Dimensions: 17" x 2" x 15"



N3010 Memory Foam P.F.® Sculptured "Comfort" Cushion w/ Rear Cut Out

- Sculpted design is anatomically correct and provides maximum comfort
- Rear cut out relieves pressure on the coccyx area
- Luxurious zippered blue velveteen cover
- Dimensions: 17" x 2½" x 14"

N2010 The Cushion™ by Essential Medical

The only cushion you'll ever need!

- Can be used as a Support Cushion, Coccyx Cushion or Invalid Ring
- Molded foam helps the user keep proper posture when sitting
- Sold one (1) per retail package



The Cushion™
by Essential Medical Supply, Inc.

N2010-Deal

- Buy six (6) cushions and get a **FREE** "Try Me" Display sample



WHEELCHAIR CUSHIONS



Coccyx Cushions

- Rear cutout relieves pressure on coccyx area, providing stability and comfort
- Navy zippered cover
- Sold one (1) per retail package

N1001 16" x 16" x 3"

N1002 18" x 16" x 3"



N1002P Plaid Coccyx Cushions with Removable Masonite Insert

- Removable masonite insert provides firm support
- Firm 3" foam base
- Attractive plaid zippered removable cover
- Dimensions: 18" x 16" x 3"
- Sold one (1) per retail package

Plaid cover may vary from photo



N1003 Sloping Coccyx Cushion

- Slope helps keep user comfortably positioned
- Rear cutout relieves pressure on coccyx area, providing stability and comfort
- Black zippered removable cover
- Dimensions: 18" x 16" x 3" x 1½"
- Sold one (1) per retail package



N1005 Sloping Bucket Seat Cushion

- Gentle slope helps keep user back and tailbone elevated in bucket seats
- Washable, removable cover
- Dimensions: 15" x 13"
- Sold one (1) per retail package



Gel Cushions

- Gel bladder with aqueous, non-toxic gel enclosed in fire retardant polyurethane foam
- Reduces pressure, supports ischial tuberosities
- 250 lb. user weight capacity
- Sold one (1) per retail package

D4001	18" x 16" x 2"	Medicare Code E2601
D4002	18" x 16" x 3"	Medicare Code E2603
D4003	20" x 16" x 3"	Medicare Code E2603
D4006	16" x 16" x 2"	Medicare Code E2601
D4007	16" x 16" x 3"	Medicare Code E2603



Fleece Covered Gel Wheelchair Cushions

- Gel bladder with aqueous, non-toxic gel enclosed in fire retardant polyurethane foam
- Reduces pressure, supports ischial tuberosities
- Nylon zippered multi cover, fluid resistant
- Removable, washable cover
- 250 lb. user weight capacity
- Sold one (1) per retail package


D4100	18" x 16" x 3"
D4101	18" x 16" x 3" with Coccyx Cutout



Wedged Cushions

- Gentle slope of cushion keeps the user back in the chair
- Navy zippered cover
- Sold one (1) per retail package


N2001	18" x 16" x 3" x 1½"
N2002	18" x 16" x 4" x 1½"



Rehab 1™ Cushions

- Basic rehab cushion made of high density polyurethane foam
- Provides firm seating and pressure relief
- Black zippered cotton cover
- Sold one (1) per retail package

N1102	18" x 16" x 2"
N1103	18" x 16" x 3"
N1104	18" x 16" x 4"



Fleece Covered Coccyx Cushions

- Excellent foam density
- Soft fleece cover for improved comfort
- Zippered, removable cover
- Rear cutout relieves pressure on coccyx area
- Sold one (1) per retail package

N1006 16" x 16" x 3"

N1008 18" x 16" x 3"



Fleece Covered Cushions

- Excellent foam density
- Soft fleece cover for improved comfort
- Zippered, removable cover
- Sold one (1) per retail package

N1106 16" x 16" x 3"

N1108 18" x 16" x 3"



N3007 Memory Foam PF® Coccyx Cushion

- Heavyweight, 6 lb. density for maximum support
- Rear cutout relieves pressure on coccyx area
- Zippered removable cover
- Dimensions: 18" x 16" x 3"
- Sold one (1) per retail package

Memory P.f.®



ESSENTIAL DONUT RINGS (INVALID)



CONVOLUTED FOAM STYLE - Navy Cover

N8000 18" Donut with Removable Cover, 4" thick

N8001 16" Donut with Removable Cover, 4" thick

MOLDED STYLE

- High density molded construction
- Choice of sizes and removable covers
- Sold one (1) per retail package



Plaid cover may vary from photo

14" Molded		16" Molded		18" Molded	
		N8002	Blue Cover	N8008	Blue Cover
N8004N	Navy Cover	N8002N	Navy Cover	N8008N	Navy Cover
		N8002P	Plaid Cover	N8008P	Plaid Cover

Designer Series Comfort Rings

- Modern designer cover is removable and washable
- Foam provides a great balance between support and comfort
- Dimensions: 16" W x 12½" D x 3" H
- Sold one (1) per retail package

N8060N Navy with Polka Dots
N8060P Pink Damask



N8100 Oversized Gel Cushion with Luxurious Velour Cover

- 4" thick for extra comfort
- Includes non aqueous gel bladder
- Reduces pressure
- Dimensions: 19" W x 16" D x 4" H
- Sold one (1) per retail package



N8101 Oversized Gel Cushion with Luxurious Velour Cover & Fleece on Top

- 4¼" thick for extra comfort
- Includes non aqueous gel bladder
- Reduces pressure
- Dimensions: 19" W x 15¾" D x 4¼" H
- Sold one (1) per retail package



INFLATABLE RING

D1001 Inflatable Vinyl Invalid Ring

- Heavy gauge vinyl
- Retail packaged
- 16" Diameter
- Sold one (1) per retail package



KNEE SEPARATORS / LEG RESTS



Knee Separators

- Includes strap with hook and loop closure
- Unique hour glass shape
- Provides ideal alignment for spine, knees and hips
- Choice of washable covers
- Dimensions: 3-6" W x 10" D x 8" H
- Sold one (1) per retail package

N6400 **Blue Cotton/Poly Cover**
N6401 **White Terry Cloth Cover**



Elevating Knee Rests

- Helps reduce back pain
- Improves circulation
- Blue, removable cover
- Sold one (1) per retail package

N6402 **17" x 10" x 7"**
N6403 **28" x 10" x 7"**



N6500 Elevating Leg Rest

- Eases stress on the spine, lower back
- Allows users to read, watch TV comfortably while in bed
- Washable removable cover
- Dimensions: 20" W x 26" D x 8" H
- Sold one (1) per retail package

LUMBAR SUPPORT CUSHIONS

Molded Lumbar Cushions



F1412 **Navy**
F1412B **Burgundy**
F1412BK **Black**
F1412G **Gray**



- Molded, high density foam
- Zippered removable cover
- Attached strap holds cushion in place
- Dimensions: 14" W x 13" H
- Sold one (1) per retail package



F1413BK Black Lumbar Cushion for Bucket Seats

- Perfect for car seats or use in home chairs
- Mesh cover for aeration
- Attached, adjustable strap keeps the cushion in place
- Dimensions: 14½" W x 13" H
- Sold one (1) per retail package

Designer Series Back Cushions

- Luxurious removable cover
- Molded high density foam
- Attached strap holds cushion in place
- Dimensions: 14" W x 13" H
- Sold one (1) per retail package

F1415N Navy with Polka Dots
F1415P Pink Damask



F1418BK Black Deluxe Lumbar Cushion

- Extra wide for added comfort and support
- Molded high density foam
- Zippered removable cover
- Attached strap holds cushion in place
- Dimensions: 17³/₄" W x 13³/₄" H
- Sold one (1) per retail package



Lumbar Rolls

- Designed to provide proper lumbar curvature while sitting or driving
- Attached strap holds the roll in place
- Sold one (1) per retail package

N6002 Half Roll 10" x 3"
N6003 Full Roll 12" x 4"



N1004 Lumbar Wheelchair Cushion

- Outstanding reimbursement level
- Straps connect it to chair back
- 16" W x 18" H x 2³/₄" D
- Sold one (1) per retail package

Medicare Code E2611



NECK SUPPORT CUSHIONS

N5001 Neck Yoke

- Provides extra firm neck support and excellent comfort while sleeping
- Blue cotton cover
- Dimensions: 14" W x 6" D x 6" H
- Sold one (1) per retail package



NECK SUPPORT CUSHIONS



N5002 Crescent Pillow

- Supports neck and head while napping or reading
- Perfect travel size pillow
- Blue cotton cover
- Sold one (1) per retail package



Fiber Filled Round Cervical Pillow (Jackson Style)

- Promotes proper cervical alignment
- Provides neck support while sleeping on the back
- Choice of removable covers
- Dimensions: 6" W x 18" H
- Sold one (1) per retail package

N5003 White Cotton Cover

N5005 Blue Satin Cover



Foam Cervical Rolls

- Made of firm but comfortable foam
- Washable, removable cover
- Use with existing pillow or alone
- Can be used under knees, back or neck
- Length: 19"
- Sold one (1) per retail package

N5006 3½" Diameter

N5007 5" Diameter

N5008 7" Diameter

E-Z SLEEP™ PILLOWS



Our Best Selling Sleep Pillow!

N7104 The Eclipse™ Pillow

- Fiber filled for comfort
- Our most versatile pillow
- Unique contour and center design provides all position relief
- The pillow that eclipses all others for stress and muscle pain relief
- Dimensions: 17" W x 20" L
- Sold one (1) per retail package



N7103 Butterfly Pillow

- Fiber filled for soft, comfortable support
- Unique butterfly design prevents neck strain
- Reduces stiff neck, shoulder and back pain
- Recommended for either back or side sleeping
- Zippered blue satin cover
- Dimensions: 17" W x 22" L
- Sold one (1) per retail package

ENDURANCE® ALTERNATING PRESSURE PUMPS AND PADS

D9000C Deluxe Pump and Standard Pad Combo

- Includes both D9000 Pump (*described below*) and D8800 Pad (*described below*)
- Sold one (1) set per retail box

D9000 Pump

- Adjustable pressure control
- Two hanging hooks to mount unit on bed frame or prop upright
- Lighted on/off switch
- Grounded cord with three prong plug
- Suction feet for quiet operation and stability
- 1 year limited warranty



D8800 Endurance® Standard Bubble Pad

- Inflated size - 78" L x 34" W
- 2 ½" bubble height meets requirements
- High grade, heavy duty vinyl to resist tearing, punctures
- Sealed edges to prevent bacterial buildup
- Bubble style cells provide maximum loft, improved user comfort
- Pad includes connecting tubing
- Fire retardant and not made with natural rubber latex
- 250 lb. weight limit
- Sold one (1) per retail package

Endurance®
Bubble Pads
will fit most brands
of APP Pumps



D8800

SHEEPETTE® BED PADS

SHEEPETTE® Synthetic Sheepskin Bed Pads

- D5001** **24" x 30"**
- D5002** **30" x 40"**
- D5003** **30" x 60"**
- D5004** **36" x 80"**

- Thick, plush Sheepette® material provides superior user comfort
- Dense, heavyweight material
- Machine washable and dryable
- **D5004 has elastic securing straps on all four corners**
- Sold one (1) per retail package



SHEEPETTE® PROTECTORS AND ACCESSORIES



Heel Protectors

- Thick, plush Sheepette® material provides superior user comfort
- Retail packs are easy to display
- Sold one (1) pair per retail package

D3001 White Polyester

D5005 Sheepette® Synthetic Sheepskin



Elbow Protectors

- Thick, plush Sheepette® material provides superior user comfort
- Retail packs are easy to display
- Sold one (1) pair per retail package

D3002 White Polyester

D5006 Sheepette® Synthetic Sheepskin



D6000 Fiber Filled Heel Protectors

- Adjustable strap so one size fits most
- Unique zippered carrying package
- Sold one (1) pair per retail package



D3004 Sheepette® Wheelchair Armrest Pads

- Thick, plush Sheepette® material provides superior user comfort
- Can be easily cut to fit desk length arms
- 12" L x 9" W
- Sold one (1) pair per retail package



D3005 Sheepette® Wheelchair Seat & Back Pad

- Thick, plush Sheepette® material provides superior user comfort
- 18" W x 40½" L
- Sold one (1) per retail package

Our Heavy Duty products are engineered for maximum strength and durability



Heavy Duty products that give you and your customers the quality you need at a price you can afford.

ENDURANCE HD HEAVY DUTY BATH SAFETY PRODUCTS

B3007 Heavy Duty Shower Bench

- Supports users up to 500 lbs
- Double cross brace for added strength
- Tool free assembly
- White seat

B3007 SPECIFICATIONS

Seat Size: 20" W x 12" D
 Height Adjustable: 16¼" - 20¼"
 User Weight Limit: 500 lbs.
 Sold in Eaches



B3008 Heavy Duty Shower Bench with Back

- Supports users up to 500 lbs
- Double cross brace for added strength
- Tool free assembly
- White seat and back

B3008 SPECIFICATIONS

Seat Size: 20" W x 12" D
 Height Adjustable: 16¼" - 20¼"
 Back Size: 15½" W x 9" H
 User Weight Limit: 500 lbs.
 Sold in Eaches



B5600 Heavy Duty Commode

- Supports users up to 500 lbs
- Large platform with 8" x 10" opening
- Padded arms

Use **B5650** Replacement Buckets

B5600 SPECIFICATIONS

Height Adjustable: 16" - 20¼"
 Overall Width: 29¼"
 Width Between Arms: 27"
 Overall Depth: 23¼" - 24⅞"
 User Weight Limit: 500 lbs.
 Sold in Case Pack of: 2





HEAVY DUTY CANES



W1305

W1306

QUAD CANES

- Supports users up to 500 lbs
- Small base weighs only 3 lbs.
- Large base weighs only 4 lbs.

Black Color

W1305 Small Base, 8" x 6"
Use T50012BL Replacement Tips

W1306 Large Base, 12" x 8"
Use T50058BL Replacement Tips

W1305 & W1306 SPECIFICATIONS

Height Adjustment: 29" - 38"
User Weight Limit: 500 lbs.
Sold in Eaches



W1307 Offset Cane - Black

- Supports users up to 500 lbs
- Made from steel, but only weighs 3 lbs.

W1307 SPECIFICATIONS

Height Adjustment: 29^{3/8}" - 37^{3/8}"
User Weight Limit: 500 lbs.
Sold in Eaches



P2750 Heavy Duty Foot Stool

- Supports users up to 500 lbs
- Rubber surface to help prevent slipping
- Double cross brace for added support

P2750 SPECIFICATIONS

Height - 9" Length - 14^{1/4}" Width - 11^{1/4}"
User Weight Limit: 500 lbs.
Sold in Eaches

ENDURANCE^{HD}
BUILT TO LAST!

W1802B-14 Four Wheel Walker

- Supports users up to 500 lbs
- Soft, extra large seat
- Carrying pouch under padded seat
- Dual front cross bar for extra support
- Packed one (1) per case

W1802B-14 SPECIFICATIONS

Overall Width:	26¾"
Seat Height From Floor:	22"
Seat Size:	18½" W x 14¼" D
Width Between Handles:	22"
Height Adjustment:	31¾" - 36"
Wheels:	7½" x 1½"
Weight:	31 lbs.
User Weight Limit:	500 lbs.
Sold in Eaches	



ENDURANCE^{HD}
BUILT TO LAST!

W1250-C Heavy Duty Folding Walker

- Supports users up to 500 lbs
- Extra wide, 20" inside width
- Tab style double button opening/closing mechanism

W1255F-C 5" Fixed Wheels
for W1250 H.D. Walker

W1250-C SPECIFICATIONS

Width:	26"
Inside Width:	20"
Weight:	6 lbs. 9 oz.
Tubing Size:	1" aluminum
Height Adjustment:	32½" - 42½"
User Weight Limit:	500 lbs.
Sold in Case Pack of:	2



Quik-Sorb[™]

INCONTINENT PRODUCTS

The comfort of cloth . . . with the economy of a reusable!

QUIK-SORB[™] REUSABLE UNDERPADS

Available in **Newly Designed** Retail Packages or Bulk Packs of three or twelve.
Choice of Birdseye 100% cotton or brushed polyester quilted top sheets

- Quilted top sheets with waterproof back
- Heavyweight 8 oz. liquidfill absorbent center
- Machine washable/dryable up to 300 times

**Quik-Sorb[™]
Underpads
Are Not Made
With Natural
Rubber Latex**

PACKAGING SPECIFICATIONS:

All Quik-Sorb[™] Underpads come packed in one of three quantities:

Retail Package: 1 | Set Pack: 3 | Bulk Pack: 12

Please note that all B-3 Set Packs and B Bulk Packs are not retail packaged. Pads that are retail packaged are noted as Retail.

Birdseye - Our Top Seller!




BIRDSEYE 100% Cotton Quilted Top Sheet

17" x 24" Chair Pad

- C2000** Retail Package
- C2000B-3** Set Pack
- C2000B** Bulk Pack

24" x 35" Bed /Sofa Pad

- C2001** Retail Package
- C2001B-3** Set Pack
- C2001B** Bulk Pack



34" x 35" Bed/Sofa Pad

- C2002** Retail Package
- C2002B-3** Set Pack
- C2002B** Bulk Pack

C2002
Bed/Sofa Pad

36" x 80" Bed Pad

- C2005** Retail Package
- C2005B-3** Set Pack
- C2005B** Bulk Pack

36" x 54" Bed/Sofa Pad

- C2006** Retail Package
- C2006B-3** Set Pack
- C2006B** Bulk Pack



Detail picture of C2005 and C2009
elastic securing straps

36" x 72" Twin Bed Pad

- C2009** Retail Package
- C2009B-3** Set Pack
- C2009B** Bulk Pack

Brushed Polyester Quilted Top Sheet

24" x 35" Bed/Sofa Pad

- C2003** Retail Package
- C2003B-3** Set Pack
- C2003B** Bulk Pack

34" x 35" Bed/Sofa Pad

- C2004** Retail Package
- C2004B-3** Set Pack
- C2004B** Bulk Pack

QUIK-SORB™ Plaid Underpads

- Plaid top sheet is designed to **hide** stains
- Quilted top sheet
- Heavyweight liquid fill center
- Waterproof back
- Machine washable, dryable



Plaid pattern may vary from photo

18" x 24" Chair Pad

C2010 Retail Package
C2010-3 Set Pack
C2010-B Bulk Pack

24" x 36" Bed/Sofa Pad

C2011 Retail Package
C2011B-3 Set Pack
C2011B Bulk Pack

34" x 36" Bed Pad

C2012 Retail Package
C2012B-3 Set Pack
C2012B Bulk Pack

QUIK-SORB™ DELUXE Reusable Underpad with Tucks

- Quilted Birdseye 100% cotton top sheet
- Heavyweight 8 oz. liquid fill absorbent center
- Waterproof back
- 36" x 36" absorbent center with two 18" x 36" flaps to secure under mattress



Ideal for restless sleepers!

36" x 72" Bed Pad

C2008 Retail Package
C2008B-3 Set Pack
C2008B Bulk Pack

ESSENTIAL Reusable Underpads

- Twill top sheet (50% cotton/50% polyester)
- Waterproof back
- Rounded corners provide extra strength and inhibit fluid runoff
- Sold in cases of 24

30" x 36" Pad

C2300B Bulk Pack



QUIK-SORB™ Reusable Underpad w/Straps

- Brushed polyester pad with heavy duty 5" long straps
- Straps help caregiver position the pad underneath the user
- Heavyweight, 8 oz. liquid fill absorbent center

34" x 35" Bed Pad

C2400 Retail Package
C2400B-3 Set Pack
C2400B Bulk Pack




QUIK-SORB™ ECONOMY Reusable Underpad

- Percale top sheet for softness and air circulation
- Liquid fill absorbent center
- Waterproof vinyl back

34" x 36" Bed Pad


C2007 Retail Package
C2007B-3 Set Pack
C2007B Bulk Pack





QUIK-SORB™ Furniture Protectors

- Perfect for all types of furniture
- Leak resistant back helps protect furniture
- 20" x 20"
- Machine washable and dryable
- Sold one (1) per retail package



C2500-DEAL

- Contains 3 each Chocolate and Tan, 2 each Blue, Green and Maroon (12 total) plus **FREE** display pad

QUIK-SORB™ REUSABLE INCONTINENT PANTS

Made from soft, medical vinyl for improved comfort, less noise when worn



- Retail package for easy display
- Pull on or snap closure styles
- Completely waterproof
- Machine washable
- Packed one (1) per retail package

Pull On Style - Elasticized waist and legs

C6000S	Small	(22" - 28")
C6000M	Medium	(30" - 36")
C6000L	Large	(38" - 44")
C6000XL	X-Large	(46" - 52")



Snap Closure Style - Elasticized waist and legs

C6001S	Small	(22" - 28")
C6001M	Medium	(30" - 36")
C6001L	Large	(38" - 44")
C6001XL	X-Large	(46" - 52")



Quik-Sorb™ Brochures help you show your customers all of their options with Essential Products!

TENS SYSTEMS

All Essential Electrodes are not made with natural rubber latex

All ESSENTIAL TENS Lead Wires Meet 62FR25497 Performance Standards

S1000 Essential TENS System

- Constant current unit
- Dual isolated channels
- Pulse amplitude of 0 to 80 mA
- Pulse frequency 2 to 150, adjustable
- Pulse width 60 to 250, adjustable
- Set contains TENS unit, lead wires, 9V battery, electrodes and carrying case
- One year warranty
- Sold one (1) per retail package



A quality unit at a great price! NOT a cheap throw away!

S2000 Digital TENS Unit

- Digital read out
- Patient compliance timer
- Dual, isolated channels
- Five mode stimulation
- Pulse adjustable 2 - 150Hz
- Pulse width 50 - 300µs
- Set contains TENS unit, lead wires, 9V battery, electrodes and carrying case
- One year warranty
- Sold one (1) per retail package



TENS ACCESSORIES / ELECTRODES / GEL

S1000E TENS Replacement Electrode Set

- Sold four (4) electrodes per retail package

S1001E Butterfly Electrode

- 4" x 6"
- Ideal for backs and other large areas
- Sold one (1) per retail package



S1000L 43" TENS Lead Wires

- Sold two (2) Lead Wires per retail package

S1000G Parker Conductive Gel

- 2 oz. Tube
- Sold twelve (12) per box



PHYSICAL AND OCCUPATIONAL THERAPY PRODUCTS



P3000 Essential Pedal Exerciser

- Easy tool free setup
- Sleek black color
- Deluxe non-slip tips for stability
- Adjustable tension knob
- Sold one (1) per retail box



P3100



Folding Pedal Exerciser

- Takes up less shelf space
- Snaps in place - totally tool free!
- Easy to transport
- Sold one (1) per retail box

**P3100
P3101**

**Folding Pedal Exerciser
Folding Pedal Exerciser
with Tracking Monitor**



P3101
with Tracking Monitor

Monitor keeps track of:
workout time, cycles,
approximate calories burned



P2100 Overdoor Traction Set

Includes: 20 lb. water bag, overdoor traction clamp, 12" spreader bar, "S" Hook, traction cord and adjustable hook and loop closure head halter that fits most.

- Latex free
- Sold one (1) set per retail box



P1102 Exercise Pulley Set

- Can be used in a sitting or standing position
- Helps restore shoulder, arm, and upper torso mobility
- Installs easily over any door
- Sold one (1) set per retail box



GAIT BELTS

Standard Gait Belts

- Cotton blend
- Choice of 54" or 72" length
- Choice of plastic or metal buckles
- Sold one (1) per retail box



54" LENGTH

- P2501 with Metal Buckle**
- P2503 with Plastic Buckle**

72" LENGTH

- P2502 with Metal Buckle**
- P2504 with Plastic Buckle**

P2500 Deluxe Gait Belt with Plastic Buckle

- Cotton/polyester blend
- Helps stabilize users who are having difficulty ambulating
- **Has hand holds on the belt for caregiver assistance**
- Fits waist sizes 30" - 45"
- Sold one (1) per retail box



L2020 Deluxe Transfer/Gait Belt

- Padded for extra comfort
- Has four webbed hand grips for safety
- Total length is 56½", width is 5½"
- Fits waist sizes 28" to 42"
- Sold one (1) per retail package



REHAB SQUEEZE BALLS

- Easy to open/close clamshell package prevents the package from becoming *shop worn*
- Available in two styles and three densities
- Latex free and lightly powdered with talc
- Sold one (1) per retail package



DIMPLED STYLE

- P2010-B Medium Density**
- P2010-P Soft Density**
- P2010-R Hard Density**

EGG SHAPED

- P2020-B**
- P2020-P**

Star Shaped Hand Exerciser

- Unique shape massages and exercises fingers, hand and forearm
- Easy to open/close clamshell package prevents the package from becoming *shop worn*
- Great impulse item for your checkout counter
- Latex free and lightly powdered with talc
- Available in two densities: Soft and Medium
- Sold one (1) per retail package



- P2030-P Soft Density**
- P2030-B Medium Density**

AS57 Counter Top Squeeze Ball, Egg and Star Display

- Compact display takes up very little room and will maximize your sales
- Contains a mix of 30 squeeze balls, eggs and stars

Display Rack Included at No Charge



For Endurance® Heavy Duty Foot Stools, see page 50



CHROME PLATED SAFETY FOOT STOOLS



P2700 Safety Foot Stool

- Chrome plated steel for strength and durability
- Rubberized top for added safety

Use **35105** Replacement Tips

P2700 SPECIFICATIONS

Height:	9"
Length:	14¼"
Width:	11¼"
Weight Capacity:	350 lbs.
Sold in Case Pack of:	1



P2701 Safety Foot Stool with Handrail

- Comes with handrail
- Chrome plated steel for strength and durability
- Rubberized top for added safety

Use **35105** Replacement Tips

P2701 SPECIFICATIONS

Height:	9"
Length:	14¼"
Width:	11¼"
Weight Capacity:	350 lbs.
Height of Handrail:	33"
Sold in Case Pack of:	2



L3045 Deluxe Shampoo Basin

- Inflatable heavy duty vinyl with two air chambers
- Long, attached drainage hose
- Extra large size: 24" L x 28" W x 6" H
- Sold one (1) per retail package



P3001 Deluxe Swivel Seat Cushion

- Helps reduce back and hip strain
 - When getting in and out of vehicles
 - When rising from chairs at home
- Large size: 15¼" x 2"
- 300 lb. weight capacity
- Sold one (1) per retail box

PATIENT BIBS / CLOTHING PROTECTORS

Choice of four unique styles

Standard Bib

- Soft, blue terrycloth bib with hook and loop closure, 18" x 30"
- 100% cotton

C3042 Retail Package

C3042B Bulk Pack

Packaging

Retail Package: 1
Bulk Pack: 12



Fashion Clothing Protector with Crumb Catcher

- Absorbent front with wet proof back
- Covers upper body and lap for full protection
- Crumb catcher easily snaps into place

C3044 Retail Package

C3044B Bulk Pack

Packaging

Retail Package: 1
Bulk Pack: 12



Deluxe Bib

- Bright, heavyweight plaid front with hook and loop closure, 18" x 30"
- Waterproof vinyl back

C3045 Retail Package

C3045B Bulk Pack

Packaging

Retail Package: 1
Bulk Pack: 12



L5050 Clothing Protector With Adjustable Crumb Catcher

- Covers upper body and lap for full protection
- Crumb catcher can be adjusted
- Easy to wipe clean
- Sold one (1) per retail package



Essential Medical Supply, Inc. provides patient gowns in **resealable style packaging** that allow your customer to examine the gown without damaging the packaging. You'll never again have **shop worn** product that doesn't sell due to customer handling.

PATIENT GOWNS



C3009

C3010

Patient Gowns

- Made from high quality print cloth material
- Traditional center tie style
- Roomy raglan sleeves
- Machine washable/dryable

Bulk Pack: All available in bulk packs of 3 or 12 gowns per poly bag. Add B-3 to code number for three (3) pack (e.g. C3009B-3) or B only for twelve (12) pack (e.g. C3009B). Bulk Packs are NOT retail packaged.

Prints

C3009	Print on Light Blue
C3010	Blue Print on White

Length

36"
41"

Pastels / Solids

C3012	Light Blue
C3013	Pink
C3014	Yellow
C3015	Mint

41"
41"
41"
41"



C3012

C3013

C3014

C3015

SPECIFICATIONS

Material:	50/50 Cotton/Polyester
Retail Package:	1
Set Pack:	3
Dozen Pack:	12

DELUXE Patient Gowns

- Made from the finest quality broadcloth
- Overlap style with snaps provides full back coverage
- Roomy raglan sleeves
- Machine washable/dryable

Bulk Pack: All available in bulk packs of 3 or 12 gowns per poly bag. Add B-3 to code number for three (3) pack (e.g. C3022B-3) or B only for twelve (12) pack (e.g. C3022B). Bulk Packs are NOT retail packaged.



C3022

C3023

Broadcloth Gowns

C3022	Blue w/snap closures
C3023	Print on White w/snap closures

Length

41"
40"

SPECIFICATIONS

Material:	50/50 Cotton/Polyester
Retail Package:	1
Set Pack:	3
Dozen Pack:	12

SOFT N' COZY™ Poly Knit Flannel Gown

- Luxurious brushed knit polyester for comfort and softness
- The warmth of flannel without the weight
- Overlap back with snap closures for full modesty
- One size fits most

Half Sleeves	Color	Length
C3072	Blue	43"

SPECIFICATIONS

Material: Knit Polyester
Retail Package: 1



C3072

The warmth of flannel . . .
The luxurious feel of knit polyester

C3200 King/Queen Size Patient Gowns

- Fits up to 3XL
- Unisex print
- Generous back overlap and elbow length sleeves for maximum modesty
- Tie back
- 42" long

SPECIFICATIONS

Material: Cotton/Polyester
Retail Package: 1



C3300 IV Patient Gowns

- Snap closures on shoulder and sleeves
- 8" overlap back for complete modesty
- Blue material
- Tie back
- One size fits most
- 41" long

SPECIFICATIONS

Material: Cotton/Polyester
Retail Package: 1



MEDICAL PLASTIC WARE

Choose high impact consumer packaging,
individual units sealed in plastic or bulk pack in cases
Color and style of products may vary from photos



Deluxe Pontoon Style Bed Pan

- Contoured for patient comfort
- Made of easy to clean, strong polypropylene

C1101 Retail Package
C1001 Sealed in Plastic Bag
C1101B Bulk Pack

PACKAGING	
Retail Boxed:	1
Plastic Bag:	1
Bulk Pack:	20



Male Urinal w/Cover

- One (1) quart - 960 cc capacity
- Wide base design resists tipping

C1102 Retail Package
C1002 Sealed in Plastic Bag
C1102B Bulk Pack

PACKAGING	
Retail Boxed:	1
Plastic Bag:	1
Bulk Pack:	50



Fracture Bed Pan

- Wedge shape is perfect for bariatric or immobile patients
- Can also be used as a female urinal

C1103 Retail Package
C1103B Bulk Pack

PACKAGING	
Retail Boxed:	1
Bulk Pack:	24



Sitz Bath

- Use for treatment of hemorrhoids and other anal/genital problems
- Fits any standard commode

C1104 Retail Package
C1004 Sealed in Plastic Bag
C1104B Bulk Pack

PACKAGING	
Retail Boxed:	1
Plastic Bag:	1
Bulk Pack:	10

Color and style of products may vary from photos

**NEW
DESIGN**


C1109 Female Urinal

- Boxed for easy display
- Anatomically correct for female users
- Packed one (1) per retail box



C1105 10" Emesis Basin

- Made of strong polypropylene
- Rounded lines for ease in cleaning
- Sealed in plastic bag
- Sold in eaches



C1106 6 qt. Rectangular Washbasin

- Ideal for bathing bed bound patients
- Sold in eaches



C1107 Commode Specimen Pan

- Fits any standard commode
- Graduated in cc's and ounces for urine and other specimen measurements
- Sold in eaches



C1108 One Ounce Medicine Cups (30cc)

- Graduated in ounce and millimeter measurements
- Made of sturdy plastic
- Sold 100 per plastic sleeve



Our Free Standing Displays Take Up a Minimum of Space and Generate a Maximum Return on Your Investment



SW1343-Display
12 Couture™ Canes with Matching Super BigFoot™ Cane Tip



SW711
Laser Cut™, Spring Garden® and Designer Offset Canes



SW726
Rolling Mobility Center



SW730-Father
Father's Day Cane Display



SW730-Holiday
Holiday Cane Display



SW730-Summer
Spring / Summer Cane Display

Maximize Your Sales and Profits With Our High Impact Planograms

All come with complete instructions and inventory control tags.




AS30 Power of Red Planogram



AS17 Comfort Zone Cushion Display



AS60 Everyday Essentials[®] Planogram



AS10 Cushion and Pillow 4ft Planogram



AS67 Bath Safety Planogram



AS20 Cushion and Pillow 8ft Planogram

INDEX

A

ADLs.....20-25
All in One Commode.....33
Alternating Pressure
Pumps and Pads.....47
Aluminum Crutches.....13
Aluminum Walkers:
Blazer®.....19
Endurance® Folding.....14
Essential Rollator.....19
Featherlight™.....18

B

Back Cushions.....44-45
Bariatric.....49-51
Baskets for
Folding Walkers.....16
Bath Benches.....27
Bath Mats.....29
Bath Safety.....27-33
Bathing Items (ADLs).....20-21
Bed Pans.....62
Bed Sheets/Sets.....36
Bed Wedges.....37
Bibs.....59
BigFoot Cane Tips.....11
Blanket Support.....35
Blazer® Walker.....19
Blind Canes.....10
Blow Molded
Shower Benches.....27
Buckets, Commodes.....33
Butterfly Electrodes.....55
Butterfly Pillows.....46

C

Cane Accessories.....11
Cane Displays.....12
Cane Seats.....9
Canes.....4-10
Carpet Glides.....15
Cat's Meow® Canes.....5
Clothing Protector with
Crumb Catcher.....59
Coccyx Cushions.....39-42
Commodes & Accessories.....33
Commode Specimen Pan.....63
Conductive Gel.....55
Convuluted Foam
Bed Pads.....37
Couture Canes.....4
Crutches and Crutch
Accessories.....13
Cup Holder.....17
Curved Handle Canes.....9
Cushions:
Comfort.....39
Lumbar.....44-45
Memory Foam.....38-39
Neck Support.....45-46
Travel.....38, 46
Wheelchair.....40-42

D

Decubitus Bed Pads.....47
Demi Featherlight™ Walker..18
Demi Spring Garden®
Folding Canes.....6
Displays, ADL's.....23
Displays, Canes.....12
Donuts.....42-43
Dressing Items (ADLs).....20-21

E

Eating / Drinking
Utensils (ADLs).....23

Eclipse™ Pillow.....46
Elbow Protectors.....48
Electrodes.....55
Emesis Basins.....63
Essential Rollator.....19
Everyday Essentials®.....20-25
Everyday Picks®.....25
Exercise Pulley Sets.....56
E-Z Sleep Pillows.....46
E-Z Walker Skis.....15

F

Fashion Canes.....7
Featherlight™ Walkers.....18
Fiber Heel Protectors.....48
Flannel Gowns.....61
Foam Bed Pads.....37
Foam Rolls.....46
Folding Canes.....4-10
Folding Reacher.....25
Foot Stools.....50, 58
Four Wheeled Walkers.....18-19
Fracture Bed Pan.....62
Free Spirit®.....3
Furniture Protectors.....54

G

Gait Belts.....57
Gel Cushions.....41
Gel for Electrodes.....55
Gentle Touch™ Canes.....4
Glide Caps.....15
Glide Wheels.....15
Gowns - Patient.....60-61
Grab Bars.....30
Grabbers.....25

H

Hair Basin.....21
Hand Bed Rails.....34
Heavy Duty.....50-51
Heavy Duty Walker.....51
Heel/Elbow Protectors.....48
Hip Kits.....25
Home Bed Pads.....37
Home Mattress Protectors.....35
Hospital Bed Sheets.....36
Household.....22
Hygiene Wipe.....21

I

Ice Attachment - Canes.....11
Incontinence Products.....52-54
Inflatable Shampooer.....21
Invalid Rings.....42-43
IV Gowns.....61

J, K

Jackson Style Pillows.....46
Jersey Knit Sheets.....36
Junior Walkers.....14
King Size Gowns.....61
Knee Separator.....38, 44
Knee Walker.....3
Knit Sheets.....36

L

Laser Cut™ Canes.....4
Lady's Fashion Canes.....6
Leg Lifter.....22
Leg Rests.....44
Leg Walker.....3

Lotion Applicator.....21
Lumbar Cushions
Back.....44-45
Bucket Seat.....44
Deluxe, Extra Wide.....45
Wheelchair.....39
Lumbar Rolls.....45

M, N

Magnifiers.....22
Male Urinals.....62
Mattress Protectors.....35
Medicine Cups.....63
Memory Foam PF®.....38-39
Nail Brush.....21
Neck Support Cushions.....45-46

O

Offset Canes, Designer.....5
Offset Canes, Standard.....4
Opener.....22
OT Products.....55-57
Overbed Tables.....35
Overdoor Traction Sets.....56

P

Padded Toilet Seats.....32
Patient Bibs.....59
Patient Gowns.....60-61
Pedal Exerciser.....56
Pillow Cases.....36
Pillow Protectors.....35
Pillows.....38-46
Planograms.....64-65
Plastic Ware.....62-63
Platform Attachment.....16
Pontoon Bed Pan.....62
Power of Red™.....24
PT Products.....55-57
Pyramid Cane.....14

Q

Quad Canes.....8
Queen Size Gowns.....61
Quik-Sorb™ Products.....52-54

R

Raised Toilet Seats.....32
Reachers.....26
Reading Aids/
Writing (ADLs).....22
REHAB 1™ Cushions.....41
Reusable Underpads.....52-54
Rings - Inflatable.....43
Rolling Mobility Center.....12
Rubber Bath Mats.....29
Rubber Replacement Tips.....13

S

Seat Canes.....9
Shampoo Basin.....21
Sheepette® Bed Pads.....47
Sheets, Hospital Bed.....36
Sheets, Waterproof.....35
Shoehorn.....20
Shower Benches.....27
Shower Mats.....29
Shower Sprays.....28
Shower Stools.....28
Sitz Bath.....62
Sleep Pillows.....46
Sock Aids.....20
Specimen Pan.....63
Splash Guards.....33

Spring Garden® Canes.....6
Squeeze Balls/Stars.....57
Steppin' Out® Canes.....7
Swivel Seats.....58
Synthetic Sheepskin.....47

T

Tablet / Pill Crusher.....22
Tablet / Pill Splitter.....22
TENS Units.....55
THE Cushion.....39
Tips - Rubber.....13
Toilet Safety Rails.....30
Toilet Seat Cushions.....31
Toilet Seat Risers.....31-32
Traction Sets.....53
Transfer Belts.....57
Transfer Benches.....28
Transfer Boards.....17
Travel Pillows.....38, 46
Trigger Release Walkers.....14
Tub Bar.....29
Tuckable Underpads.....53

U, V

Underpads-Reusable.....52-54
Urinal - Female.....63
Urinal - Male.....62

Vinyl Handle Cane.....10
Vinyl Incontinent Pants.....54
Vinyl Mattress Protectors.....35

W, X

Walkers
Heavy Duty.....47
Folding.....14
4 Wheeled.....18-19
Walker Accessories
Baskets.....16
Folding.....14
Platform Attachment.....16
Pouches.....16-17
Sheepette® Grip Covers...16
Skis.....15
Tray.....15
Wheels.....15
Washbasins.....63
Wedged Cushions.....41
Wedges - Bed.....37
Wheelchair
Accessories.....17
Cushions.....38-39
Pouches.....15, 17
Seat Covers.....44
Trays.....17
Wheels - Walker.....15
Wide Commodes.....31
Wood Canes.....8
3XL Patient Gowns.....57

Medicare Codes









Codes are correct at the time of publication. Care should be taken, however, to double check codes for current validity



Let Our Literature Work For You!

LITERATURE ORDER FORM

Quantity Requested

<input type="checkbox"/>		Quik-Sorb™ Trifold	_____
<input type="checkbox"/>		Everyday Essentials® Brochure	_____
<input type="checkbox"/>		Power of Red™ Trifold	_____
<input type="checkbox"/>		Cushion Brochure	_____
<input type="checkbox"/>		Free Spirit® Trifold	_____
<input type="checkbox"/>		Essential Bath Safety™ Brochure	_____
<input type="checkbox"/>		Essential Medical Supply Catalog	_____
<input type="checkbox"/>		Essential Medical Supply Price List	_____

Three easy ways to request literature

First, copy/scan and print this sheet. Then:

- Give this sheet to your Essential Rep
- Fax it to: 407-770-0624
- E-mail it to: sales@essentialmedicalsupply.com

Dealer Name

Attention

Address

City, State, Zip

Phone Number



**Brands that You can Trust for
Quality and Performance!**

For Product Related Questions, Please Contact:

Toll-Free:.....(800) 826-8423

Fax:.....(407) 770-0624

Email:.....sales@essentialmedicalsupply.com

Website:.....www.essentialmedicalsupply.com

July 2019