

EFFICIENCY BY DESIGN



IQ
spiro





In the United States, a patient dies from COPD every four minutes. It's the fourth leading cause of death, with 24 million people estimated to have the disease and 12 million of those undiagnosed.

The IQspiro™ Digital Spirometry System makes detection at the primary care level simple and effective.

SIMPLE & EFFECTIVE IN-OFFICE TESTING

IQspiro™ Digital Spirometry System

The IQspiro™ from Midmark combines a digital spirometer with a Windows® based software solution that makes it fast and easy to capture, interpret and review patient test information.

- Accurate and reliable – meets or exceeds ATS 2005 Standards
- Easy grasp, lightweight design makes testing comfortable for adult or pediatric patients
- Ergonomic, T-shaped design promotes correct patient positioning during testing
- Unique eject lever alleviates handling of used mouthpieces
- Easy-latch door design improves mouthpiece stability and protection
- USB connection requires no battery to operate
- View real-time Flow/Volume Loop and Volume/Time graphs with fully customizable parameters
- Pre- and post-bronchodilator testing with automatic interpretation
- Trend patient's test measurements over time in graphical and tabular formats
- Easy to print full page reports on standard paper
- Use with today's top EHR/EMR systems or as a standalone device with IQmanager™



Save test results directly to patients' electronic records. Visit the Technical Support page on midmark.com for a complete list of our EHR/EMR partners.

Increase efficiency and improve patient care within your existing workflow

The speed, accuracy and ease of use make the IQspiro™ a cost-effective device for detecting and managing pulmonary disease – but the benefits don't end there. Talk to a Midmark representative today to see how effective and economical it can be to add spirometry testing to your patient care strategy whether you're a single practitioner or part of a larger physician group.

1 Preparing the test is as easy as retrieving the patient data from the EHR system or IQmanager™ software.



2 Patient coaching is made easier with a choice of easy-to-see incentives on full screen (image: candles).



	Actual	Predicted	Diff. Actual
FEV1	1.14	1.14	0.00%
FVC	2.15	2.15	0.00%
FEV1/FVC	53%	53%	0.00%
IRV	1.01	1.01	0.00%
TLC	3.12	3.12	0.00%
RV	2.11	2.11	0.00%
RV/TLC	68%	68%	0.00%
RVIC	0.00	0.00	0.00%
RVIC/TLC	0.00%	0.00%	0.00%

Age: 55 years Lung Age: 55 years PRED FVC: 2.15 L

Sex: Female PRED FEV1: 1.14 L

Height: 165cm Gender: Female

Race: Caucasian CURVA: 0.00

Interpretation: Normal spirometry



4 Detailed test results with auto interpretation help clinician diagnose patient's lung health accurately.



3 Avoid contact of used mouthpieces during disposal while test results are being saved directly to the patient's records.

Make your practice more efficient with an IQ Digital Diagnostic System

Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us at 1-800-MIDMARK today and we'll put together a customized suite of solutions to grow as your practice grows.

Diagnostic Products				
	EHR/EMR Connectivity	Single Platform (IQmanager™)	Standalone	Thin - Client (IQpath™)
IQecg™	■	■		■
IQholter™	■	■		■
IQholter™EX	■	■		■
IQholter™EP	■	■		■
IQstress™	■	■		
IQspiro™	■	■		■
IQvitals®	■	■	■	■
IQvitals® PC	■	■		■

Minimum System Requirements

- Windows® Vista, Windows® XP , Windows® 7, (32-bit and 64-bit)
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB port
- Mouse, keyboard and high resolution display monitor

Contact us or visit midmark.com for more detailed system requirements.

Windows is a registered trademark of Microsoft Corporation.

Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.

Models

4-000-0026 IQspiro™ USB Version with Syringe

4-000-0025 IQspiro™ USB Version

For information, pricing or to place an order, contact us at:

1-800-624-8950 or

1-310-516-6050 or

midmark.com

Midmark Diagnostics Group is an ISO 13485 Certified Company.

Midmark Diagnostics Group
1125 W. 190th Street, Gardena, CA 90248
midmark.com
1-800-624-8950

© 2011 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0915-00 Rev. A1 (8/12)



RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working together with you to redefine the future of the clinical space.

midmarkclinicalsolutions.com

Accessories



IQcal™ 3L Calibration Syringe
1-100-0007



IQmanager™ Software
4-100-1200



Disposable Nose Clips (Box of 25)
2-244-0001



Disposable Mouthpieces
2-100-1205 (Box of 25)
2-100-1206 (4 Boxes of 25)
2-100-1207 (Box of 500)



IQcart™
3-004-1000



Because we care.