

PATIENT ROOM SHARPS CONTAINERS



MDS705153
(Shown with cabinet MDS707953)

These sharps disposal containers are designed for use in patient rooms, on med carts, or in high traffic areas.

- » Counter-balanced lid creates a torturous path to facilitate safe disposal
- » Lid design facilitates one-handed operation
- » Horizontal drop design for maximum capacity to fill line
- » Easy-to-see fill levels
- » Multiple cabinet and bracket options

Item No.	Description	Accessories	Pkg
MDS705153	5 Quart, Translucent Red, Counter-Balance Lid	A C E F	20/cs
MDS705154	5 Quart, Translucent, Counter-Balance Lid	A C E F	20/cs
MDS705201F	1 Gallon, Red, Counter-Balance Lid	G H	32/cs
MDS701201F	1 Gallon, Translucent Red, Counter-Balance Lid	G H	32/cs
MDS705202F	2 Gallon, Red, Counter-Balance Lid	G H	20/cs
MDS701202F	2 Gallon, Translucent Red, Counter-Balance Lid	G H	20/cs
MDS705203	3 Gallon, Translucent Red, Counter-Balance Lid	B C I J	12/cs

PHLEBOTOMY SHARPS CONTAINERS

Portable sharps containers are designed to fit into blood drawing trays for point-of-use disposal. Square base for added stability.

- » Translucent top allows for easy identification of contents
- » Choice of single or dual openings to accommodate a variety of sharps
- » Containers lock for final disposal and safe transport
- » Brackets available to mount containers to wall and free up counter space

Item No.	Description	Accessories	Pkg
MDS705110	1 Quart, Red, Phlebotomy Container	D	100/cs
MDS705115	1.5 Quart, Red, Transportable Container		20/cs



MDS705110



MDS705115

**SEE BACK
PAGE FOR
ACCESSORIES
KEY.**