


# PROBASICS





# Table of Contents

<b>Wheelchairs</b>	
ProBasics® Transport Chairs	4
ProBasics® Value K0001 Wheelchair	6
ProBasics® Lightweight Wheelchair	7
ProBasics® High Performance Lightweight Wheelchair	8
ProBasics® K0004 Transformer™ Wheelchair	9
ProBasics® Heavy-Duty Wheelchairs	10
ProBasics® Full Reclining Wheelchair	11
<b>Walking Aids</b>	
ProBasics® Rollators	12
ProBasics® Knee Walker	15
ProBasics® Walkers	16
ProBasics® Canes and Crutches	18
<b>Bath Safety</b>	
ProBasics® Raised Toilet Seats	20
ProBasics® Commodes	21
ProBasics® Shower Chairs	22
ProBasics® Transfer Benches	23
<b>Therapeutic Support Surfaces</b>	
ProBasics® SatinAir™ APM 2 Mattress	24
ProBasics® SatinAir™ Bariatric Mattresses	25
<b>Patient Room</b>	
ProBasics® IV Stands	26

# WHEELCHAIR

## ProBasics® Transport Chairs

Featuring a fold down back for easy storage, swingaway removable footrests, durable nylon upholstery and the added safety of push-to-lock rear wheel locks and a seat positioning strap, users will feel safe and comfortable in the ProBasics Transport Wheelchair line.



### A. ProBasics® Transport Wheelchair

Weighing only 27 lb. and sporting a weight capacity of 300 lb., the Steel Transport Wheelchair is available in 17" or 19" seat widths in a stylish silver vein finish. HCPCS Code: E1038

Model #	Description
0005	17" Transport Chair
9105	19" Transport Chair

### B. ProBasics® Lightweight Aluminum Transport Wheelchair

Weighing only 19 lb. and sporting a weight capacity of 300 lb., the Lightweight Aluminum Transport Wheelchair is available in 17" or 19" seat widths and stylish burgundy, blue green or black finish. HCPCS Code: E1038

Model #	Description
9201	19" Aluminum Transport Chair - Burgundy
9201BL	19" Aluminum Transport Chair - Blue
9201BK	19" Aluminum Transport Chair - Black
9201GR	19" Aluminum Transport Chair - Green
9200BL	17" Aluminum Transport Chair - Blue



### C. ProBasics® Heavy-Duty Transport Wheelchair

- Swingaway footrests with metal footplates
- Removable padded desk arms
- Wipe-clean black vinyl upholstery
- Dual crossbraces - double gusseted
- Push-to-lock rear wheel locks
- 12" Flat free rear wheels
- HCPCS Code: E0143, E0156

Model #	Description
PB9500	Heavy-Duty Transport Chair - Burgundy
PB9500BK	Heavy-Duty Transport Chair - Black



Model #	Folded Width	Seat Dimensions	Seat-to-Floor Height	Caster Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
0005	10"	17"W x 16"D	19"	8"	300 lb.	33 lb.	30" x 23" x 10"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
9105	10"	19"W x 16"D	19"	8"	300 lb.	33 lb.	31" x 23" x 10"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
9200BL	9"	17"W x 16"D	19"	8"	300 lb.	29 lb.	31" x 23" x 9.5"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
9201 9201BL 9201BK 9201GR	9"	19"W x 16"D	19"	8"	300 lb.	29 lb.	31" x 23" x 9.5"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
PB9500 PB9500BK	N/A	22"W x 17.5"D	20"	8"	400 lb.	39.5 lb.	35.25" x 29" x 10"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Value K0001 Wheelchair

HCPCS CODE: K0001 /K0002

The ProBasics Value K0001 Wheelchair features swingaway footrests or padded, elevating legrests for user comfort, and flip-back, desk-length arms to make transfers in and out of the wheelchair easier. Composite footplates, handrims and wheels provide durability.

- Removable, padded desk-length arms
- Durable, nylon upholstery is easy to clean
- Composite footplates and handrims for durability
- Available in black finish

Model #	Description
PBECFR	Footrest
PBECLR	Legrest

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
PBEC03FR	16"W x 16"D	19"	17"	250 lb.	44 lb.	31.5" x 34.5" x 11.25"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PBEC04LR	16"W x 16"D	19"	17"	250 lb.	45.9 lb.	31.5" x 34.5" x 11.25"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PBEC09FR	18"W x 16"D	19"	17"	250 lb.	44 lb.	31.5" x 34.5" x 11.25"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PBEC10LR	18"W x 16"D	19"	17"	250 lb.	45.9 lb.	31.5" x 34.5" x 11.25"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PBEC11FR	20"W x 16"D	19"	17"	250 lb.	44.5 lb.	31.5" x 34.5" x 11.25"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PBEC12LR	20"W x 16"D	19"	17"	250 lb.	46.5 lb.	31.5" x 34.5" x 11.25"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



MODEL # PBEC09FR

# ProBasics® Lightweight Wheelchair

HCPCS CODE: K0003

The ProBasics Lightweight K0003 Wheelchair features swingaway footrest or elevating legrest models, and comes equipped with adjustable height, flip-back desk-length arms to make transfers in and out of the wheelchair easier, composite footplates and handrims for durability, and a dual axle for easy conversion to hemi height.

- Dual axle - can convert to hemi
- Aluminum multi-position front forks
- Flip-back, adjustable height arms
- Heavy-duty composite side panels
- Composite footplates
- Composite handrims
- Silver vein finish

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
2010AH	16"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2011AH	16"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
2012AH	18"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2013AH	18"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
2014AH	20"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2015AH	20"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



MODEL # 2012AH

# ProBasics® High Performance Lightweight Wheelchair

HCPCS CODE: K0004

The ProBasics High Performance Lightweight K0004 Wheelchair features swingaway footrest or elevating legrest models and comes equipped with flip-back, desk-length, space saver arms to make transfers in and out of the wheelchair easier, anti-tip devices for added safety, composite footplates and handrims for durability, and a dual axle for easy conversion to hemi height.

- Flip-back, padded, desk-length space saver arms
- Composite footplates
- Composite handrims
- Heel loops included with footrest
- Durable nylon upholstery
- Black finish

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
PB1610	16"W x 16"D	18"-20"	16"	300 lb.	39.2 lb.	32" x 31" x 12"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1611	16"W x 16"D	18"-20"	16"	300 lb.	42.5 lb.	32" x 31" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1810	18"W x 16"D	18"-20"	16"	300 lb.	39.2 lb.	32" x 31" x 12"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1811	18"W x 16"D	18"-20"	16"	300 lb.	42.5 lb.	32" x 31" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1812	20"W x 16"D	18"-20"	16"	300 lb.	41.8 lb.	32" x 32" x 12"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1813	20"W x 16"D	18"-20"	16"	300 lb.	43.1 lb.	32" x 32" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



MODEL # PB1811



# ProBasics® K0004 Transformer™ Wheelchair

HCPCS CODE: K0004

The ProBasics K0004 Transformer Wheelchair is a dual-purpose, high-strength, lightweight wheelchair that can quickly change from a self-propelled wheelchair to a compact transport wheelchair. It comes complete with two sets of lightweight wheel locks for both modes and folds for easy storage.

- Flip-back, desk-length arms with padded armrests
- Lightweight, black, woven nylon upholstery with accessory pouch sewn on back
- Footrests include heel loops and composite footplates
- Fold-down back with easy-release paddles
- Detachable, flip-up rear anti-tippers
- Quick-release rear wheels for easy conversion to transport chair
- Composite handrims for durability
- Stylish blue marble finish

Model # Description  
PB191ELR Legrest Conversion Kit

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
PB1916	16"W x 16"D	19"	16"	250 lb.	31.9 lb.	31" x 32" x 14"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1918	18"W x 16"D	19"	16"	250 lb.	33 lb.	31" x 32" x 14"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1920	20"W x 16"D	19"	16"	300 lb.	36.3 lb.	31" x 32" x 14"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



Transport Chair Configuration



MODEL # PB1916

# ProBasics® Heavy-Duty Wheelchairs

HCPCS CODE: K0007

## A. ProBasics® Heavy-Duty Wheelchair

The ProBasics Heavy-Duty K0007 Wheelchair features swingaway footrest or elevating legrest models.

- Dual axle converts to hemi
- Available with hammertone finish
- Heavy-duty reinforced frame has gusseted side frames, caster journals, dual crossbraces and 5/8" rear mounted axle
- Mid-position slide post adds strength and allows easy folding
- 24" Composite mag wheel with urethane tires
- Heavy-duty reinforced upholstery with extra heavy-duty liner



MODEL # PB1013



MODEL # PB1326

## B. ProBasics® Extra Wide Wheelchair

The reinforced frame with double crossbrace offers a rugged, heavy-duty wheelchair for bariatric users. Two (26" and 28") sizes of seat width are provided to suit the individual needs of users.

- Armrests and front riggings are Invacare compatible
- Detachable desk arms
- Flat-free tires
- Double crossbraces with link bars
- 24" x 1-3/8" spoke wheels with high strength M15 axles
- Swingaway footrests with large aluminum footplates optional
- Elevating legrests optional
- Available in black finish

Model #	Description
PB132ELR	Elevating Legrest
PB132FR	Swingaway Footrest

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
PB1013	22"W x 18"D	19" - 21"	16"	450 lb.	70.4 lb.	37.5" x 33.5" x 12"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1014	22"W x 18"D	19" - 21"	16"	450 lb.	72 lb.	37.5" x 33.5" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1124	24"W x 18"D	19" - 21"	16"	450 lb.	72.6 lb.	37.5" x 33.5" x 12"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
PB1224	24"W x 18"D	19" - 21"	16"	450 lb.	74 lb.	37.5" x 33.5" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1326	26"W x 20"D	21"	17"	700 lb.	57.5 lb.	38" x 35" x 13"	Sold Separately	5 Years Frame/Crossbraces 1 Year Components
PB1328	28"W x 20"D	21"	17"	700 lb.	59.4 lb.	38" x 35" x 13"	Sold Separately	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Full Reclining Wheelchair

HCPCS CODE: K0006

The ProBasics K0006 Reclining Wheelchair features removable desk-length arms and padded elevating legrests. An easy-to-operate reclining mechanism allows for the chair's back to recline smoothly, and rear anti-tippers and a headrest extension add safety and comfort.

- Composite mag wheels with dual seal precision bearings
- Rear anti-tippers & headrest extension
- Padded removeable desk-length armrests
- Padded elevating legrests come standard
- Finger trigger release mechanism allows for smooth reclining
- Enhanced reclining mechanism
- Warranty: 5 years on weld and crossbrace, 1 year on upholstery and wheel locks
- Available with silver vein finish



Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
PB1816	16"W x 16"D	19.5"	31"	300 lb.	62 lb.	40" x 36" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1818	18"W x 16"D	19.5"	31"	300 lb.	63 lb.	40" x 36" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1820	20"W x 16"D	19.5"	31"	300 lb.	63 lb.	40" x 36" x 12"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
PB1822	22"W x 16"D	19.5"	31"	400 lb.	64 lb.	39.5" x 36" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

**MODEL # PB1816**

## Wheelchair Accessories

Model #	Description
PB1811RAT	Anti-Tippers for K0003 and K0004 Wheelchairs - Not Shown
PB0535	IV Pole
PB016	Oxygen Holder
PB017	IV Holder - Not Shown
PB3065	Seatbelts



PB0535



PB016



PB3065

# ProBasics® Rollators



## A. ProBasics® Junior Rollator

The ProBasics Junior Aluminum Rollator features an easy-squeeze loop brake system, adjustable handle height, comfortable, padded backrest and seat and lightweight aluminum frame. Comes standard with a storage pouch under seat for added privacy. HCPCS Code: E0143, E0156

Model #	Description
PB1028	Junior Rollator - Flame Burgundy
PB1028BL	Junior Rollator - Flame Blue



## B. ProBasics® Steel Rollator

The ProBasics Steel Rollator provides comfort and durability. Featuring adjustable-height handles, molded hand grips, removable storage pouch and curved padded backrest, users can customize their comfort. HCPCS Code: E0143, E0156

Model #	Description
PB1026	Steel Rollator - Burgundy
PB1026BL	Steel Rollator - Blue
PB1026BK	Steel Rollator - Black



## C. ProBasics® Value Rollator

The ProBasics Value Rollator provides adjustability and comfort with a durable, lightweight aluminum frame. Featuring a padded seat and a removable storage pouch and curved padded backrest, users will feel comfortable using the Value Rollator. HCPCS Code: E0143, E0156

Model #	Description
PB1024-BK	Value Rollator - Black
PB1024-BL	Value Rollator - Blue
PB1024-RD	Value Rollator - Red



### D. ProBasics® Low Profile Aluminum Rollator

The ProBasics Low Profile Aluminum Rollator features an easy-squeeze loop brake system and a two-button, removable back rest all on a lightweight aluminum frame. The low 19" seat-to-floor height provides a comfortable seat height for petite users. HCPCS Code: E0143, E0156

Model #	Description
PB1025LP	Aluminum Rollator - Metallic Burgundy
PB1025LPBK	Aluminum Rollator - Black
PB1025LPBL	Aluminum Rollator - Blue



### E. ProBasics® Flame Finish Aluminum Rollator

Featuring a lightweight, aluminum welded frame, full swivel front casters and easy-squeeze brakes for progressive braking control, this rollator is easy to maneuver. The flip-up, padded seat reveals a quick-access storage pouch with an articulating backrest for convenient storage and transportation. HCPCS Code: E0143, E0156

Model #	Description
1037	Aluminum Rollator - Flame Red
1037BL	Aluminum Rollator - Flame Blue
1037BK	Aluminum Rollator - Flame Silver
1037F	Aluminum Rollator - Flame Fuchsia



### F. ProBasics® Deluxe Aluminum Rollator

Ideal for those who spend more time outdoors, the Deluxe Aluminum Rollator features a durable, lightweight aluminum frame, large 8" wheels and a rugged in-line braking system. This rollator offers a seamless, padded polyurethane seat for improved comfort and ease of care. The textured seat is waterproof, anti-bacterial and flips up to reveal a large, zippered compartment for secure storage. HCPCS Code: E0143, E0156

Model #	Description
1032MAR	Deluxe Aluminum Rollator - Flame Red
1032BLMAR	Deluxe Aluminum Rollator - Flame Blue

**G. ProBasics® Bariatric Aluminum Rollator**

The ProBasics Bariatric Aluminum Rollator features a durable, welded frame, locking loop brakes, large 8" wheels and a 400 lb. weight capacity to provide users with a heavy-duty, reliable product. With a low seat-to-floor height and five ribs on the heavy-duty front and rear wheels, users can feel more stable. HCPCS Code: E0149



Model #	Description
1033MAR	Bariatric Aluminum Rollator - Flame Red
1033MARBL	Bariatric Aluminum Rollator - Flame Blue

Model #	Dimensions Inside Hand Grips	Seat Dimensions	Seat-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Weight Capacity	Limited Warranty
PB1028 PB1028BL	14"	13.5"W 13"D	19"	29"-34"	6"	250 lb.	3 Years Frame 1 Year Seat/Brakes
PB1026 PB1026BL PB1026BK	18"	14"W 13"D	20"	30.5"-36.75"	6"	300 lb.	3 Years Frame 1 Year Seat/Brakes
PB1024-BK PB1024-BL PB1024-RD	18"	14"W 12"D	21"	31"-37"	6"	300 lb.	3 Years Frame 1 Year Seat/Brakes
PB1025LP PB1025LPBK PB1025LPBL	18"	14"W 13"D	19"	31.5"-35.5"	6"	300 lb.	3 Years Frame 1 Year Seat/Brakes
1037 1037BL 1037BK 1037F	18"	14"W 12"D	21"	31"-37"	6"	300 lb.	3 Years Frame 1 Year Seat/Brakes
1032MAR 1032BLMAR	18"	14"W 12"D	21"	31.5"-35.5"	8"	300 lb.	3 Years Frame 1 Year Seat/Brakes
1033MAR 1033MARBL	23.5"	19"W 13"D	22"	33.5"-38.5"	8"	400 lb.	3 Years Frame 1 Year Seat/Brakes

NOTE: Rollators are not intended to be used as wheelchairs

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Knee Walker

HCPCS CODE: E0118

The perfect aid for those recovering from ankle or foot conditions such as surgery, breaks, sprains and amputations, the ProBasics Knee Walker is comfortable to use and easy to maneuver. The handle and seat height can be easily adjusted without tools to meet a variety of user needs. The Knee Walker is simple to fold for storage and comes with a removable wire storage basket.

Model #	Overall Dimensions	Kneepad Dimensions	Kneepad-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Weight Capacity	Units Per Carton	Limited Warranty
1030KW	16"W 31"L	6"W 14"L	17"-23"	32"-41.75"	7.5"	300 lb.	1	3 Years Frame 1 Year Seat and Brakes

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.





## A. ProBasics® Value Walker - Not Shown

- Easy-to-use, two-button release allows compact storage
- Each side operates independently
- Weight Capacity: 300 lb.
- Warranty: Limited lifetime against manufacturer's defects on the frame
- HCPCS Code: E0143

Model #	Description
PB1064	Value Walker with 5" Wheels
PB1064Y	Value Walker with 5" Wheels - Youth

## B. ProBasics® Dual Release Walker

- Easy-to-use, red button release allows simple folding by fingers, palm, or side of hand
- Lightweight 1" anodized aluminum frame construction provides superior strength
- Dual release folding mechanism allows walker sides to fold independently
- Weighs under 7 lb.
- HCPCS Code: E0143 with wheels, E0135 without wheels

Model #	Description
PB1050	Dual Release Walker
PB1050Y	Dual Release Walker - Youth
PB1060MRBL	Dual Release Walker with Wheels - Flame Blue



## C. ProBasics® Deluxe Two Button Folding Walker with Wheels Installed

- Easy to use two button release allow compact storage
- Reinforced with saddle formed riveted side supports
- Weight Capacity: 300 lb.
- Warranty: Limited lifetime against manufacturer's defects on the frame
- HCPCS Code: E0143

Model #	Description
PB1062	2 Button Walker with Wheels Installed - Adult
PB1062Y	2 Button Walker with Wheels Installed - Youth



**D. ProBasics® Sure Lever Release™ Folding Walker**

- Lightweight 1" anodized aluminum construction
- Dual lever release offers users easy folding while holding onto walker
- Compatible with all ProBasics walker accessories
- HCPCS Code: E0143 with wheels, E0135 without wheels

Model #	Description
PB1080	Sure Lever Release Walker
PB1080Y	Sure Lever Release Walker - Youth
PB1085	Sure Lever Release Walker with 5" Wheels - Hammetone
PB1085BL	Sure Lever Release Walker with 5" Wheels - Flame Blue



**E. ProBasics® Bariatric Walker**

- Easy-to-use, push button mechanisms may be operated by fingers, palm, or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability when standing up
- Front crossbraces are located near top of walker to allow full stride
- Sturdy, 1" diameter aluminum tubing provides stability while remaining lightweight
- Folds to a slim 4.5" for easy transportation and storage
- Height adjusts from 32" to 39" in 1" increments
- HCPCS Code: E0143 with wheels, E0135 without wheels

Model #	Description
PB1090	Bariatric Walker
PB1091B	Bariatric Walker with 5" Wheels



Model #	Overall Width	Width Inside Hand Grips	Height Adjustment	Weight Capacity	Limited Warranty
PB1064 PB1064Y	23"	17"	32"-39" 25"-32"	300 lb.	Limited Lifetime
PB1050 PB1050Y PB1060MRBL	23"	17"	32"-39" 25"-32"	300 lb.	Limited Lifetime
PB1062 PB1062Y	23"	17"	32"-39" 25"-32"	300 lb.	Limited Lifetime
PB1080 PB1080Y	22"	18"	32"-39" 25"-32"	300 lb.	Limited Lifetime
PB1085 PB1085BL	22"	18"	32"-39" 25"-32"	300 lb.	Limited Lifetime
PB1090 PB1091B	23"	19"	32"-39"	500 lb.	Limited Lifetime

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Canes and Crutches

## A. ProBasics® Value Curved Handle Cane

- Anodized aluminum finish with hypolon grip
- Curved handle allows cane to be retained on arm
- Adjustable height with safety locking silencer
- Adjusts from 31" to 40" in 1" increments
- HCPCS Code: E0100

Model #	Description
PB1057	Value Curved Handle Cane

## B. ProBasics® Offset Canes

- Anodized aluminum finish
- Soft cushioned hand grip
- Locking height adjustment silencer
- Height adjustment from 31" to 40" in 1" increments
- Available in bronze, black, or silver
- HCPCS Code: E0100

### Bariatric Offset Canes - Not Shown

- Solid brass, push-button adjustment for extra stability
- Secured with a locking ring for added safety
- Comfortable, black, PVC handgrip
- Handle height adjusts in 1" increments from 29" to 38"
- 500 lb. weight capacity
- HCPCS Code: E0105

Model #	Description
PB2057	Offset Cane with Strap - Silver
PB2057K	Offset Cane with Strap - Black
PB2057B	Offset Cane with Strap - Bronze
PB2058	Bariatric Offset Cane

## C. ProBasics® Quad Canes

- Lightweight aluminum upper shaft
- Welded steel base, chrome finish
- Push-button height adjustment from 31" to 40" in 1" increments
- Offset handle provides balance and stability
- Handle rotates for left or right hand use
- Locking ring system
- HCPCS Code: E0105

### Bariatric Quad Canes - Not Shown

- All steel construction
- Double-plated base with extremely strong welds
- Solid brass, push button adjustment for extra stability
- Handle height adjusts in 1" increments from 29" to 38"
- Secured with a locking ring for added safety
- Comfortable, black, PVC handgrip
- Dual tubing insert for added strength and safety
- 500 lb. weight capacity
- HCPCS Code: E0105

Model #	Description
PB4057	Quad Cane - Large Base
PB4058	Bariatric Quad Cane - Large Base
PB3057	Quad Cane - Small Base
PB3058	Bariatric Quad Cane - Small Base



### D. ProBasics® Folding Canes

- Black or bronze anodized aluminum finish
- "T" Shaped handle
- Folds easily for travel or storage
- Adjustable height 32.5" to 36.5"
- Convenient carry and display pouch
- HCPCS Code: E0100

Model #	Description
PB2055B	Folding Cane - Bronze
PB2055K	Folding Cane - Black



### E. ProBasics® Folding Walk Cane

- Anodized aluminum frame
- Lightweight compact design folds easily
- Special angled footpieces provide stability
- Adjustable height from 28" to 33"
- HCPCS Code: E0135


Model #	Description
PB6055	Folding Walk Cane



### F. ProBasics® Aluminum Crutches

- Solid one piece lightweight anodized aluminum frame
- High impact plastic armrests and hand grips
- Fully assembled with cushions, handgrips, and rubber tips
- Push button adjustment allows for safe, accurate and quick adjustment of height
- Height adjustment scale on center tube
- HCPCS Code: E0114

Model #	Description
PB1055A	Aluminum Crutches - Adult
PB1055T	Aluminum Crutches - Tall
PB1055Y	Aluminum Crutches - Youth



Model #	Height Adjustment	Weight Capacity	Quantity Per Carton	Limited Warranty
PB1057	31"-40"	250 lb.	10	2 Years Frame
PB2057 PB2057K PB2057B	31"-40"	250 lb.	10	2 Years Frame
PB2058	29"-38"	500 lb.	1	2 Years Frame
PB4057 PB3057	31"-40"	250 lb.	6	2 Years Frame
PB4058 PB3058	29"-38"	500 lb.	2	2 Years Frame
PB6055	28"-33"	250 lb.	4	2 Years Frame
PB2055B PB2055K	32.5"-36.5"	250 lb.	1	2 Years Frame
PB1055A	45"-53"	300 lb.	8	2 Years Frame
PB1055T	53"-61.5"	300 lb.	8	2 Years Frame
PB1055Y	37"-46"	300 lb.	8	2 Years Frame

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

## ProBasics® Raised Toilet Seats



### A. ProBasics® Clamp-On Raised Toilet Seat

- Front clamping mechanism ensures secure, easy locking onto toilet
- Prevents shifting of seat while transferring
- Comfortable support arms on PB408 add security when using toilet seat
- Seat elevation 5.5"
- HCPCS Code: E0243

Model #	Description
PB308	Clamp-On Raised Toilet Seat
PB408	Clamp-On Raised Toilet Seat with Arms



### B. ProBasics® Raised Toilet Seat

- Recommended for use with Toilet Safety Frame
- Attaches without tools and fits on most round bowl configurations\*
- Rear seat wings for more secure attachment
- Contoured seat for internal rotation of hips
- HCPCS Code: E0243

Model #	Description
PB307	Raised Toilet Seat



### C. ProBasics® Toilet Safety Frame

- Arms fold back for cleaning and are easily removed for transfers
- Armrests provide secure hand support
- Frame is made of anodized, non-corrosive 1" aluminum tubing
- HCPCS Code: E0243

Model #	Description
PB416	Toilet Safety Frame

Model #	Seat Dimensions	Quantity Per Carton	Weight Capacity	Limited Warranty
PB308	17" x 17"	1	250 lb.	2 Years Frame
PB408	17" x 17"	1	250 lb.	2 Years Frame
PB307	17" x 17"	1	250 lb.	2 Years Frame
PB416	N/A	2	250 lb.	2 Years Frame

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

\* ALWAYS make sure that the Raised Toilet Seat is securely in place and stable before each use. Although the Raised Toilet Seat is designed for use on most popular toilets, it may not work on certain unusual, elongated or artistic bowls.

# ProBasics® Commodes

## A. ProBasics® Deluxe All-in-One Commode

- Elongated seat
- 1" gray powder-coated steel
- Snap-on seat and lid
- Three versatile functions: bedside, toilet safety frame or raised toilet seat
- Includes bucket, lid and splash guard
- Adjustable height
- Fully assembled
- HCPCS Code: E0163

Model #	Description
PB412GPC	Deluxe All-in-One Commode



## B. ProBasics® Folding Commodes

- Fast, easy assembly and disassembly
- Snap-on, elongated seat and lid
- Includes bucket, lid and splash guard
- HCPCS Code: E0163

Model #	Description
PB420	Value Folding Commode
PB419	Deluxe Folding Commode with Elongated Seat



## C. ProBasics® Drop Arm Commode

- 1" Gray powder-coated steel
- EZ Release drop arm allows for safe side transfer
- Bucket with lid and splash guard
- Removable back rest
- Soft hypolon armrests
- Snap-on seat with lid
- Tool-free assembly
- HCPCS Code: E0165

Model #	Description
PB413	Drop Arm Commode



## D. ProBasics® Bariatric Commodes

- 450 lb. weight capacity
- Designed for bedside use
- Rectangular pan and pan holder
- Armrests for added comfort
- HCPCS Code: E0168

Model #	Description
PB411	Bariatric Commode
PB310	Bariatric Commode - Extra Wide Seat

## E. ProBasics® Bariatric Drop Arm Commode

- 1" Heavy duty steel construction
- Large comfortable snap-on seat
- 650 lb. weight capacity
- Drop arm for easy transfer
- HCPCS Code: E0168

Model #	Description
PB413BAR	Bariatric Drop Arm Commode

Model #	Width Between Arms	Seat Dimensions		Seat-to-Floor Height	Quantity Per Carton	Weight Capacity	Limited Warranty
PB412GPC	19"	14" x 16"	16"-22"	4	350 lb.	2 Years Frame	
PB420 PB419	18"	14" x 16.75"	23"	4	300 lb.	2 Years Frame	
PB413	20"	14" x 16"	17"-22"	1	300 lb.	2 Years Frame	
PB411	21"	15.5" x 15.5"	17"-23"	2	450 lb.	2 Years Frame	
PB310	24"	24" x 18.5"	17"-23"	2	650 lb.	2 Years Frame	
PB413BAR	24.5"	24" x 18"	19.5"-23.5"	1	650 lb.	2 Years Frame	

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Shower Chairs

## A. ProBasics® Bath Chairs

- 1" Anodized aluminum frame
- Legs angle outward to resist tipping
- Blow molded plastic seat and back with drainage holes & built-in handles
- Available with or without back
- Tool-free assembly
- HCPCS Code: E0240

### Bariatric Bath Chair - Not Shown

- Comfortable, blow-molded seat with integrated handles for additional support
- Dual cross brace for extra stability with extremely durable welds
- High-strength, aluminum frame to prevent rusting
- HCPCS Code: E0240

Model #	Description
PB102	Bath Chair with Back
PB202	Bath Bench without Back
PB401	Bariatric Bath Chair with Back



## B. ProBasics® Bath Chair with Padded Arms

The ProBasics Bath Chair with Padded Arms features adjustable height legs to accommodate multiple user heights. The durable aluminum frame secures to surfaces with large, suction cup feet for added security while drainage holes in the bench aid in reducing slipping. HCPCS Code: E0240

Model #	Description
PB302	Bath Chair with Arms



Model #	Base Dimensions	Seat-to-Floor Height	Quantity Per Carton	Weight Capacity	Limited Warranty
PB202	17" x 17"	14" - 19"	4	300 lb.	2 Years Frame
PB102	17" x 17"	14" - 19"	4	300 lb.	2 Years Frame
PB302	17" x 14.5"	15" - 20"	4	300 lb.	2 Years Frame
PB401	21" x 17.5"	16.5" - 20.5"	2	550 lb.	2 Years Frame

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® Transfer Benches



## A. ProBasics® Transfer Bench

- Tool-free assembly
- Textured seat with drain holes
- Hand hole in back for easy transport
- Lightweight aluminum frame
- Large suction cup feet
- HCPCS Code: E0247

Model #	Description
PB417	Transfer Bench



## B. ProBasics® Padded Transfer Bench

- Tool-free assembly
- Soft cushioned seat and back for comfort
- Lightweight aluminum frame
- Large suction cup feet
- HCPCS Code: E0247

Model #	Description
PB415	Padded Transfer Bench



## C. ProBasics® Bariatric Transfer Bench

- 1" Aluminum tube construction
- Added crossbraces for stability
- Tool-free assembly
- 500 lb. weight capacity
- HCPCS Code: E0248

Model #	Description
PB418	Bariatric Transfer Bench

Model #	Seat Dimensions	Seat-to-Floor Height	Quantity Per Carton	Weight Capacity	Limited Warranty
PB417	27" x 16"	18"-23"	2	300 lb.	2 Years Frame
PB415	24" x 16"	21"-25"	2	300 lb.	2 Years Frame
PB418	26" x 18"	18"-24"	2	500 lb.	2 Years Frame

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# ProBasics® SatinAir™ APM 2 Mattress

HCPCS CODE: E0277

The ProBasics SatinAir APM 2 Mattress is ideal for those seeking a budget friendly, quality product. Combining low air loss and alternating pressure therapies in one, the durable and user-friendly SatinAir APM 2 Mattress is designed for the prevention and treatment of stage IV pressure ulcers. This quiet and powerful mattress system features 20 8" bladders, allowing for deep immersion and excellent pressure redistribution. The laser air holes improve air circulation and minimize heat build-up to keep the user cool and dry.

## Mattress and Cover Features

- Removable, water-resistant, non-shear, antimicrobial nylon cover is quilted for extra comfort
- 20, 8" deep air bladders alternate at 10-minute intervals
- 16 laser air holes keep the user cool and dry
- CPR release allows for rapid deflation

## Pump Features

- Powerful yet quiet 7.5 LPM pump
- Low air loss and alternating pressure therapies in one system
- Air flow dial allows customization for individual user needs
- 10-minute alternating pressure intervals
- Low pressure alarm warns the user and caregiver during power outage or low pressure



Model #	Overall Dimensions	Weight Capacity	Compressor	Alternating Time Cycle	Alarms	Regulation	Limited Warranty
SATINAPM2QC	80"L 35"W 8"H	350lb.	7.5 LPM	10-Minute Intervals	Audio and Visual	UL Approved	18 Months

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# THERAPEUTIC



# ProBasics® SatinAir™ Bariatric Mattresses

HCPCS CODE: E0277

The SatinAir Bariatric Mattress Series is the first choice of clinicians to treat and prevent bariatric users from pressure ulcers.

- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution
- Low air loss function controls and manages moisture to prevent skin maceration
- Cell-on-cell mattress design prevents user from bottoming out in case of power failure
- 10" air cells provide comfort and support
- Audible low pressure alarms alert clinicians when preset pressure levels decrease
- Static mode allows for easy transfer in and out of bed
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode
- CPR pull tab quickly releases air from mattress



Model #	Overall Dimensions	Weight Capacity	Compressor	Alternating Time Cycle	Alarms	Regulation	Limited Warranty
PBAPMBARI42	80"L 42"W 10"H	650lb.	15 LPM	10, 15, 20, 25 Minute Intervals	Audio	UL Approved	2 Years
PBAPMBARI48	80"L 48"W 10"H	650lb.	15 LPM	10, 15, 20, 25 Minute Intervals	Audio	UL Approved	2 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# C SUPPORT


## ProBasics® IV Stands

HCPCS CODE: E0776

### A. ProBasics® 5 Leg IV Stands

- Chrome plated steel pole
- Aluminum base provides more stability
- Twin wheel swivel casters
- Twist or knob lock provides a full range of height positions
- Detachable head for easy on/off
- Height adjustments 45.5" to 79"

Model #	Description
PB2135	5 Leg 4 Hook - Not Shown
PB4035	5 Leg 2 Hook



### B. ProBasics® 4 Leg IV Stand

- Chrome plated and hammertone steel pole
- Offset low profile base
- Twin wheel swivel casters
- Twist lock provides full range of height positions
- Two pronged removable top for easy storage and transport
- Height adjustment 45.5" to 79"
- Base size 23"

Model #	Description
PB5035	4 Leg 2 Hook





**Invacare Corporation**

www.invacare.com

**Homecare**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Continuing Care**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 668-2337

**Canada**

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2015

Invacare Corporation.  
All rights reserved.  
Trademarks are  
identified by the symbols  
™ and ®. All trademarks  
are owned by or licensed  
to Invacare Corporation  
unless otherwise  
noted. ProBasics is a  
registered trademark  
and Transformer, Sure  
Lever Release and  
SatinAir are trademarks  
of Garden City Medical  
Inc. Specifications  
are subject to change  
without notification.  
Form No. 14-387  
Rev. 02/15  
150201



**Yes, you can.®**