


Posey® Protective Caps

The Newest Design for Patient Fall Protection

6040, 6041, 6042


Falls with serious injury are estimated to cost hospitals \$27,000 per incident¹ and fallers typically stay in the hospital 6.3 days longer than non-fallers.²

6040



Leading Causes of TBI

One of the more serious fall related injuries is a fall leading to a Traumatic Brain Injury (TBI). In 2010, about 2.2 million emergency department (ED) visits were associated with TBI and it was a diagnosis in more than 280,000 hospitalizations. TBIs contribute to the deaths of more than 50,000 people each year.



6041



Effective fall protection for your patients

- Available in three unique styles to better accommodate a wide variety of patients/residents.
- Unique design of Caps with built in protective pads offers a more discreet and low profile look than many commercially available protective head protectors, helping to preserve dignity for patient/residents.
- Inner lining of Caps offers comfort and reduced heat buildup, helping to improve patient/resident compliance and adoption.
- Lightweight design is available in 3 sizes to fit a wide population of patients and residents.
- All three caps are machine washable.

6042



NOTE: Due to the random possibility of falls, the Posey Company makes no guarantee, express or implied, that the user is protected from Traumatic Brain Injury.

¹ S et al. A cost-effectiveness analysis of proposed national falls prevention program. *Can Geriatr Med.* 2010;26(4):751-66.



² Catherine A. Wong, M.P.H., M.T.(A.S.C.P.); Angela J. Recktenwald, M.P.H.; Marilyn L. Jones, B.S.N., M.P.H.; Brian M. Waterman, M.P.H.; Mara L. Bollini, B.S.N., M.H.A.; Wm. Claiborne Dunagan, M.D. The Cost of Serious Fall-Related Injuries at Three Midwestern Hospitals, *The Joint Commission Journal on Quality and Patient Safety*, February 2011 Volume 37 Number 2.




Posey®

The Newest Design for Patient Fall Protection

The new Posey Protective Caps help protect patients from a serious head injury by providing a layer of protection around the head. The Posey Protective Caps are designed with integrated D3O foam construction, which is tested to provide high level impact absorption and protection of the head in the event of a fall.



Impact Reduction, Posey Baseline



SIZING CHART

Size	Circumference
S	21.0" (53.3cm)
M	22.5" (57.2cm)
L	24.0" (61.0cm)

All styles are available in Small, Medium, and Large sizes



Free 30-day trial. Call 1.800.447.6739

*The warranty is for the life of the product against all defects in materials and workmanship.

To request a Posey catalog or speak with an associate, call 1.800.44.POSEY (76739) or visit us online at Posey.com

