



Sprinter XL4

The Sprinter XL4 is a heavy duty four-wheel scooter with a 350 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The 12in. solid wheels give the unit a 3in. ground clearance for operation on multiple terrains.

- + 350lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



Burgundy

Model Number: 889B-4
Medicare Billing Code: K0899

HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.



Enduro XL3⁺

The Enduro XL3⁺ is a heavy duty three-wheel scooter with a 500 lb. weight capacity. It has a delta style tiller which houses the controls for the lighting system, turn signals, wig wag and horn. It has a dual purpose wig wag which allows the user to push forward or pull back on the drive lever depending on dexterity for ease of driving. The unit comes standard with a 20in. seat that can accommodate most riders, and has swivel and recline features for added comfort and ease of use.

- + 500lb Weight Capacity
- + Slide & Swivel Seating
- + Full Suspension System
- + Full Lighting Package
- + Adjustable Tiller



Burgundy

Model Number: 778XLSBN
Medicare Billing Code: K0899

HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.



Technical Specifications

Overall Dimensions (L x W x H)	in	52" x 23" x 30-42"
Number, Size of Tires	Front, Rear	2, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 350
Weight of Heaviest Part ²	lb	108
Weight ³	lb	181
Each Battery Module Weight	lb	32
Battery Size	Type x PCS	12V33-40AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(^)	14(8)
Range (Per Charge) ⁵	Mile	25
Turning Radius	in	54"
Ground Clearance	in	3"
Handle Bar	Type	T-Bar
Model #		889B-4
Standard Seat	W x D x BH	SV191821B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.

Technical Specifications

Overall Dimensions (L x W x H)	in	53" x 28" x 52"
Number, Size of Tires	Front, Rear	1, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 500
Weight of Heaviest Part ²	lb	104
Weight ³	lb	222
Each Battery Module Weight	lb	50
Battery Size	Type x PCS	12V50-75AH x 2
Charger	Amp	Off board (8A)
Output Power of Motor	HP	2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(^)	14(8)
Range (Per Charge) ⁵	Mile	25-39
Turning Radius	in	49"
Ground Clearance	in	4"
Handle Bar	Type	Delta
Model #		778XLSBN
Standard Seat	W x D x BH	SV201820B