Clean!

## **CPAP Concentrated Cleanser**

## A Fresh Approach for Cleaning Tubing and Accessories

Citrus II<sup>®</sup> CPAP Concentrated Cleanser contains all natural solvents effective in cleaning all types of CPAP machine tubes and accessories.

Use weekly to ensure a clean and healthy night's sleep. CPAP Concentrated Cleanser is made from vegetable derived cleaners and leaves no residue! Simply soak tubes and accessories, with the exclusion of the mask, in a mixture of 4 oz. CPAP Concentrated Cleaner to 1 gallon of water. Then, simply rinse and let air dry. Also available, Citrus II CPAP Mask Wipes and Citrus II CPAP Mask Cleaner in two

convenient ready-to-use

spray bottles.



Citrus II CPAP Mask Wipes Contains 62 Wipes 5" x 8" (12.7 cm x 20.3 cm) # 635871639



# 635871164 8 oz. Spray # 635871165



20-Piece Display Contains: 9 - CPAP Mask Wipes, 7-1.5 oz. Spray Cans, 4-CPAP Concetrated Cleansers (16.5"H x 17.5"W x 10.75"D) # 236371333



Citrus II CPAP Concentrated Cleanser 32 oz. Bottle Makes 8 Gallons # 635871849

SS #63026b REV 0508