

SAFETY DATA SHEET

SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

PRODUCT NAMECitrus II CPAP Mask Wipes

PRIMARY APPLICATION CPAP mask wipe. See product label for proper usage instructions.

COMPANY IDENTITY

Beaumont Products, Inc.

COMPANY ADDRESS

1560 Big Shanty Drive

COMPANY CITY, STATE, ZIP

Kennesaw, GA 30144

COMPANY PHONE

1-800-451-7096

EMERGENCY PHONE CHEMTEL: 1-800-255-3924, Int'l: 1-813-248-0585

SECTION 2: HAZARDS IDENTIFICATION

Warning

HAZARD CLASSIFICATION Acute toxicity, oral (chapter 3.1). Hazard Category 5.

Skin corrosion, irritation (chapter 3.2). Hazard Category 3.

Serious eye damage/eye irritation (chapter 3.3). Hazard Category 2B.

HAZARD STATEMENTS H303: May be harmful if swallowed.

H316: Causes mild skin irritation.

H320: Causes eye irritation.

PRECAUTIONARY STATEMENTSP101: If medical advice is needed, have product container or label at hand.

P102: Keep out of reach of children.

P103: Read label before use.

P264: Wash hands thoroughly after handling.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

CHEMICAL NAME CAS# WT%

There are no OSHA Hazardous Substances in this product. The exact percentages of ingredients in this product have been withheld as a trade secret. Any hazards associated with the finished product are listed in Section 2 of this SDS.

SECTION 4: FIRST AID MEASURES

RELEVANT ROUTES OF EXPOSURE Skin Contact, Eye Contact, Ingestion

SKIN CONTACT FIRST AID If skin irritation occurs: Get medical advice/attention (P332+P313).

EYE CONTACT FIRST AID IF IN EYES: Rinse cautiously with water for several minutes. Remove contact

lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: Get medical advice/attention (P337+P313).

INGESTION FIRST AID Call a POISON CENTER/doctor/physician if you feel unwell (P312).

Revision Number: 3.0

Citrus II CPAP Wipes Safety Data Sheet

SECTION 5: FIRE-FIGHTING MEASURES

EXTINGUISHING MEDIAThis product is non-flammable. Use fire-extinguishing media appropriate for

surrounding materials.

SPECIFIC HAZARDS Decomposition may product carbon dioxide and carbon monoxide.

PROTECTIVE MEASURES IN FIREUse protective equipment appropriate for firefighting surrounding materials.

SECTION 6: ACCIDENTAL RELEASE MEASURES

PERSONAL PRECAUTIONS Safety Data Sheets are intended for occupational-use products; this is a

consumer-use product. Wear protective clothing as described in Section 8 of

Revision Number: 3.0

this Safety Data Sheet.

CONTAINMENT AND CLEAN-UP MEASURESMop or wipe up any spilled material. For larger spills, dike around spill and

pump to intact container. Dispose of material and container in accordance

with federal, state, and local regulations.

SECTION 7: HANDLING AND STORAGE

USAGE PRECAUTIONS Isolate from strong oxidizing agents and strong acids, including acidic clays,

peroxides, halogens, vinyl chloride, and iodine pentafluoride..

Keep product from freezing.

Keep out of reach of children and pets.
When handling, DO NOT eat, drink, or smoke.
Keep container closed when not in use.

Wash hands with soap and water after handling.

STORAGE PRECAUTIONS Keep separated from strong oxidizing agents and strong acids, including acidic

clays, peroxides, halogens, vinyl chloride, and iodine pentafluoride.

Keep product from freezing. Store in a cool, dry place.

Keep away from food, drink, and animal feed. Keep out of reach of children and pets.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

PROTECTIVE EQUIPMENTWhen using this product as directed on the product label, no protective

equipment is necessary. Take normal precautions necessary when using a

CPAP mask wipe product.

RESPIRATORY EQUIPMENTNone required during normal consumer use.

Industrial use: Normal ventilation is adequate.

SKIN PROTECTION None required during normal consumer use.

Industrial use: Wear impervious gloves and suitable protective clothing.

EYE PROTECTION None required during normal consumer use.

Industrial use: Wear safety glasses or chemical goggles.

ENGINEERING MEASURES None.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE Colorless to slight yellow transparent liquid wicked in white wipes.

Not determined.

pH 8.0 – 8.5.

ODOR Characteristic.

ODOR THRESHOLD Not determined.

FREEZING POINT 32°F/0°C.

INITIAL BOILING POINT 213°F/100°C.

FLASH POINT None.

BOILING POINT RANGE

EVAPORATION RATE

FLAMMABILITY

UPPER FLAMMABILITY LIMIT

LOWER FLAMMABILITY LIMIT

VAPOR PRESSURE

VAPOR DENSITY

Not determined.

Not determined.

SPECIFIC GRAVITY 1.050 – 1.070 @ 77°F/25°C. **SOLUBILITY** Completely soluble in water.

PARTITION COEFFICIENT

AUTO IGNITION TEMPERATURE

DECOMPOSITION TEMPERATURE

VISCOSITY

Water-thin.

0.007

SECTION 10: STABILITY AND REACTIVITY

STABILITY This product is stable under normal conditions.

Do not freeze this product.

Do not store with strong oxidizers and strong acids, including acidic clays,

peroxides, halogens, vinyl chloride, and iodine pentafluoride.

Product may darken over time. This color change does not present a safety

issue.

OTHER Hazardous polymerization will not occur.

Hazardous decomposition of by-products: carbon dioxide (CO₂), carbon

monoxide (CO).

Incompatible materials: strong oxidizers and strong acids, including acidic

clays, peroxides, halogens, vinyl chloride, and iodine pentafluoride.

SECTION 11: TOXICOLOGICAL INFORMATION

ACUTE TOXICITY

INGESTIONUnknown.SKIN CONTACTUnknown.EYE CONTACTUnknown.

OTHER INFORMATION ON TOXICOLOGICAL EFFECTS

TOXIC DOSE LD50 Unknown.

Revision Number: 3.0

Citrus II CPAP Wipes Safety Data Sheet

SECTION 12: ECOLOGICAL INFORMATION

ECOTOXICITY DATANo data available.

MOBILITY This product is completely soluble in water.

DEGRADABILITY No data available.

SECTION 13: DISPOSAL CONSIDERATIONS

PRODUCT DISPOSALDispose of in accordance with federal, state, and local regulations.

SECTION 14: TRANSPORT INFORMATION

US DOT This product is not regulated for transport.

SECTION 15: REGULATORY INFORMATION

OSHA HAZARD COMMUNICATION RULE, 29 CFR 1910.1200:

The Occupational Safety and Health Administration requires Safety Data Sheets to provide any hazards that may be associated with the product and make this information available in the workplace. Since the use pattern and exposure in the workplace are generally not consistent with those experienced by consumers, this SDS may contain additional health hazard information not pertinent to consumer use.

FDA/EPA: This product is not regulated by the FDA or EPA.

TSCA: All chemicals used in this product are either listed on or exempt from the TSCA Inventory.

OSHA Hazardous Substance: This product contains no chemicals listed as OSHA Hazardous Substances.

Proposition 65: None.

SECTION 16: ADDITIONAL INFORMATION

SDS REVISION NUMBER 3.0

SDS PREPARATION DATE March 17, 2014 **SUPERCEDES** 2.0; February 8, 2013

This information is supplied under the OSHA Hazard Communication Standard, 29 CFR 1910.1200, and is offered in good faith based on data available that is believed to be true and accurate. The recommended hygiene and safe handling procedures are believed to be generally applicable to this product. However, each user should review these recommendations in a specific context of the intended use and determine whether they are appropriate for that use. No warranty, express or implied, is made regarding the accuracy of this data, the hazards connected with the use of the product, or the results to be obtained from the use thereof. We assume no responsibility for damage or injury from the use of the product described herein.

Revision Number: 3.0