
AeroChamber miniTM Holding Chamber (HC)

Small Design, Big Results

The **AeroChamber miniTM** HC is a unique device designed to deliver pMDI aerosols to patients in conjunction with:

- mechanical ventilation
- manual resuscitation bag
- standard aerosol resuscitation mask



Anti-Static Chamber

The **AeroChamber miniTM** HC is comprised of an anti-static, dual chamber design. The patented exhalation channel has a one way valve that preserves aerosol dose and minimizes dead space.

Designed with the Patient in Mind

The **AeroDockTM** canister port makes the **AeroChamber miniTM** HC a device designed to accept standard pMDIs as well as those with integrated dose counters. The counters will activate with this feature.

Effective Treatment

The **AeroChamber miniTM** HC is ergonomically correct and fits comfortably in your hand. It can be easily connected to the endotracheal tube in front of a wye or HME in a ventilator circuit. While in a ventilator circuit it is designed to maintain positive end expiratory pressure. It can be used with a standard aerosol resuscitation mask and is effective for use with low tidal volume patients.

Of course it's from Monaghan.TM

Copyright © Monaghan Medical Corporation 1997-2013. All rights reserved. The use of our site is governed by our Terms of Use and Privacy Policy.

TM/_® Trade marks and registered trade marks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation.

† Trade marks and registered trade marks of the respective companies.

Source URL (retrieved on Mar 17 2014 - 1:38pm):

<http://www.monaghanmed.com/products/clinician/aerochamber-mini-holding-chamber-hc>