

**NEW...**

# AirSep External Power Cartridge and Charger

Compatible with AirSep Focus™, FreeStyle™, and FreeStyle™ 5  
Portable Oxygen Concentrators

The new AirSep External Power Cartridge provides extended battery life for all AirSep POCs.

- **NEW LIGHTWEIGHT** alternative to the AirSep AirBelt
- External charger can charge 2 batteries simultaneously
- Built-in belt clip for easy wearability
- Fully charges battery in 5 hours or less
- Easy-access battery gauge displays remaining charge

## Battery Durations:

### Focus

Setting	External Power Cartridge Duration
2	3.5 Hours

Focus adapter cable (CD035-1) must be purchased separately to use the external power cartridge with the Focus.

### FreeStyle

Setting	Duration of FreeStyle Internal Battery + External Power Cartridge
1	10 Hours
2	6 Hours
3	5 Hours

### FreeStyle 5

Setting	Duration of FreeStyle 5 Internal Battery + External Power Cartridge
1	7 Hours
2	4 Hours
3	3 Hours
4	2.5 Hours
5	2 Hours



Charger contains an easy-to-read light system:

**Flashing Green Light:** Charging

**Solid Green Light:** Fully Charged

**Flashing Red Light:** Error

**Solid Red Light:** No battery inserted

## ORDERING INFORMATION

BT024-1	AirSep External Power Cartridge
BT025-1	Charger
BT026-1	External Power Cartridge & Charger Combination (US)
BT026-2	External Power Cartridge & Charger Combination (EUROPEAN)
BT026-3	External Power Cartridge & Charger Combination (AUSTRALIAN)
BT026-4	External Power Cartridge & Charger Combination (UK)

Additional batteries may be purchased independently.