AL835 High Capacity Truck Lift





Direct from Door to Truck Bed

Higher capacity model for standard and extended cab trucks provides lifting capacity to handle manual chairs, scooters and power chairs.

Specifications	AL835
Accommodates	Power chairs & scooters
Lifting Capacity	350 lbs. (160 kg)
Installed Weight	116 lbs. (53 kg)
Powered	Lifting & rotation
Adjustable Boom Height	34" - 41"
Adjustable Boom Length	28.75" – 46.5"
Securement	Docking device included
Available Options	Battery pack, wireless remote
Warranty	3-year transferable

Features & Benefits

- Delivers chair close to driver or passenger side door
- Powered lifting and rotation
- Offset base mounting fits tightly to sidewall and fender walls
- Driver or passenger-side mounting bases available
- Only two motors to control with remote hand control
- Dashboard-mounted on/off switch
- Dual-strapped hook prevents mobility device from spinning during load/unload
- Remote hand control included; wireless option available
- Docking device included



800-833-0478 | www.harmar.com

