

		Procedure	Typical Incision Size	Overall Island Dressing Size*	Medipore™ +Pad Soft Cloth Dressing	Tegaderm™ +Pad Film Dressing
CARDIOVASCULAR/ THORACIC		<b>Coronary Artery Bypass Graft (CABG)</b>	Traditional sternal: Avg. 6–8" (15–20cm) up to 12" (31 cm)	3½" x 10"	3571	3591
				3½" x 13¾"	3573	3593
			Minimally invasive: 3–4" (8–10 cm)	3½" x 6"	3569	3589
			Vein harvesting: Several 1–½" (4 cm) incisions	3½" x 8"	3570	3590
			2¾" x 4"	3564	3584	
		<b>Peripheral Stent</b>	Very small (approx. 1–3cm), sometimes in groin	2" x 2¾"	3562	3582
				2¾" x 4"	3564	3584
GENERAL SURGERY		<b>Appendectomy</b>	Lap: 1" (2.5 cm)	2" x 2¾"	3562	3582
				2¾" x 4"	3564	3584
				3½" x 4"	3566	3587(oval)
			Open incision: 3–6" (8–15 cm)	3½" x 6"	3569	3589
				3½" x 8"	3570	3590
		<b>Breast Biopsy</b>	2" (5 cm)	2¾" x 4"	3564	3584
				3½" x 4"	3566	3586 3587 (oval)
		<b>Cholecystectomy</b>	Lap: 4 to 5 incisions ½" (1 cm) in size	2" x 2¾"	3562	3582
				2¾" x 4"	3564	3584
			Open incision: 3–6" (8–15 cm)	2¾" x 4"	3564	3584
				3½" x 4"	3566	3586
				3½" x 6"	3569	3589
3½" x 8"				3570	3590	
	<b>Gastric Bypass</b>	Lap procedure: Four 1" (2.5 cm)	2" x 2¾"	3562	3582	
			2¾" x 4"	3564	3584	
		Midline: 6–10" (15–25 cm)	3½" x 4"	3566	3586	
			3½" x 8"	3570	3590	
			3½" x 10"	3571	3591	
			3½" x 13¾"	3573	3593	
	<b>Inguinal Hernia</b>	Lap: 1" (2.5 cm)	2" x 2¾"	3562	3582	
			2¾" x 4"	3564	3584	
		Open: up to 4" (10 cm)	3½" x 4"	3566	3586	
			3½" x 6"	3569	3589	
			3½" x 8"	3570	3590	
	<b>Mastectomy</b>	Avg. 6" (15 cm), depending on the procedure	3½" x 8"	3570	3590	
			3½" x 10"	3571	3591	
			6" x 6"	3568	3588	
	<b>Partial Colectomy</b>	4" (10 cm) transverse abdominal or 7–8" (18–20 cm) midline	3½" x 6"	3569	3589	
			3½" x 8"	3570	3590	
			3½" x 10"	3571	3591	
OB-GYN		<b>Cesarean Section</b>	6–8" (15–20 cm) bikini cut, pubic; lower abdominal	3½" x 10"	3571	3591
				3½" x 13¾"	3573	3593
		<b>Hysterectomy</b>	Abdominal: 5–7" (13–18 cm)	3½" x 8"	3570	3590
				3½" x 10"	3571	3591
Laparoscopic: Three ½" (1 cm) LAVH: Two ½" (1 cm)			2" x 2¾"	3562	3582	
			2¾" x 4"	3564	3584	
			3½" x 4"	3566	3586	
			3½" x 6"	3569	3589	
			3½" x 8"	3570	3590	
			3½" x 10"	3571	3591	
ORTHOPEDIC		<b>Carpal Tunnel</b>	1–3" (2.5–8 cm)	2¾" x 4"	3564	3584
				3½" x 4"	3566	3586
				3½" x 6"	3569	3589
		<b>Fracture Fixation</b>	Varies, depending on location of fracture. Often approx. 4" (10 cm) long. May have 2 incisions.	3½" x 6"	3569	3589
				3½" x 8"	3570	3590
		<b>Hip Replacement</b>	8" (20 cm) over incision; some 8–12" (20–31 cm), or as small as 4–6" (10–15 cm)	3½" x 8"	3570	3590
				3½" x 10"	3571	3591
				3½" x 13¾"	3573	3593
		<b>Knee Replacement</b>	4–6" (10–15 cm) total knee vertical incision, arthroplasty	3½" x 6"	3569	3589
				3½" x 8"	3570	3590
3½" x 10"				3571	3591	
	<b>Spine</b>	Laminectomy spinal infusion 2–6" (5–15 cm)	2¾" x 4"	3564	3584	
			3½" x 4"	3566	3586	
			3½" x 6"	3569	3589	
			3½" x 8"	3570	3590	
PLASTIC SURGERY		<b>Breast Augmentation or Reduction</b>	3" (8cm) vertical and/or 3" 'smile' under breast reduction	2¾" x 4"	3564	3584
				3½" x 4"	3566	3586
				3½" x 6"	3569	3589
		<b>Liposuction</b>	Laparoscopic procedure: Can be up to 12 incisions, ¼"–½" (½–1 cm)	2" x 2¾"	3562	3582
			2¾" x 4"	3564	3584	

# dress for SUCCESS

## Suggested Island Dressings for Surgical Procedures

**3M™ Tegaderm™ and 3M™ Medipore™ +Pad dressings are all-in-one sterile dressings designed:**

- for ease of use
- to reduce application time
- to conform to body contours, providing comfort
- with breathable film to allow skin to function normally with good exchange of moisture vapor and oxygen
- to conform to body contours, and to flex with movement or swelling

**Both feature non-adherent, absorbent pads:**

- to promote normal healing
- with less pain and trauma
- for light to moderately draining wounds
- won't stick to the wound

## 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad

**Film Barrier** – provides both a viral\*\* and a bacterial barrier. Requires fewer dressing changes, while helping to reduce the risk of surgical site infections

**Transparent Border** – allows monitoring of peri-wound skin without removing dressing

**Absorbent, Non-Adherent Pad** – for light to moderately draining wounds; won't stick to the wound, allowing normal wound healing process

**Waterproof** – patient may sponge bathe or shower without removing dressing

**Versatile** – available in 9 different sizes

## 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

**Gentle to skin** – provides excellent adhesion, yet is gentle to the skin

**All-in-one design** – reduces application time

**Absorbent, Non-Adherent Pad** – for light to moderately draining wounds; won't stick to the wound, allowing normal wound healing process

**Breathable** – better exchange of moisture vapor and oxygen, and reduced risk of skin maceration

**Versatile** – available in 8 different sizes

\*Choose a dressing size large enough to provide a margin that adheres to dry healthy skin around the incision site.

\*\*In vitro testing shows that the transparent film dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.