



**RELIEF FOR PAIN SYMPTOMS** 


## PLACEMENT GUIDE FOR 2 WIRELESS ELECTRODE CONTROLLERS


The placement of electrodes can be one of the most important parameters in achieving success with therapy. Apply electrodes to the exact site indicated by your physician or healthcare professional. The skin at the application site must be clean and dry. Press electrodes firmly to the skin, making good contact between the skin and the electrodes.


Never remove the self-adhesive electrodes from the skin while the device is turned on.






















**RELIEF FOR PAIN SYMPTOMS** 

## PLACEMENT GUIDE FOR 1 WIRELESS ELECTRODE CONTROLLER

The placement of electrodes can be one of the most important parameters in achieving success with therapy. Apply electrodes to the exact site indicated by your physician or healthcare professional. The skin at the application site must be clean and dry. Press electrodes firmly to the skin, making good contact between the skin and the electrodes.

Never remove the self-adhesive electrodes from the skin while the device is turned on.

