



VIP Solutions, LLC



Phone: 800-726-8620 / 701-293-9175

Fax: 701-297-9702

3309 Fiechtner Drive unit 3 Fargo ND 58103

Website: aquatrek2.com

AQ-1000 – Beach & Trail- Wheelchair

Specifications

Height- 49"

Seating width- 21"

Length-51"

Foot-rest-12"x 11"

Width- 46"

Seat Height- 25"

Seat-backrest-27" tall

Weight Capacity-350 lbs.

Distance between Front wheels- 38"

Distance between arm rest-21 1/4"

Brake-has hitch-pins to lock it in place when parked

Wheels Can withstand surface temps 25°F to 165°F

The AQ-1000 wheelchair-Comes equipped with a Brake & a cushion for the seat & back

Additional Options-

Seat-belt or Seat-belt with Shoulder harness

Will not fit through standard doors-Requires double-doors for entry

AQ-1000-
Showing brakes extended



Shown with Brake down

