



Heavy-Duty Aspirator

The **Roscoe Heavy-Duty Aspirator** is the perfect solution for all your suction needs. It is simple to use and has a durable design to withstand the rigors of homecare use. It produces a maximum vacuum of 560 mmHg and features a high-quality vacuum regulator with anti-vibration setting, an anti-bacteria filter design helps prevent contamination, oil-less, permanent lubrication, and a collection bottle with overflow protection. Designed for easy portability at just 12 pounds, its easy to carry and easy to clean.

PRODUCT FEATURES:

- Filter design helps prevent contamination
- High quality vacuum regulator and anti-vibration vacuum gauge setting
- Oil-free, permanent lubrication
- Collection bottle with overflow protection
- Superior performance and high reliability for continuous operation
- Easy carry and cleaning
- HCPCS Code: E0600*

PRODUCT SPECIFICATIONS:

Dimensions:	14.8" L x 6.8" W x 9.8" H
Weight:	12.1 lbs
Maximum Vacuum:	560 +/- 3% mm Hg
Flow Range:	>40 LPM
Sound Level:	<58 dBA
Electrical Requirements:	AC 115V 60 Hz
Power Consumption:	70 Watts
Filter Pore Size:	0.45 microns
Collection Capacity:	800 cc
Operating Temperature:	32° - 104°F
Operating Humidity:	95% RH
Warranty:	Limited Lifetime

OPTIONAL ACCESSORIES:



Item #	Description	Qty
ASP-HD	Roscoe Heavy-Duty Aspirator	each
90386	Suction Kit (Includes elbow, filter, collection jar, 6' tubing and two 7" tubes)	each
Replacement Parts:		
50038	800 cc Disposable Collection Jar	12/cs
ILB-ROS	Filter, 3/8" to 1/2" barbed	each
COMP-GAUGE	Replacement Gauge	each
ROS-610E	Elbow	20/cs

* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.