


Free U.S. Shipping on the purchase of an Eagle Health Bench/Chair ([/pages/shipping](#))  
 PLUS 10% OFF ANY BENCH/CHAIR (<https://www.eaglehealth.com/pages/10-discount-for-benches-and-chairs>)  
 OR  
 20% OFF OUR 78668 ERGO SWIVEL TRANSFER BENCH ([/pages/78668-special](#))  
 (Limited Time Only)



**\*\*Disclaimer: Armrest is not included with bench.  
 Must be purchased separately\*\***



**EAGLE HEALTH SUPPLIES INC.**  
[www.eaglehealth.com](http://www.eaglehealth.com)



**Clamp Specifications**  
 Fits Tub Wall Widths of 2.75" - 5.75"  
 \*\*Must have flat-top tub wall\*\*  
 Fits Tub Wall Heights of 13" - 17.5"  
 (measure from inside tub)

**Tub Requirements**  
 Min. Width Required Inside Tub (At Top) 23" or wider  
 Min. Width Required Inside Tub (At Bottom) 16.5" or wider  
 Min. Width Required from Outer Tub Wall to Inner Tub Wall (At Bottom) 25"  
 Min. Width Required from Outer Tub Wall to Handlebar (At Top) 29.5"



**Item #77762: TUB-MOUNT SWIVEL SLIDING TRANSFER BENCH**