

PRODUCT DETAILS

| | | | |
|------------------------------------|-----------------------------------|--------------------|---|
| Thickness | Micro | Cuff Length | Standard |
| Material | Natural Rubber Latex | Freedom From Holes | • In process control before packing: 0.65 (Inspection level GI) • Manufacturing final release: 0.65 (Inspection level GI) |
| Not Made From Natural Rubber Latex | No | Allergy Prevention | None |
| Powder Content | Powder-Free | Finger Thickness | 0.19 ± 0.02 / 7.48 ± 0.79 |
| External Glove Surface | Smooth with Micro-Textured Finish | Palm Thickness | 0.19 ± 0.03 / 7.48 ± 1.18 |

PRIMARY INDUSTRIES

- Healthcare

STANDARDS AND CERTIFICATIONS

EN ISO 374-1:2016
Type B



KPT

EN ISO 374-5



VITUS

EN ISO 21420:2020

ASTM D3577;

ASTM D7160

EN 16523-1;

EN 374-2 and -4

ISO 10282

EN 455 1-4