




Product Catalog

2021/2022


MY COLD
THERAPY
Authorized Distributor

Call (704) 953-2545

Order online at: mycoldtherapy.com



COLD THERAPY



Polar Care Wave Shown

Polar Care Wave™

Polar Care Wave combines motorized cold therapy with active compression in a simple and compact system, making it ideal for both facility and home use. Polar Care Wave is easily transported from the hospital, physical therapy clinic or athletic training room to the patient's home, extending the benefits of cold and compression.

Polar Care Wave is designed to reduce pain and edema to optimize patient recovery. Federal law restricts this device to sale by or on the order of a licensed health care practitioner.

Simple

- Single patient use facility-to-home unit
- Remarkably easy to operate
- No additional barrier needed
- Ice bag accessory available for an easier way to refill system

Customizable

- Varied levels of cold + compression settings
- Y-adapter for bilateral usage
- Multiple pad configurations

Convenient

- Compact design facilitates transport
- System alone and system/pad combos available
- Hose-management storage system
- Whisper quiet pump technology encourages rest and recovery

Quality

- High-quality Breg Cold Compression Pads
- Class II 510(k) clearance



Description	Part #
Polar Care Wave	100577-000
Polar Care Wave Combo Cold Compression Knee	100578-000
Polar Care Wave Combo Cold Compression Shoulder	100579-000
Polar Care Wave Combo Cold Compression Universal	100580-000

Description	Part #
Accessories	
Breg Carrying Bag	C00015
PC Wave Knee Y Adapter	100581-000
Breg Ice Bags (kit of 12)	100582-000
PC Wave Replacement Power Supply	100584-000

Description	Part #
Pads	
Breg Cold/Compression Pad, Knee, L	C00003
Breg Cold/Compression Pad, Shoulder	C00004
Breg Cold/Compression Pad, Foot/Ankle	C00005
Breg Cold/Compression Pad, Hip	C00013
Breg Cold/Compression Pad, Universal	C00016
Breg Cold/Compression Pad, Knee	C00017
Breg Cold/Compression Pad, Back	C00020



Standard Knee Pad



Hip Pad



Shoulder Pad



Back Pad



Large Knee Pad



Foot/Ankle Pad



Universal Pad

Breg Cold/Compression Pad Sizes				
Description	Max Circumference	Min Circumference	Dimensions	Strap Lengths
Standard Knee	27" (68.58 cm)	8" (20.32 cm)	16" x 12" (40.64 x 30.48 cm)	11" (27.94 cm)
Shoulder Pad	50" (127 cm) shoulder	40" (101.6 cm) shoulder	15" x 13" (38.1 x 33.02cm)	34" (86.36 cm)
	20" (50.8 cm) bicep	12" (30.48 cm) bicep		9.5" (24.13 cm)
Foot/Ankle Pad	16" (40.64 cm) ankle	9" (22.86 cm) ankle	17"x12" (43.18 x 30.48 cm)	
	12" (30.48 cm) foot	7" (17.78 cm) foot		
Back Pad	53" (134.62 cm) waist	32.5" (82.55 cm) waist	13.5"x14.5" (34.29 x 36.83 cm)	21" (53.34 cm) on each side
Hip Pad	53" (134.62 cm) waist	33" (83.82 cm) waist	16"x12" (40.64x30.48 cm)	20" (50.8 cm) each (2 straps)
	31" (78.74 cm) thigh	8" (20.32 cm) thigh		25" (63.5 cm)
Large Knee Pad	30.5" (77.47 cm) thigh	16" (40.64 cm) thigh	16"x17.5" (40.64 x 44.45 cm)	14.5" & 11.5" (36.83 & 29.21 cm)
Universal Pad	49" (124.46 cm) chest/shoulder	8" (20.32 cm) at the knee	13"x12.5" (33.02 x 31.75 cm)	
	20" (50.8 cm) bicep/calf strap	31" (78.74 cm) at the chest/shoulder		
	47" (119.38 cm) for the thigh strap			

Polar Care® Kodiak®

The Polar Care Kodiak Cold Therapy System is a convenient and versatile offering in Breg's Polar Care line. Its easy-to-use, compact design makes it great for clinic, hospital, and home use. With the addition of a little ice and water, your patients will enjoy 6-8 hours of effortless cold therapy. Offering the only battery powered option on the market, Breg ensures your patients can enjoy the benefits of cold therapy from anywhere, on the sidelines or in the backyard. Each battery pack comes with four replaceable AA batteries ready to power 10-14 hours of motorized cold therapy. The revolutionary Intelli-Flo® pads specify temperature to each treatment area. Proper use requires an insulation barrier between the Intelli-Flo pad and the patient's skin. Federal law restricts this device to sale by or on the order of a licensed health care practitioner.

Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post trauma
- Chronic pain
- Physical therapy

Features

- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Whisper-quiet pump technology allows patients to rest and recover in comfort
- Durable, powerful construction provides peace of mind
- Battery Pack option for patients on the go
- Intelli-Flo Pad technology meets the unique demands of the treatment area
- Sterile Dressings available for each part



Polar Care Kodiak Battery Pack (optional)

Description	Part #
Polar Care Kodiak*	10601
Kodiak Combo Units (includes Kodiak and 1 single-patient use pad)	
PC Kodiak with Intelli-Flo Multi-Use Pad <i>Sterile Intelli-Flo Dressing, Multi-Use</i>	10602 10630
PC Kodiak with Intelli-Flo Multi-Use, Long Stem Pad (16" long) <i>Sterile Intelli-Flo Dressing, Multi-Use</i>	10603 10630
PC Kodiak with Intelli-Flo Knee Pad <i>Sterile Intelli-Flo Dressing, Multi-Use</i>	10604 10630
PC Kodiak with Intelli-Flo Knee, Compression Pad <i>Sterile Intelli-Flo Dressing, Multi-Use</i>	10605 10630
PC Kodiak with Intelli-Flo Shoulder Pad <i>Sterile Intelli-Flo Dressing, Multi-Use</i>	10606 10630
PC Kodiak with Intelli-Flo XL Shoulder Pad <i>Sterile Intelli-Flo Dressing, XL Shoulder</i>	10607 10650
PC Kodiak with Intelli-Flo Ankle Pad <i>Sterile Intelli-Flo Dressing, Ankle</i>	10608 10640
PC Kodiak with Intelli-Flo Ankle, Compression Pad <i>Sterile Intelli-Flo Dressing, Ankle</i>	10609 10640
PC Kodiak with Intelli-Flo Back Pad <i>Sterile Intelli-Flo Dressing, Back</i>	10611 09810
PC Kodiak with Intelli-Flo Hip Pad <i>Sterile Intelli-Flo Dressing, Hip</i>	10621 10680

*Contact customer care for CE Marking international part numbers.

Note: Matching Sterile Intelli-Flo Dressing listed under each combo unit.

Description	Part #
Bulk 6-Packs	
PC Kodiak, Bulk 6-Pack	28027
PC Kodiak Combo, Intelli-Flo Multi-Use, Bulk 6-Pack	28028
PC Kodiak Combo, Intelli-Flo Knee, Bulk 6-Pack	28029
PC Kodiak Combo, Intelli-Flo Knee Compression, Bulk 6-Pack	28030
PC Kodiak Combo, Intelli-Flo Shoulder, Bulk 6-Pack	28031
Accessories	
Kodiak Battery Pack, includes 4 AA batteries	97050
Replacement Power Supply	10698
Intelli-Flo Pad Y Adaptor	10699
Compression Strap (60" long)	10670
Intelli-Flo Hip Pad Strap	10690
Replacement Kit	10690
Polar Care Bed Hanger	000038-000

Intelli-Flo® Pads

These revolutionary pads offer temperature specific to each treatment area, allowing for consistent cold delivery. The ergonomic design provides exceptional coverage, static compression, and patient comfort. Intelli-Flo Pads are compatible with Polar Care Kodiak only. Proper use requires an insulation barrier between the Intelli-Flo Pads and the patient's skin. The water impermeable Sterile Polar Dressings offered by Breg provide an appropriate and complete barrier between the pad and the patient's skin.



Intelli-Flo Multi-Use Pad
(13.5" x 14.25")



Intelli-Flo Knee Pad
(13.5" x 14")



Intelli-Flo Shoulder Pad
(13.5" x 14")
X-Large (13" x 21")



Intelli-Flo Ankle Pad
(12.75" x 18.5")
Small (8.25" x 9.75")
Part number 10210 shown



Intelli-Flo Back Pad
(8.25" x 11.25")



Intelli-Flo Hip Pad
(13" x 20")



Intelli-Flo Hand / Wrist Pad
(8.25" x 9.75")



Shown on Ankle
Intelli-Flo 3x5 Pad
(3" x 5")



Shown on Hand
Intelli-Flo 3x5 Pad
(3" x 5")



Shown on Foot
Intelli-Flo 3x5 Pad
(3" x 5")

Description	Part #
Intelli-Flo Pad, Multi-Use	10240
Sterile Intelli-Flo Dressing, Multi-Use	10630
Intelli-Flo Pad, Multi-Use, Long Stem (16" long)	10245
Sterile Intelli-Flo Dressing, Multi-Use	10630
Intelli-Flo Pad, Knee	10230
Sterile Intelli-Flo Dressing, Multi-Use	10630
Intelli-Flo Pad, Knee, Compression	10235
Sterile Intelli-Flo Dressing, Multi-Use	10630
Intelli-Flo Pad, Shoulder	10220
Sterile Intelli-Flo Dressing, Multi-Use	10630
Intelli-Flo Pad, XL Shoulder	10225
Sterile Intelli-Flo Dressing, XL Shoulder	10650
Intelli-Flo, Ankle*	10210
Sterile Intelli-Flo Dressing, Ankle	10640

Description	Part #
Intelli-Flo, Small Ankle*	10211
Sterile Intelli-Flo Dressing, Back, Ankle, Wrist	09810
Intelli-Flo Pad, Back	10250
Sterile Intelli-Flo Dressing, Back, Ankle, Wrist	09810
Intelli-Flo Pad, Hip	10280
Sterile Intelli-Flo Dressing, Hip	10680
Intelli-Flo Pad, Hand / Wrist*	10260
Sterile Intelli-Flo Dressing, Back, Ankle, Wrist	09810
Intelli-Flo Pad 3x5*	10205
Sterile Intelli-Flo Dressing, 3x5	10660

*Dressing is included

Note: Matching Sterile Intelli-Flo Dressing listed under each Intelli-Flo Pad.

Contact customer care for CE Marking international part numbers.

Polar Care Cube™

The Polar Care Cube Cold Therapy System is simple and reliable. Its easy-to-use, compact design makes it great for clinic, hospital, and home use. With the addition of ice and water, your patients will enjoy 6-8 hours of effortless cold therapy. Proper use requires an insulation barrier between the pad and the patient's skin. Federal law restricts this device to sale by or on the order of a medical professional.

Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post trauma
- Chronic pain
- Physical therapy

Features

- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Whisper-quiet pump technology allows patients to rest and recover in comfort
- Durable, powerful construction provides peace of mind
- Broad selection of premium WrapOn Pads conform to the affected area
- Sterile Polar Dressings available for each part



Description	Part #
Polar Care Cube	10701
Cube Combo Units with WrapOn Pad (includes Cube and 1 single-patient use pad)	
PC Cube with WrapOn Knee Pad <i>Sterile Polar Dressing, Knee L</i>	10705 04708
PC Cube with WrapOn L Knee Pad <i>Sterile Polar Dressing, Knee L</i>	10706 04708
PC Cube with WrapOn XL Knee Pad <i>Sterile Polar Dressing, Shoulder</i>	10707 04908
PC Cube with WrapOn Ankle Pad <i>Sterile Polar Dressing, Back, Ankle, Wrist</i>	10708 09810
PC Cube with WrapOn Multi-Use XL Pad with Long Stem (16" Long) <i>Sterile Polar Dressing, XL Multi-Use</i>	10709 02344
PC Cube with WrapOn XL Multi-Use Pad <i>Sterile Polar Dressing, XL Multi-Use</i>	10710 02344
PC Cube with WrapOn Shoulder Pad <i>Sterile Polar Dressing, Shoulder</i>	10711 04908
PC Cube with WrapOn XL Shoulder Pad <i>Sterile Polar Dressing, XL Shoulder</i>	10712 04918
PC Cube with WrapOn Hip Pad <i>Sterile Intelli-Flo Dressing, Hip</i>	10713 10680

Note: Matching Sterile Polar Dressing listed under each pad. Body part indicated for Sterile Polar Dressing may or may not correlate with body part in pad name, but the sizes have been paired appropriately.

Description	Part #
Bulk Packs	
Polar Care Cube, 4-Pack	28600
PC Cube Combo, WrapOn XL Multi-Use, 4-Pack	28601
PC Cube Combo, WrapOn Knee, 4-Pack	28602
PC Cube Combo, WrapOn Shoulder, 4-Pack	28603
PC Cube Combo Units with Polar Pad (includes Cube and 1 single-patient use pad)	
PC Cube with Polar Knee / Shoulder Pad <i>Sterile Polar Dressing, Knee / Shoulder</i>	10702 02328
PC Cube with Polar Multi-Use Pad <i>Sterile Polar Dressing, Multi-Use</i>	10704 02348
Accessories	
Replacement Power Supply	10698
Compression Strap (60" long)	10670
PC Y Adaptor	07640
Polar Care Bed Hanger	000038-000
4ft Extension Tube	09915

Polar Care Glacier™

The Polar Care Glacier Cold Therapy System balances the needs of the patient and the clinician. Its easy-to-use, feature-rich design makes it a great choice for hospital use. With the addition of ice and water, your patients will enjoy up to 13 hours of cold therapy according to their own comfort level. Proper use requires an insulation barrier between the pad and the patient's skin. Federal law restricts this device to sale by or on the order of a medical professional.

Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post trauma
- Chronic pain
- Physical therapy

Features

- Once-a-shift fill requirement increases a hospital's efficiency— 8hr shifts? 12hr shifts? The Glacier has them covered
- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Whisper-quiet pump technology allows patients to rest and recover in comfort
- Temperature control dial allows physicians to individualize protocols
- Broad selection of premium WrapOn Pads conform to the affected area
- Sterile Polar Dressings available for each part



Description	Part #
Polar Care Glacier	10901
PC Glacier Combo Units with WrapOn Pad (includes Glacier and 1 single-patient use pad)	
PC Glacier with WrapOn Knee Pad <i>Sterile Polar Dressing, Knee L</i>	10905 04708
PC Glacier with WrapOn Knee L Pad <i>Sterile Polar Dressing, Knee L</i>	10906 04708
PC Glacier with WrapOn Knee XL Pad <i>Sterile Polar Dressing, Shoulder</i>	10907 04908
PC Glacier with WrapOn Multi-Use XL Pad <i>Sterile Polar Dressing, XL Multi-Use</i>	10908 02344
PC Glacier with WrapOn Multi-Use XL Pad with Long Stem <i>Sterile Polar Dressing, XL Multi-Use</i>	10909 02344
PC Glacier with WrapOn Shoulder Pad <i>Sterile Polar Dressing, Shoulder</i>	10911 04908
PC Glacier with WrapOn Shoulder XL Pad <i>Sterile Polar Dressing, XL Shoulder</i>	10912 04918
PC Glacier with WrapOn Ankle Pad <i>Sterile Polar Dressing, Back, Ankle, Wrist</i>	10913 09810
PC Glacier with WrapOn Hip Pad <i>Sterile Intelli-Flo Dressing, Hip</i>	10914 10680

Description	Part #
PC Glacier Combo Units with Polar Pad (includes PC Glacier and 1 single-patient use pad)	
PC Glacier with Knee / Shoulder Pad <i>Sterile Polar Dressing, Knee / Shoulder</i>	10902 02328
PC Glacier with Multi-Use XL Pad <i>Sterile Polar Dressing, XL Multi-Use</i>	10904 02348
Accessories	
Replacement Power Supply	10698
Compression Strap (60" long)	10670
PC Y Adapter	07640
4ft Extension Tube	09915
Polar Care Bed Hanger	000038-000
<i>Note: Matching Sterile Polar Dressing listed under each Pad. Body part indicated for Sterile Polar Dressing may or may not correlate with body part in pad name, but the sizes have been paired appropriately.</i>	

Polar Care Cub

The manually powered Polar Care Cub Cold Therapy System is a simple and cost-effective solution for patients and practitioners looking for an effective alternative to rudimentary ice therapy methods. The Cub uses a manually powered hand pump to deliver cold therapy to the treatment area. With the addition of ice and water, hours of pain relief are just a few squeezes of the pump away.

Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post trauma
- Chronic pain
- Physical therapy

Features

- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Small, lightweight, portable and manually powered design ensures the Cub can go anywhere it's needed—on the sidelines, or in the backyard
- Simple, hand-powered design is intuitive
- Broad selection of WrapOn and Polar Pads to cover the affected area (not compatible with WrapOn or Polar Ankle Pads)
- Sterile Polar Dressings available for each part



Description	Part #
Polar Care Cub	04009
Polar Care Cub Combo Units with WrapOn Pad (includes Cub and 1 single-patient use pad)	
PC Cub with WrapOn XL Multi-Use Pad <i>Sterile Polar Dressing, XL Multi-Use</i>	04690 02344
PC Cub with WrapOn XL Multi-Use Pad Long Stem (16" long) <i>Sterile Polar Dressing, XL Multi-Use</i>	04640 02344
PC Cub with WrapOn Knee Pad <i>Sterile Polar Dressing, Knee L</i>	04600 04708
PC Cub with WrapOn Shoulder Pad (universal) <i>Sterile Polar Dressing, Shoulder</i>	11040 04908
Accessories	
PC Cub Pump Assembly	04010

Note: Matching Sterile Polar Dressing listed under each combo unit.

WrapOn Polar Pads

The WrapOn Polar Pads' ergonomic design provides exceptional coverage and patient comfort. Elastic straps offer static compression while holding the pad firmly in place. WrapOn Pads may be used with the Polar Care Glacier, Cube, and Cub only. Proper use requires an insulation barrier between the pad and the patient's skin. The water impermeable Sterile Polar Dressings offered by Breg provide a complete barrier between the pad and the patient's skin.



Multi-Use XL WrapOn Pad
(11.25" x 11.25")



Knee WrapOn Pad
(10.25" x 11.25")
Large (11.5" x 12") XL (12.25" x 19")



Shoulder WrapOn Pad
(10.25" x 11.5")
X-Large (13" x 21.5") (10" x 11.75")



WrapOn Hip Pad
(11.5" x 12")
Long strap length = 44"
Short strap length = 26"



Ankle WrapOn Pad
(8.25" x 9.75")



Back WrapOn Pad
(8.25" x 11")



Hand / Wrist WrapOn Pad
(8.5" x 9.75")

Description	Part #
WrapOn Pad, XL Multi-Use	04790
Sterile Polar Dressing, XL Multi-Use	02344
WrapOn Pad, XL Multi-Use Long Stem (16" long)	04740
Sterile Polar Dressing, XL Multi-Use	02344
WrapOn Pad, Knee	04700
Sterile Polar Dressing, Knee L	04708
WrapOn Pad, L Knee	04703
Sterile Polar Dressing, Knee L	04708
WrapOn Pad, XL Knee	04705
Sterile Polar Dressing, Shoulder L	04908
WrapOn Pad, Shoulder (Universal)	04900
Sterile Polar Dressing, Shoulder L	04908
WrapOn Pad, Shoulder, w/ Extra Long Straps (9", 12" and 29" long)	04903
Sterile Polar Dressing, Shoulder L	04908
WrapOn Pad, XL Shoulder (Universal)	04905
Sterile Polar Dressing, XL Shoulder	04918
WrapOn Pad, Hip	04750
Sterile Intelli-Flo Dressing, Hip	10680

Description	Part #
WrapOn Pad, Ankle*	04730
Sterile Polar Dressing, Back, Ankle, Wrist	09810
WrapOn Polar Pad, Back**	09800
Sterile Polar Dressing, Back, Ankle, Wrist	09810
WrapOn Polar Pad, Back w/ Long Straps (two 17.5" x 5" straps)	09805
Sterile Polar Dressing, Back, Ankle, Wrist	09810
WrapOn Polar Pad, Hand / Wrist***	10750
Sterile Polar Dressing, Back, Ankle, Wrist	09810

Note: Matching Sterile Polar Dressing listed under each WrapOn Pad. Body part indicated for Sterile Polar Dressing may or may not correlate with body part in pad name, but the sizes have been paired appropriately.

*Not Compatible with Polar Care Cub.

**Must be combined with Back Brace PN 1015X or 0733X.

***Dressing is included.

Polar Pads

Polar Pads' ergonomic design provides exceptional coverage and patient comfort. Polar Pads may be used with the Polar Care Glacier, Cube, and Cub only. Proper use requires an insulation barrier between the pad and the patient's skin.



Multi-Use Polar Pad
(9.75" x 11.25")



Multi-Use Polar Pad XL
(11.5" x 11.75")



Rectangle Polar Pad
Small (4.5" x 10.75")
Large (10" x 14.25")



Knee / Shoulder Polar Pad
(9.75" x 11.25")



Ankle Polar Pad
(8.25" x 9.75")



Back Polar Pad
(6.25" x 11.25")



TMJ Polar Pad
(4.25" x 16")

Description	Part #
Polar Pad, Multi-Use	02340
Sterile Polar Dressing, Multi-Use	02348
Polar Pad, Multi-Use, Long Stem (16" long)	02480
Sterile Polar Dressing, Multi-Use	02348
Polar Pad, XL Multi-Use	02346
Sterile Polar Dressing, XL Multi-Use	02344
Polar Pad, XL Multi-Use, Long Stem (16" long)	02486
Sterile Polar Dressing, XL Multi-Use	02344
Polar Pad, Small Rectangle	02400
Sterile Polar Dressing, Large Rectangle	02428
Polar Pad, Large Rectangle	02420
Sterile Polar Dressing, Large Rectangle	02428
Polar Pad, Large Rectangle, Long Stem (16" long)	02520
Sterile Polar Dressing, Large Rectangle	02428

Description	Part #
Polar Pad, Knee / Shoulder	02320
Sterile Polar Dressing, Knee / Shoulder	02328
Polar Pad, Ankle*	02500
Sterile Polar Dressing, Back, Ankle, Wrist	09810
Polar Pad, Back	09900
Sterile Polar Dressing, Back, Ankle, Wrist	09810
Polar Pad, TMJ**	02460
Sterile Polar Dressing, TMJ	02468

Note: Matching Sterile Polar Dressing listed under each Polar Pad.

*Not Compatible with Polar Care Cub.

**Dressing is included.

Polar Wraps

Polar Wraps can be wrapped around each of the Polar Pads, allowing patients to conveniently secure the pads to the affected area without using additional wraps or bandages.

Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post-trauma
- Chronic pain
- Physical therapy

Description	Part #
Polar Wrap, Back M	02703
Polar Wrap, Back L	02704
Polar Wrap, Hip M	02753
Polar Wrap, Hip L	02754

Description	Part #
Polar Wrap, Ankle	02780
Polar Wrap, Knee / Thigh	02600
Polar Wrap, Shoulder S/M	02653
Polar Wrap, Shoulder L/XL	02655



Polar Care Packs

Breg's Polar Care Pack line is a complete line of gel pack products with wraps for the ankle, knee, hip, and shoulder. The premium gel packs are vacuum packed, so the gel stays evenly distributed and pliable at standard freezer temperatures. These unique gel packs are specifically formulated to provide cold therapy for more than two hours. Each wrap features elastic straps that provide static compression while holding the gel packs firmly in place. The Polar Care Pack line is an excellent addition to your motorized cold therapy protocol.

Description	Part #
Knee Wrap Only	02873
Knee Wrap plus (2) Gel Packs*	02874
Knee Wrap plus (4) Gel Packs*	02875
Shoulder Wrap Only	02876
Shoulder Wrap plus (2) Gel Packs*	02877
Shoulder Wrap plus (4) Gel Packs*	02878
Hip Wrap Only	02879
Hip Wrap plus (2) Gel Packs*	02880
Hip Wrap plus (4) Gel Packs*	02881
Ankle Wrap Only	02882
Ankle Wrap plus (2) Gel Packs*	02883
Ankle Wrap plus (4) Gel Packs*	02884
Universal Wrap Only	100419-000

Description	Part #
Universal Wrap plus (2) Gel Packs*	100420-000
Universal Wrap plus (4) Gel Packs*	100421-000
Gel Pack Single	02885
Gel Pack Case (12)	100418-000
Gel Wrap Patient Bag	02887

*Includes drawstring bag.



Soft Stuff

Soft Stuff gel wraps stay cold for 45-60 minutes. Times will vary depending upon size of the patient, the environment and temperature of the freezer.

Features

- Always pliable regardless of the application
- Uses latex-free elastic straps
- All the straps are 3" (8cm) wide, spreading the pressure over a wide area
- The straps that come with the universal knee/shoulder allow the specific pad to be used comfortably on both the knee and shoulder
- All ingredients inside Soft Stuff are food grade products. If there ever was a leak in the nylon bag that holds the cold solution, it would not be a safety or health hazard

Description	Part #
Soft Stuff, Univ, Knee / Shoulder (2 inserts)	SS800
Soft Stuff, Univ, Knee / Shoulder (1 insert)	SS801
Soft Stuff, Univ, Foot / Hand / Elbow (1 insert)	SS811
Soft Stuff, Univ, Foot / Hand / Elbow (2 inserts)	SS812

