

AquaPak® Prefilled Sterile Systems



Designed to add humidity to inspired gas and mitigate drying and irritation



HUD00640

HUD00540

HUDRHP340U

AquaPak® Prefilled Bubble Humidifiers

- Features a micro-diffuser that produces small bubbles and generates surface agitation
- Quiet operation for use at a flow rate of 1-5 LPM
- Available individually or as a kit, which includes water and adaptor
- Reservoir bottles have built-in oxygen tubing connectors
- Prefilled sterile reservoirs come in 3 convenient sizes: 340 mL, 540 mL, 640 mL

Ordering information

Item No.	Size	Sterile Water	Adaptor	Pkg.
HUDRHP340U	340 mL	•	•	20/cs
HUD00540	540 mL	•	•	12/cs
HUD00640	640 mL	•	•	10/cs
HUD00301	340 mL	•		20/cs
HUD00501	540 mL	•		12/cs
HUD00601	640 mL	•		10/cs





HUD03700

AquaPak® Prefilled Large Volume Nebulizers

- Wide variety of reservoirs, adaptors, and solutions to meet specific patient requirements
- Features a return tube that recycles large particles back to the reservoir to be aerosolized, designed to help reduce rainout in tubing
- Separate heater available to provide heated aerosol therapy
- Prefilled Aquapak sterile reservoirs are available in 3 convenient sizes: 440 mL, 760 mL, 1,070 mL

Ordering information

Item No.	Size	Sterile Water	Saline	HUD03128 (028) Adaptor	HUD3133 (033) Adaptor	Pkg.
HUD00400	440 mL	•				20/cs
HUD00405	440 mL		0.45%			20/cs
HUD00428	440 mL	•		•	•	20/cs
HUD00433	440 mL	•				20/cs
HUD4435	440 mL		0.45%	•		20/cs
HUD00470	440 mL		0.45%		•	20/cs
HUD03700	760 mL	•				10/cs
HUD3705	760 mL		0.45%			10/cs
HUD3709	760 mL		0.90%			10/cs
HUDRHP728U	760 mL	•		•		10/cs
HUD03733	760 mL	•			•	10/cs
HUD3735	760 mL		0.45%	•		10/cs
HUD3739	760 mL		0.90%	•		10/cs
HUD03770	760 mL		0.45%		•	10/cs
HUD04000	1,070 mL	•				10/cs
HUD04045	1,070 mL		0.45%			10/cs
HUD04128	1,070 mL	•		•		10/cs
HUD4133	1,070 mL	•			•	10/cs
HUD4135	1,070 mL		0.45%	•		10/cs
HUD03128	1,070 mL			•		10/cs
HUD03133	1,070 mL				•	10/cs

Accessories



HUD00040

Humidifier Adaptor

- Audible alarm system alerts clinicians of flow restriction or occlusion of oxygen cannula or mask tubing
- Use with items HUD00301, HUD00501, HUD00601
- Individually packaged, 120/cs



HUD03128

HUD03128 Nebulizer Adaptor Yellow

- Venturi-style adaptor allows the selection of oxygen concentrations from 28% to 98%
- Individually packaged, 10/cs
- Compatible only with AquaPak Prefilled Large Volume Nebulizer bottles



HUD3133

HUD3133 Nebulizer Adaptor Blue

- Venturi-style adaptor allows the selection of oxygen concentrations from 33% to 98%
- Individually packaged, 10/cs
- Compatible only with AquaPak Prefilled Large Volume Nebulizer bottles



HUD05014

External Electric Heater

- Prefilled nebulizer system provides heated aerosol
- Adjustable aerosol temperature helps meet varying patient needs
- Illuminated power switch for easy on and off
- Individually packaged, 1 ea

Learn more. Contact your Medline Representative, visit medline.com or call 1-800-MEDLINE to learn how we can help with all your humidifier needs.