

## LAB-X SERIES, BARIATRIC, BLOOD DRAWING CHAIRS


Built tough enough to support 700 lbs. under normal use, the Clinton Lab X Series, Bariatric, Blood Drawing Chairs are in a class by themselves. Sturdy, steel construction is combined with thick foam padding to provide an extra measure of comfort and security for larger patients.

### Lab X Bariatric Series Chairs feature:

- A full 35" wide seat (89 cm)
- 1 1/4" square, heavy duty, all-welded, tubular steel frame (3 cm)
- Extra steel center support
- Adjustable height and depth flip arm(s)
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- Optional ClintonClean™ arms available
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray powder-coated frame
- Adjustable foot levelers
- Easy-clean, durable vinyl in a selection of colors
- 700 lbs. combined load capacity under normal use (181.4 kg)



Armrests are adjusted by two star knobs



Adjustable depth flip-arm



66000B



-665 ClintonClean™ Arms option

### Lab X Series, Bariatric, Blood Drawing Chair with Padded Arms

- Padded seat and backrest
- Dual, padded flip arms and armrests
- Option -665, ClintonClean™ solid plastic arms available

#### 66000B

	dimensions	seat height	
seat	35" x 21" 89 cm x 53 cm	20" 51 cm	
overall	width	depth	height
	46" 117 cm	29" 74 cm	36" 91 cm
arms	arms height range		
	6" 15 cm		

**IMPORTANT NOTE:** Clinton Blood Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

### Lab X Series, Bariatric, Blood Drawing Chair with Padded Flip Arm and Drawer

- Padded seat and backrest
- Padded flip arm and armrest
- ClintonClean™ drawer top
- One-piece, all-plastic drawer
- Option -665, ClintonClean™ solid plastic arms available

#### 66002B

	dimensions	seat height	
seat	35" x 21" 89 cm x 53 cm	20" 51 cm	
overall	width	depth	height
	53" 135 cm	29" 74 cm	36" 91 cm
arms	arms height range		
	6" 15 cm		

### Lab X Series, Bariatric, Extra-Tall, Blood Drawing Chair with Padded Arms

- Padded seat and backrest
- Dual, padded flip arms and armrests
- Extra height reduces bending for phlebotomist
- Option -665, ClintonClean™ solid plastic arms available

#### 66099B

	dimensions	seat height	
seat	35" x 21" 89 cm x 53 cm	26" 66 cm	
overall	width	depth	height
	46" 117 cm	29" 74 cm	42 1/2" 108 cm
arms	arms height range		
	6" 15 cm		

### Lab X Series, Bariatric, Extra-Tall, Blood Drawing Chair with Padded Flip Arm and Drawer

- Padded seat and backrest
- Padded flip arm and armrest
- ClintonClean™ drawer top
- One-piece, all-plastic drawer
- Extra height reduces bending for phlebotomist
- Option -665, ClintonClean™ solid plastic arms available

#### 66092B

	dimensions	seat height	
seat	35" x 21" 89 cm x 53 cm	26" 66 cm	
overall	width	depth	height
	53" 135 cm	29" 74 cm	42 1/2" 108 cm
arms	arms height range		
	6" 15 cm		



66002B



66099B



66092B

**WARNING:** Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).