

INVACARE®
9000 SL
WHEELCHAIR



The **Invacare® 9000 SL Wheelchair** is the perfect high-strength lightweight rental product for those who want quality and durability at a moderate price. This low-maintenance wheelchair also has a full complement of options and accessories. Its contemporary design and high-performance frame make the 9000 SL wheelchair a superior value.

Features

- Stylish, contemporary carbon steel frame is long-lasting
- Multiple frame color variations for customizing
- Dual axle positions allow variations of seat-to-floor heights
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Space-saver arms reduce overall frame width without limiting seat width
- Meets RESNA standards for manual wheelchairs in product performance

HCPCS Code: K0004



- Dual axle positions
- Urethane casters
- Space saver arms
- Seating options available
- Builder options available



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-226
111241



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width (2" increments)	Seat Depth	Seat-to-Floor Height (1" increments)
9SL	23", 25"	36"	30"	15"	16"-20"	16", 18"	15.5"-19.5"
9SL	15"-17", 17"-19"	6", 8"	20", 22", 24"	250 lb.	43 lb.	32"x13"x 37"	5 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.