

Compressor nebulizer system

The McKesson compressor nebulizer system is designed to convert liquid medicine into a mist by using compressed air technology. The compressor nebulizer system helps with the treatment of asthma, allergies and other respiratory disorders. It's suitable for children and adults who need nebulized medicine and want a compact, easy-to-use nebulizer.

Features and benefits

- Compressor forces air to the nebulizer
- Converts the prescribed medication into an aerosol of microscopic droplets that can easily be inhaled
- Lightweight (weighs less than 3 lbs.)
- Compact and durable making it easy to store and transport
- Able to deliver a consistent particle size of less than 5 microns
- Treatment times are fast and effective
- Includes five filters
- Can be used with McKesson Nebulizer Mask
Adult: #32642, #32643, #32644, Infant: #32641, Child: #32645
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	CONTENTS	QTY
16-3655LT	Compressor Nebulizer System, Compact, Lightweight, Accessories Included	Compressor Nebulizer Unit, Nebulizer Kit, Inhalation Tube, Corrugated Pipe, T-Type Connector, Inhalation Mouthpiece	1 ea, 6 ea/cs
16-3655LTF	Nebulizer Filters	N/A	5/pack, 20 packs/cs

Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.