Unna Boot/Unna Boot with Calamine

Unna Boot: Gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide to provide venous ulcer compression therapy. Contains no gelatin. Non-raveling, thread-locked gauze helps prevent loose threads that can contaminate wounds.

Unna Boot with Calamine: Gauze bandage that has been evenly impregnated with a non-hardening pasteof zinc oxide and calamine to provide venous ulcer compression therapy and help ease itching.

Indications: Provides venous ulcer compression therapy in actively ambulatory patients. Dries the oozing and weeping of poison ivy, poison oak, or poison sumac. Ideal for managing extremities with conditions such as venous ulcers that require protection and support.

Contraindications: Patients with a known hypersensitivity to the product itself or to its components, including zinc oxide or calamine. Unna boots are not intended to provide compression therapy in non-ambulatory patients.

Change Frequency: Change dressing as directed by a healthcare professional, typically once or twice a week.

Recommended Secondary Dressings: Cohesive bandage.



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY	HCPCS	EXUDATE LEVEL*
Unna Boot 3", with Zinc Oxide, Non-Sterile, 3" x 10 yds	2066	1/bx, 12bx/cs	A6456	N/A
Unna Boot 4", with Zinc Oxide, Non-Sterile, 4" x 10 yds	2067	1/bx, 12bx/cs	A6456	N/A
Unna Boot 3", with Zinc Oxide & Calamine, Non-Sterile, 3" x 10 yds	2068	1/bx, 12bx/cs	A6456	N/A
Unna Boot 4", with Zinc Oxide & Calamine, Non-Sterile, 4" x 10 yds	2069	1/bx, 12bx/cs	A6456	N/A

EXUDATE LEVEL KEY: $\Diamond \Diamond \Diamond = \text{Light}$ $\Diamond \Diamond \Diamond = \text{Moderate}$ $\Diamond \Diamond \Diamond = \text{Heavy}$

*Exudate level indicates highest level of exudate that the dressing can accommodate

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2014 McKesson Medical-Surgical Inc.