

# Tranquility ThinLiner Absorbent Sheets

Super thin and super-absorbent. Ideal for light bowel or urinary incontinence, moisture management, skin-fold, ostomy care and more! Use in regular underwear or as a booster pad - diaper doubler.



25 Sheets per pack

## Product Description:

Soft, supple absorbent sheets, used for the collection and retention of skin-fold moisture. ThinLiners have TQ-A/P (Tranquility's Antimicrobial and Antifungal Protection) treatment bonded to both sides of the sheet. The treatment inhibits the growth of microorganisms while preserving an ideal environment on the product for skin integrity. Ideal for tucking, folding, wrapping in or around skin folds, groin area, under breasts or where light-to-moderate moisture is present. Ultra-thin sheets wick moisture away from the skin and lock it into the absorbent core. This action creates a dry environment that can prevent significant skin breakdown. Latex-free. Non-sterile. These cost-effective, multipurpose absorbent sheets can also be used for urinary and fecal incontinence management. Latex free. Non-sterile.

**Size:** 6 x10 inches

**Capacity:** 5 ounces

## Product Features and Descriptions:



### **BR – Breathable**

Allows air to circulate around the body, promoting skin integrity.

### **LF – Latex-Free**