

3.5 Gram Nitrile Exam Gloves, Blue

(Tested for Use with Chemotherapy Drugs)

Features:

- **Colors:** Available in Ocean Blue and Violet Blue (when ocean blue stock has been depleted, violet blue gloves will be shipped instead)
 - Textured at fingertip surface
 - Larger quantity per box to help save time and reduce waste
 - Non-sterile
 - Ambidextrous
 - Powder-free exam gloves
 - ASTM D6319 compliant
 - Manufactured in compliance with ISO 9001, ISO 13485
 - Not made with natural rubber latex
 - Chemotherapy drugs tested: 5-Fluorouracil, Carmustine (BCNU), Cisplatinum, Cyclophosphamide (Cytosan), Dacarbazine (DTIC), Doxorubicin Hydrochloride, Etoposide (Toposar), Metotrexate, Paclitaxel (Taxol), Thiotepa, and Vincristine Sulfate
 - **County of Origin:** Thailand
- **Storage:** Keep storage area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Do not store above 100°F (38°C).



Ocean Blue



Violet Blue



| MFR # | E1 Item # | Description | Size | Quantity |
|---------|-----------|--------------------------------|---------|-----------------|
| 7025428 | 1175360 | 3.5g Nitrile Exam Gloves, Blue | X-Small | 250/bx, 10bx/cs |
| 7025429 | 1175361 | 3.5g Nitrile Exam Gloves, Blue | Small | 250/bx, 10bx/cs |
| 7025430 | 1175362 | 3.5g Nitrile Exam Gloves, Blue | Medium | 250/bx, 10bx/cs |
| 7025431 | 1175363 | 3.5g Nitrile Exam Gloves, Blue | Large | 250/bx, 10bx/cs |
| 7025432 | 1175364 | 3.5g Nitrile Exam Gloves, Blue | X-Large | 225/bx, 10bx/cs |



DO NOT REUSE

McKesson Medical-Surgical Inc.
 Richmond, VA 23233
 Questions? Call 1-800-777-4908
mms.mckesson.com

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information. The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.