

## Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters

HCPCS: E0143



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F

- Comes standard with basket, tray and pouch (Figure A)
- Basket conveniently folds for easy storage (Figure B)
- Ships in ultra-compact carton Comes with frame folded and compact Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob
- 7.5" casters for indoor or outdoor use
- Handles easily adjust in height with self-threading knob (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure D)
- Easy, one-hand folding

PRODUCTS	Description
10289BL	3 Whl Alum Rollator,Blue,1/ea
10289RD	3 Whl Alum Rollator,Red,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Handle Height	31" - 38"
Weight Capacity	300 lb
Carton Shipping Weight	16 lb
Weight	14 lb

PARTS	Description
9503E17101-L	Hand Brake with Cable, 10289, left 1/ea, EA
9503E17101-R	Hand Brake with Cable, 10289, Rt, 1/ea, EA
816T	Tray, 10289, 1/ea, EA
815B	Folding Basket, RTL199/10289, 1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
9503E17106	Rear Axle with Nut,10289, 1/ea, EA
9503E17108	Wheel Assembly,10289, 1/ea, EA
9503E17112	Folding Hinge,10289, 1/ea, EA
9503E17103	Height Adj.Knobs, 10289/RTL199, 1/pr, PR
9503E17104	Wheel Brake,Left,10289, 1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
950330049901	Hand Grip, Right, 3-wheel Rollator, 1/ea, EA
950330039901	Hand Grip, Left, 3-wheel Rollator, 1/ea, EA