


Preserve

Safeguard granulating tissue



CUR250383



CUR253180

CURAD® Oil Emulsion Dressings

Compare to **Adaptic® non-adhering dressing***

- Knitted, high porosity cellulose acetate is evenly impregnated with U.S.P. white petrolatum
- Emulsion blend permits the flow of exudates without adhering to granulating tissue
- Minimizes skin maceration; not made with natural rubber latex
- Use directly on wounds beneath an absorbent secondary dressing
- Ideal for most draining wounds including: minor burns, suture lines, lacerations, skin ulcers, episiotomy incisions, abrasions

Item No.	Description	Pkg.
CUR250330	3" x 3" (7.6 x 7.62 cm)	50/bx, 600/cs
CUR250381	3" x 8" (7.6 x 20.3 cm)	1/pk, 24/bx, 144/cs
CUR250383	3" x 8" (7.6 x 20.3 cm)	3/pk, 108/bx, 648/cs
CUR250317	3" x 16" (7.6 x 40.6 cm)	36/bx, 216/cs

*This product is not manufactured or distributed by Johnson & Johnson Corporation, distributors of Adaptic® non-adherent dressings.

CURAD® Xeroform Gauze Dressings

- Petrolatum gauze with the addition of 3% Xeroform (Bismuth Tribromophenate), a deodorizing agent
- Contraindicated for heavily draining wounds
- Occlusive dressing; not made with natural rubber latex
- Ideal for non-draining to lightly draining wounds including: donor sites, lacerations, minor or partial thickness burns, abrasions, skin graft sites, suture wounds

Item No.	Description	Pkg.
Overwrap Gauze Dressing, Sterile		
CUR253590W	5" x 9" (12.7 x 22.9 cm)	50/bx, 200/cs
CUR253180W	1" x 8" (2.5 x 20.3 cm)	50/bx, 200/cs
Gauze Dressing, Sterile		
CUR253180	1" x 8" (2.5 x 20.3 cm)	50/bx, 200/cs
CUR253590	5" x 9" (12.7 x 22.9 cm)	50/bx, 200/cs
CUR253220	2" x 2" (5.1 x 5.1 cm)	25/bx, 150/cs
CUR253440	4" x 4" (10.2 x 10.2 cm)	25/bx, 150/cs
CUR253200	4" x 3 yd (10.2 cm x 2.7 m)	6/bx, 36/cs