

# Enteral Feeding Pumps and Pump Sets

## Kangaroo™ ePump, Joey and Connect



382400

### Kangaroo™ ePump Enteral Feeding Pump

Features feed-and-flush technology that automates flushing for tube patency and delivery of hydration.

- Fully programmable variable feed and flush pump
- Flow rate: 1 to 400 mL/hr in 1 mL increments
- Flush rate: 1960 mL/hr
- Flush volume: 10 to 500 mL/hr in 1 mL increments
- Flush interval: 1 to 24 hr in 1 hr increments
- Multilingual: 15 languages
- Battery life: 15 hours
- Weight without pole clamp: 2.4 lbs
- Weight with pole clamp: 2.7 lbs
- Feeding history up to 72 hours
- Attitude-independent
- Intuitive user interface
- Auto-prime



773656

Cat. no.	Description	Ship case
<b>Pump</b>		
<b>382400</b>	Kangaroo™ ePump with pole clamp	1
<b>Pump accessories</b>		
<b>382491</b>	Kangaroo™ ePump Power Adapter	1
<b>382492</b>	Kangaroo™ ePump Pole Clamp	1
<b>1054715</b>	Kangaroo™ ePump Battery	1



773662

### Kangaroo™ ePump Enteral Feeding Sets

- Anti-free flow design
- DEHP-free
- ENFit™ Connection

Cat. no.	Description	Ship case
<b>Pump sets</b>		
<b>772055</b>	Kangaroo™ ePump 500 mL Pump Set	30
<b>773656</b>	Kangaroo™ ePump 1000 mL Pump Set	30
<b>674668</b>	Kangaroo™ ePump 500 mL Feed Set with 500 mL flush bag	30
<b>773662</b>	Kangaroo™ ePump 1000 mL Feed Set with 1000 mL flush bag	30
<b>716154</b>	Kangaroo™ ePump Sterile 100 mL Burette Set	30
<b>775659</b>	Kangaroo™ ePump ENPlus Spike Set	30
<b>775100</b>	Kangaroo™ ePump ENPlus Spike with 1000 mL flush bag	30
<b>776150</b>	Kangaroo™ ePump 100 mL Burette Recertification Set (not ENFit™)	5



716154

