

## ORDERING INFORMATION

## Part Number Description

## ColPac Chilling Unit Model C-5

3102	C-5 ColPac Chilling Unit Includes 6 Standard (1500) and 6 Half Size (1506) ColPacs
------	--

## Standard Accessories

4010	Side Table Rack <i>See page 4 for Description and image</i>
4230	Extra Shelf for Side Table Rack ( <i>not shown</i> )

## ColPac Freezer

90910	ColPac Freezer 21 1/4" x 26" x 33 1/2" (53 cm x 66 cm x 86 cm) Includes 12 Standard ColPacs (1500) Shipping weight 150 lbs (68 kg)
-------	---

## TECHNICAL SPECIFICATIONS

## ColPac Chilling Unit Model C-5

Mains Power:	115 VAC ~ 60 Hz
Current:	4A
Weight:	165 lb (74 kg) Shipping Weight - 190 lb (86 kg)
Dimensions*:	27" x 16" x 33" (69 x 41 x 84 cm)
Electrical Safety Class:	Type B
Safety Tests:	UL 544 (U.S. Only)

\* Dimensions are listed as L x W x H

## COLPAC® CHILLING UNIT



Large capacity and fast chilling time means your patients will never have to wait for a ColPac. The ColPac Chilling Units help you keep pace throughout a busy day. No plumbing required, simply plug into a standard electrical outlet.



## FEATURES

## ColPac Chilling Unit

- Closed-cell foam insulation for efficiency
- Stainless steel, welded construction for long life
- Heavy-duty compressor for more efficient cooling
- Specially designed interior chilling coils cool packs faster
- 3-inch (8 cm) swivel-type rubber casters for silent, friction-free movement of mobile units
- Drain valve for easy cleaning and defrosting
- One-year warranty

## ColPac Freezer

- Five cubic feet of spacious capacity
- Four durable stainless steel shelves
- Comes complete with trays and cup holders
- Shipped with 12 standard size ColPacs
- UL listed and cUL listed
- One-year warranty



## COLPAC®, BOO-BOO PAC™, FLEXIPAC®, SPORT-PAC™

Seven different sizes of blue vinyl and four sizes of black polyurethane ensure effective cold therapy for any body area or shape.



### FEATURES

#### The Original ColPac – There's a ColPac® for Every Application!

- Designed for heavy clinical use.
- Effective pain relief
- Filled with a non-toxic silica gel that remains pliable throughout treatment
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Latex free
- One-year warranty

#### boo-boo pac

- Fun, bear-shaped design especially made for kids
- Ideal for pediatric applications
- Filled with non-toxic silica gel
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Shipped in a display box containing 12 packs

#### FlexiPac

- Freezable to provide penetrating cold therapy
- Microwavable for soothing heat
- Flexible gel conforms to any body contour at any temperature

#### Sport-Pac

- Soccer ball design is ideal for youth sport applications
- Easily transferred from the freezer to the team cooler for sideline use
- Delivers up to 30 minutes of effective cold therapy for acute pain due to strains, sprains and bruises
- Shipped in a display box containing 10 packs

### ORDERING INFORMATION

Part Number	Description
<b>Blue Vinyl ColPac</b>	
A 1512	Override 11" x 21" (28 cm x 53 cm)
B 1500	Standard Size 11" x 14" (28 cm x 36 cm)
C 1506	Half Size 7½" x 11" (19 cm x 28 cm)
D 1508	Neck Contour 23" long (58 cm)
E 1502	3" x 11" (8 cm x 28 cm)
F 1504	Quarter Size 5½" x 7½" (14 cm x 19 cm)
G 1510	Eye (Specially shaped to cover closed eyelids)
<b>Black Polyurethane ColPac</b>	
H 1554	Neck Contour 21" long (53 cm)
I 1556	Override 12½" x 18½" (32 cm x 47 cm)
J 1552	Standard Size 10" x 13½" (25 cm x 34 cm)
K 1562	Half Size 6½" x 11" (17 cm x 28 cm)
<b>L boo-boo pac</b>	
Shipped in a display box containing 12 packs. NOTE: Must be purchased in quantities of 12 each	
1534	Royal Blue boo-boo pac (12 per case)
<b>M FlexiPac Hot &amp; Cold Compresses</b>	
4026	5" x 6" (13 cm x 15 cm) 48/case
4020	5" x 10" (13 cm x 25 cm) 24/case
4029	8" x 14" (20 cm x 35 cm) 12/case
4032	5" x 10" (13 cm x 25 cm) 25/case Chiropractically Yours (not shown)
<b>N Sport-Pac</b>	
Shipped in a display box containing 10 packs. NOTE: Must be purchased in quantities of 10 each	
1517	Sport-Pac (10 per case)