

ComfortGlide® TRP Patient Repositioning System

A safer, easier way to move patients



**Take the
pain out of
repositioning**

How are staff injuries impacting you?

Caregivers suffer some of the highest rates of musculoskeletal injuries from overexertion resulting in missed workdays, decreased retention and billions of dollars in direct and indirect cost.^{1,2,3}

Studies show that patient repositioning using a low-friction glide sheet may help reduce the risk of musculoskeletal strain and injuries to caregivers.⁴

Do you know how much staff injuries are costing you?

\$15,600

Average worker's compensation claim related to patient handling⁵

\$20B

Estimated direct and indirect costs associated with only back injuries in the healthcare industry⁶

**We surveyed 149 nurses
Here's how they responded.⁷**

36%

have thought about leaving the profession due to pain or back injury

75%

of respondents have strained themselves while repositioning a patient

ComfortGlide is your solution

ComfortGlide provides caregivers the important tools to safely boost, turn, reposition and transfer patients based on your facility's safe patient handling requirements.

Safe patient handling is equally as important to patients. Bed-ridden patients are at risk for pressure injuries, which can develop in as little as one hour.²² Frequent repositioning using transfer aids to reduce friction and shear is recommended to help prevent the development of pressure injuries.⁷

Caregiver Benefits


- Provides the proper tools for caregivers to follow facility protocols

Patient Benefits

- Wedges provide the NPIAP recommended angle for offloading
- Moisture management pad keeps the patient's skin dry
- Allows for better patient experience and outcomes



ComfortGlide Air Pg. 4



ComfortGlide Sling Pg. 6



ComfortGlide Repositioning Sheets Pg. 8



ComfortGlide Wedges Pg. 10



ComfortGlide Drypads Pg. 10



ComfortGlide Fitted Sheets Pg. 10



Moving patients is physically demanding

Caregivers suffer back and other overexertion injuries severe enough that they have to miss work.^{8,9,10} The number one reason why: manually lifting, transferring, boosting and repositioning patients multiple times a day.⁸

ComfortGlide Air

Inflatable sheet gives caregivers a safe and physically easy way to move patients.

- Always ready to use (remains under patient)
- Significantly reduces strain of boosting for staff
- Compatible with low air-loss mattresses
- Patient remains below bed rail when inflated
- Safety rated for patients up to 1,000 lbs.²¹





Scan the QR code
to watch the
ComfortGlide Sling
instructional video

Move patients more safely and efficiently


Lifting patients manually can cause staff injuries. 20% of nurses leave direct care positions as a result of low back pain,¹ creating high costs for facilities.

ComfortGlide Sling

2-in-1 repositioning sheet and patient sling

Every time your staff lifts a patient to place a sling, there are back injury risks. Our ComfortGlide Sling can be left under a patient throughout their stay and used as both a glide sheet and a sling.

- **Reduce manual patient handling.** Once the glide sheet is in place, there's no need to move the patient to use a lift.
- **Save staff time.** Fewer steps are needed to lift and move patients.
- **Optimize care.** Use with ComfortGlide system of wedges and dry pads.
- Low-friction nylon requires less exertion than a typical draw sheet
- Multiple handles provide leverage for several caregivers
- Compatible with low air-loss mattresses
- 1,000-lb. lifting capacity and large size accommodates the majority of your patients
- Sturdy loops for easy attachment to a lift



Use as a sling



Use as a glide sheet





78"
X-Large
Length

55"
Standard
Length

Reduce the effort needed to turn patients

Intuitive and easy to use, ComfortGlide Repositioning Sheets are safety rated to accommodate virtually all patients—575 lb. weight capacity for standard, 1,000 lb. weight capacity for X-Large.¹¹

ComfortGlide Repositioning Sheets

These sheets are designed and tested to help reduce the exertion needed to turn patients.


- Low-friction nylon bottom surface is shown to require 35 percent less effort than a standard draw sheet to initiate and maintain sliding¹²
- Multiple handles help optimize weight distribution and leverage
- Air-permeable materials allow sheet to remain under patients for an extended period of time and are always available to aid repositioning¹³
- Soft, breathable microfiber covers the entire top surface for patient comfort and is compatible with low air-loss mattresses¹³
- Optional retention straps anchor the sheet to the bed, helping keep the patient in the desired position



The Evidence Base

Average force required to initiate slide


This graph shows the force required to initiate the slide a 200-pound object across a hospital bed.¹⁴ The lower the number, the easier it is for healthcare workers to initiate sliding a patient.



lbf = pounds of force.

WVTR test for breathability

Water Vapor Transmission Rate (WVTR in units $g/(m^2 \cdot day)$) measures the amount of water vapor that passes through a material. The higher the number, the more breathable the material. The ComfortGlide sheet performed as well as the leading competitor.¹⁵



g = grams. m^2 = square meters, day = days



Scan the QR code to watch the ComfortGlide Repositioning Sheet instructional video



ComfortGlide Wedges

High-density foam wedges help keep patients in the 30-degree tilt position recommended by NPIAP¹⁶, helping to reduce pressure on bony prominences.

- Non-slip base grips the surface of any mattress or linen
- “Down” arrows ensure proper placement
- Soft, composite fabric covering can be wiped clean
- We offer a variety of wedges including a standard, XL, XXL and even a 15 degree wedge that can be used for microturns



Scan the QR code to watch the ComfortGlide Wedges instructional video

ComfortGlide Drypad


Ultra-absorbent drypads are designed to manage moisture over an extended period of time and keep patients 26 times drier than the leading competitor.¹⁷

- Draws in moisture, locks it away from skin and feels dry to the touch within minutes¹⁷
- Air-permeable for optimal patient comfort and zero strikethrough^{17,18}
- Compatible with low air-loss mattresses¹⁷
- Both standard (36" x 57") and XL size (40" x 80") available; each size designed to completely cover the repositioning sheets
- Testing available upon request

ComfortGlide Fitted Sheets

Fitted sheets feature a low-friction,¹⁹ breathable¹⁸ polyester fabric that helps reduce the effort to slide a patient.

- Contoured elastic perimeter hugs the mattress for a close fit without awkward buckles
- Compatible with low air-loss mattresses¹⁸
- Can be laundered up to 50 times²⁰
- Available in four different sizes to fit a wide range of mattresses



Customized ordering

Purchase only what you need

Item No.	Description	Pkg
ComfortGlide Air, Single-Patient-Use		
MSC60011	ComfortGlide Air Inflatable Sheet, 1,000-lb. weight capacity	12/cs
MSC60012	ComfortGlide Air Inflatable Sheet w/ Drypad, 1,000 lb. weight capacity	12/cs
MSC60090	ComfortGlide Air Blower	1/ea
MSC60099	ComfortGlide Air Blower Mobile Cart with Basket	1/ea
ComfortGlide Sling, Single-Patient-Use		
MDSC40018	2-in-1 Sling and Repositioning Sheet, 1,000-lb. weight capacity, 42" x 78"	12/cs
MDSC40017	2-in-1 Sling and Repositioning Sheet, 1,000-lb. weight capacity, 42" x 65"	12/cs
ComfortGlide Sheets, Single-Patient-Use		
CGLIDESHEET	Glide Sheet, 35" x 55", with straps, 575-lb. weight capacity	12/cs
CGLIDESHEETNS	Glide Sheet, 35" x 55", 575-lb. weight capacity	12/cs
CGSHEETXLD	Glide Sheet, X-Large, 42" x 78", 2 optional straps, 1,000-lb. weight capacity	12/cs
MSC50013	Glide Sheet, 35" x 55", Tacky Top, 450-lb weight capacity	12/cs
CGLIDESHEETM	Glide Sheet, Mini, 30" x 35", 500-lb weight capacity	12/cs
ComfortGlide Foam Wedges, Single-Patient-Use		
MSC50023	30-Degree Foam Wedges, non-slip base: 11" x 5" x 16"	2/pr, 6 pr/cs
MSC50024	30-Degree Foam Wedges, non-slip base: 15" x 6.½" x 16"	2/pr, 2 pr/cs
MSC50023A	30-Degree Foam Wedges, non-slip base: 11" x 5" x 16", Made in USA	2/pr, 6 pr/cs
MSC50024A	30-Degree Foam Wedges, non-slip base: 15" x 6.5" x 16", Made in USA	2/pr, 2 pr/cs
MSC50025A	30-Degree Foam Wedges, non-slip base: 24" x 12" x 7", Made in USA	2 ea/cs
MSC50070	15-Degree Foam Wedges, non-slip base: 16" x 11" x 3", Made in USA	2/pr, 6 pr/cs
ComfortGlide Drypads, Single Use Only		
CGLIDEPAD	Drypads, 36" x 57"	5/bg, 60/cs
CGLIDEPADXL	X-Large Drypads, 40" x 80"	5/bg, 40/cs
ComfortGlide Fitted Base Sheets, Reusable to 50 Launderings²⁰		
CGLIDEBASE36	Fitted Base Sheet, 36" x 80", Blue	12/cs
CGLIDEBASE42	Fitted Base Sheet, 42" x 80", Purple	12/cs
CGLIDEBASE50	Fitted Base Sheet, 50" x 80", Black	12/cs
CGLIDEBASE54	Fitted Base Sheet, 54" x 80", Black	12/cs



We're here to support you



Learn more. Contact your Medline Representative today or call 1-800-MEDLINE.



We make healthcare run better™

Medline Industries, LP

Three Lakes Drive, Northfield, IL 60093
 Medline United States | 1-800-MEDLINE (633-5463)
 medline.com | info@medline.com

Medline Canada

1-800-268-2848 | medline.ca | canada@medline.com

Medline México

01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

References: **1.** de Castro AB. Handle With Care®: The American Nurses Association's Campaign to Address Work-Related Musculoskeletal Disorders. The Online Journal of Issues in Nursing. 2004;9(3). Available at <http://www.nursingworld.org/MainMenuCategories/ANAMarketplace/ANAPeriodicals/OJIN/TableofContents/Volume92004/No3Sept04/HandleWithCare.html> Accessed October 2, 2015. **2.** Safe Patient Handling. Occupational Safety and Health Administration (OSHA). Available at <https://www.osha.gov/SLTC/healthcarefacilities/safepatienthandling.html> Accessed October 2, 2015. **3.** Safe Patient Handling -- Preventing Musculoskeletal Disorders in Nursing Homes. Occupational Safety and Health Administration (OSHA). Available at <https://www.osha.gov/Publications/OSHA3708.pdf> Accessed October 2, 2015. **4.** Theou O, Soon Z, Filek S, et al. Changing the Sheets: A New System to Reduce Strain During Patient Repositioning. Nurs Res. 2011; 60(5):302-8. Available at <http://www.ncbi.nlm.nih.gov/pubmed/21873921#> Accessed October 2, 2015. **5.** Safe Patient Handling Programs: Effectiveness and Cost Savings. Safety and Health Topics. Occupational Safety and Health Administration (OSHA). Available at <https://www.osha.gov/Publications/OSHA3279.pdf> Accessed on October 2, 2015. **6.** Shekelle PG, Wachter RM, Pronovost PJ, et al. Making Health Care Safer II: An Updated Critical Analysis of the Evidence for Patient Safety Practices. Comparative Effectiveness Review No. 211. (Prepared by the Southern California-RAND Evidence-based Practice Center.) AHRQ Publication No. 13-E001-EF. Rockville, MD: Agency for Healthcare Research and Quality. March 2013. Available at www.ahrq.gov/research/findings/evidence-based-reports/ptsafetyupt.html Accessed October 2, 2015. **8.** Safe Patient Handling and Mobility (SPHM). Centers for Disease Control and Prevention website. Available at: <https://www.cdc.gov/niosh/topics/safepatient/#source1> Accessed May 16, 2017. **9.** Hospitals fail to protect nursing staff from becoming patients—part 1 [transcript]. All Things Considered. National Public Radio. February 4, 2015. Available at: <http://www.npr.org/2015/02/04/382639199/hospitals-fail-to-protect-nursing-staff-from-becoming-patients> Accessed May 16, 2017. **10.** Bureau of Labor Statistics. Case and Demographic Incidence Rates. In this figure, "hospitals" represents NAICS 622, which covers all types of hospitals. "Construction" represents NAICS supersector GP1CON, "manufacturing" represents NAICS supersector GP1MFG, and "professional and business services" represents NAICS supersector SP1PBS. Data are limited to private industry. **11.** Safe Patient Handling. Occupational Safety and Health Administration (OSHA). Available at: www.osha.gov/SLTC/healthcarefacilities/safepatienthandling.html Accessed April 26, 2017. **7,12-21.** Data on file. Available upon request. **22.** OSHA. Occupational Safety and Health Administration. Safety and Health Topics. Healthcare. Safe Patient Handling. Available at: <https://www.osha.gov/SLTC/healthcarefacilities/safepatienthandling.html>. Accessed June 17, 2019.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

©2021 Medline Industries, LP. All rights reserved. ComfortGlide and Medline are registered trademarks of Medline Industries, LP. MKT19W2157002 / LIT424R / 2.5M / SG / 3