

Folding Steel Commode

HCPCS: E0163



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Easily opens and folds
- Folds flat for convenient storage and transportation
- Tool-free set-up (Figure A - F)
- Durable plastic snap on seat and lid
- Grey powder coated steel welded construction increases strength and durability
- Plastic arms for added comfort
- Comes complete with 7.5 qt commode bucket, cover and splash shield
- Limited Lifetime Warranty
- Piggy back label supplied on product hangtag
- Made using recycled materials

PRODUCTS

Description

11148-1	Commode, Folding, 1/cs
---------	------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	63 lbs
Seat Dimensions	13.5" (W) x 15" (D)
Seat Height	16.6" x 22.5" (H)
Warranty	Limited Lifetime
Weight	15 lbs
Weight Capacity	350 lbs
Width	18" between arms 22.5" outside