

## Deluxe, Trigger Release Folding Walker with 5" Wheels

HCPCS: E0143



FIGURE A



FIGURE B




FIGURE C

- Durable, Composite Trigger Release resists cracking and breaking
- Sturdy aluminum construction ensures maximum strength while remaining lightweight
- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity
- Durable, composite trigger release can be reached with hands of all sizes and resists cracking and breaking, ensuring safety and durability
- Ribbed, contoured handgrips relieve pressure on users hands
- Using the marked peg holes, easily adjust the height in 1" increments for a customized fit
- 5" wheels are ideal for maneuvering both indoors and outdoors
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Folds for easy transport and the sides are detachable for storing in tight spaces
- Quick and tool-free assembly

### PRODUCTS

#### Description

10226-1	Walker, Trigger Rls, Adlt, w/5"Whl, 1/ea
10227-4	Walker, 2 Trigger Rel. Junr. w/5"Whls, 4/cs

### SPECIFICATIONS

#### Description

Carton Shipping Weight	10226-1: 29 lbs   10227-4: 26 lbs
Depth Opened at Base	10226-1: 21.5"   10227-4: 17"
Handle Height Range	10226-1: 32.5" - 39.5"   10227-4: 25.5" - 32.5"
Weight	10226-1: 7 lbs   10227-4: 6 lbs
Weight Capacity	10226-1 and 10227-4: 350 lbs
Wheel Size	10226-1 and 10227-4: 5"
Width inside back legs	10226-1: 27"   10227-4: 24"

### ACCESSORIES

#### Description

10200B	Walker Basket, 1/ea, RTL, EA
RTL10131	Walker Caddy, 1EA, 3/CS, CS
10121	Tennis Ball Glides, 1pr./bx, RTL, BX
10110	Glide Ski for Walker, 1pr./bx, RTL, BX
RTL10254BFL	Univ Mobility Bag, Med, Blue Flrl, 1EA, EA
RTL10254GR	Univ Mobility Bag, Med, Army Green, 1EA, EA
RTL10254TFL	Univ Mobility Bag, Med, Trop Flrl, 1EA, EA
RTL10255NV	Univ Mobility Bag, Small, Navy, 1EA, EA

**ACCESSORIES****Description**

RTL10255PK	Univ Mobility Bag, Small, Pink, 1EA, EA
RTL10256NV	Univ MobilityBag, Med Flat, Navy, 1EA, EA
RTL10256PK	Univ MobilityBag, Med Flat, Pink, 1EA, EA
RTL10256TL	Univ MobilityBag, Med Flat, Teal, 1EA, EA