

P903

Item Number	Item Description	UOM
Standard Hi/Low Beds		
P903	PrimeCare 903 Full Elect Hi/Low Bed with Integrated 42" width expansion and one-step central brake lock 600 lb. weight capacity. 7" low height - 30" high frame height.	EA
Assist Device		
PLTCAB	Rotating Assist Bar for P503/P703/P903 Beds (sold as pairs)	PR
PLTCAR	Rotating Assist Rail for P503/P703/P903 Beds (sold as pairs)	PR
PLTCDR	DELUXE "soft" Rotating Assist Rail for P703/P903 Beds (sold as pairs)	PR
PLTCABR	RIGHT side, Rotating Assist Bar for P503/P703/P903 Beds	EA
PLTCABL	LEFT side, Rotating Assist Bar for P503/P703/P903 Beds	EA
PLTCARR	RIGHT side, Rotating Assist Rail for P503/P703/P903 Beds	EA
PLTCARL	LEFT side, Rotating Assist Rail for P503/P703/P903 Beds	EA
PLTCDRR	RIGHT side, DELUXE "soft" Rotating Assist Rail for P703/P903 Beds	EA
PLTCDRL	LEFT side, DELUXE "soft" Rotating Assist Rail for P703/P903 Beds	EA
Embedded Staff Control		
P903SC03	Embedded Staff Control for P903 bed, panel Only.	EA
35" Bed Ends		
35SHBE-A	35" Smokey Hollow bed end in Oslo Cherry	EA
35SHBE-B	35" Smokey Hollow bed end in Oslo Oak	EA
35SHBE-FM	35" Smokey Hollow bed end in Oslo Figured Mahogany	EA
35MCBE-A	35" Mill Creek bed end in Oslo Cherry	EA
35MCBE-B	35" Mill Creek bed end in Oslo Oak	EA
35MCBE-FM	35" Mill Creek bed end in Oslo Figured Mahogany	EA
35BWBE-A	35" Brownlee Woods bed end in Oslo Cherry	EA
35BWBE-B	35" Brownlee Woods bed end in Oslo Oak	EA
35BWBE-FM	35" Brownlee Woods bed end in Oslo Figured Mahogany	EA
35" Bed Ends with Embedded Staff Control		
35SHBE-A-SC03	35" Smokey Hollow bed end in Oslo Cherry w/ Embedded Staff Control	EA
35SHBE-B-SC03	35" Smokey Hollow bed end in Oslo Oak w/ Embedded Staff Control	EA
35SHBE-FM-SC03	35" Smokey Hollow bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
35MCBE-A-SC03	35" Mill Creek bed end in Oslo Cherry w/ Embedded Staff Control	EA
35MCBE-B-SC03	35" Mill Creek bed end in Oslo Oak w/ Embedded Staff Control	EA
35MCBE-FM-SC03	35" Mill Creek bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
35BWBE-A-SC03	35" Brownlee Woods bed end in Oslo Cherry w/ Embedded Staff Control	EA
35BWBE-B-SC03	35" Brownlee Woods bed end in Oslo Oak w/ Embedded Staff Control	EA
35BWBE-FM-SC03	35" Brownlee Woods bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
42" Bed Ends		
42SHBE-A	42" Smokey Hollow bed end in Oslo Cherry	EA
42SHBE-B	42" Smokey Hollow bed end in Oslo Oak	EA
42SHBE-FM	42" Smokey Hollow bed end in Oslo Figured Mahogany	EA
42MCBE-A	42" Mill Creek bed end in Oslo Cherry	EA
42MCBE-B	42" Mill Creek bed end in Oslo Oak	EA
42MCBE-FM	42" Mill Creek bed end in Oslo Figured Mahogany	EA
42BWBE-A	42" Brownlee Woods bed end in Oslo Cherry	EA
42BWBE-B	42" Brownlee Woods bed end in Oslo Oak	EA
42BWBE-FM	42" Brownlee Woods bed end in Oslo Figured Mahogany	EA
42" Bed Ends with Embedded Staff Control		
42SHBE-A-SC03	42" Smokey Hollow bed end in Oslo Cherry w/ Embedded Staff Control	EA
42SHBE-B-SC03	42" Smokey Hollow bed end in Oslo Oak w/ Embedded Staff Control	EA
42SHBE-FM-SC03	42" Smokey Hollow bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
42MCBE-A-SC03	42" Mill Creek bed end in Oslo Cherry w/ Embedded Staff Control	EA
42MCBE-B-SC03	42" Mill Creek bed end in Oslo Oak w/ Embedded Staff Control	EA
42MCBE-FM-SC03	42" Mill Creek bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
42BWBE-A-SC03	42" Brownlee Woods bed end in Oslo Cherry w/ Embedded Staff Control	EA
42BWBE-B-SC03	42" Brownlee Woods bed end in Oslo Oak w/ Embedded Staff Control	EA
42BWBE-FM-SC03	42" Brownlee Woods bed end in Oslo Figured Mahogany w/ Embedded Staff Control	EA
CUSTOM Bed Ends		
*** Drop Ship Only from WI, 4-6 week lead time***		
35SHBE-CUSTOM	Custom 35" Smokey Hollow bed end	EA
35MCBE-CUSTOM	Custom 35" Mill Creek bed end	EA


35BWBE-CUSTOM	Custom 35" Brownlee Woods bed end	EA
35BHBE-CUSTOM	Custom 35" Brier Hill bed end	EA
35BLHB-ARCH	Custom 35" Bronze Lake Head Board, ARCH, bed end	EA
35BLFB-ARCH	Custom 35" Bronze Lake Foot Board, ARCH, bed end	EA
35BLHB-CROWN	Custom 35" Bronze Lake Head Board, CROWN, bed end	EA
35BLFB-CROWN	Custom 35" Bronze Lake Foot Board, CROWN, bed end	EA
35SHBE-CUSTOM-SC03	Custom 35" Smokey Hollow bed end w/ Embedded Staff Control	EA
35MCBE-CUSTOM-SC03	Custom 35" Mill Creek bed end w/ Embedded Staff Control	EA
35BWBE-CUSTOM-SC03	Custom 35" Brownlee Woods bed end w/ Embedded Staff Control	EA
35BHBE-CUSTOM-SC03	Custom 35" Brier Hill bed end w/ Embedded Staff Control	EA
35BLHB-ARCH-SC03	Custom 35" Bronze Lake Head Board, ARCH, bed end w/ Embedded Staff Control	EA
35BLFB-ARCH-SC03	Custom 35" Bronze Lake Foot Board, ARCH, bed end w/ Embedded Staff Control	EA
35BLHB-CROWN-SC03	Custom 35" Bronze Lake Head Board, CROWN, bed end w/ Embedded Staff Control	EA
35BLFB-CROWN-SC03	Custom 35" Bronze Lake Foot Board, CROWN, bed end w/ Embedded Staff Control	EA
42SHBE-CUSTOM	Custom 42" Smokey Hollow bed end	EA
42MCBE-CUSTOM	Custom 42" Mill Creek bed end	EA
42BWBE-CUSTOM	Custom 42" Brownlee Woods bed end	EA
42BHBE-CUSTOM	Custom 42" Brier Hill bed end	EA
42BLHB-ARCH	Custom 42" Bronze Lake Head Board, ARCH, bed end	EA
42BLFB-ARCH	Custom 42" Bronze Lake Foot Board, ARCH, bed end	EA
42BLHB-CROWN	Custom 42" Bronze Lake Head Board, CROWN, bed end	EA
42BLFB-CROWN	Custom 42" Bronze Lake Foot Board, CROWN, bed end	EA
42SHBE-CUSTOM-SC03	Custom 42" Smokey Hollow bed end w/ Embedded Staff Control	EA
42MCBE-CUSTOM-SC03	Custom 42" Mill Creek bed end w/ Embedded Staff Control	EA
42BWBE-CUSTOM-SC03	Custom 42" Brownlee Woods bed end w/ Embedded Staff Control	EA
42BHBE-CUSTOM-SC03	Custom 42" Brier Hill bed end w/ Embedded Staff Control	EA
42BLHB-ARCH-SC03	Custom 42" Bronze Lake Head Board, ARCH, bed end w/ Embedded Staff Control	EA
42BLFB-ARCH-SC03	Custom 42" Bronze Lake Foot Board, ARCH, bed end w/ Embedded Staff Control	EA
42BLHB-CROWN-SC03	Custom 42" Bronze Lake Head Board, CROWN, bed end w/ Embedded Staff Control	EA
42BLFB-CROWN-SC03	Custom 42" Bronze Lake Foot Board, CROWN, bed end w/ Embedded Staff Control	EA
Trapeze		
P903TPZ	Trapeze for P903	EA



Model #PLTCAB



Model #PLTCDR



Model #PLTCAR