



ENCORE<sup>®</sup>

Powder-Free Latex Surgical Gloves

**Ansell**

Everything you touch...we touch.

# Encore® Powder-Free Latex Surgical Gloves

At Ansell, we understand how much responsibility rests in your hands. That's why we work hard to develop innovative surgical gloves like our powder-free and latex Encore® line. The Encore® family includes a glove for every need and provides the superior barrier protection of latex while decreasing the risk associated with powder. All our Encore® family gloves provide the following:

## Safety and Protection

### Infection Control

**Viral Penetration.** Although not required by the FDA, Ansell's premium line of surgical gloves are tested to ensure protection against viral permeation. The Encore® family of gloves has passed the viral penetration requirements according to ASTM F1671.

**Cuff Roll-Down.** Ansell's SureFit Technology™ overcomes the "cuff roll-down" typically experienced while wearing impervious surgical gowns designed to prevent fluid penetration. By maintaining the surgical glove securely around the arm, SureFit Technology™ protects the integrity of the sterile field by shielding the gown cuff and ensures the efficacy of the surgical glove as a barrier protection.

### Skin Protection

**Low Protein Content.** All Encore® gloves contain 50µg/g or less of water/extractable protein per gram, to help eliminate powder-related post-operative complications and sensitivities.



## Comfort & Fit

### Donning

**Alcohol Resistant.** Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants and adds another layer of safety between your skin and the latex.

**Damp-Hand Donnability.** Polymer coating for dry and damp-hand donnability.

### Hand Fatigue

**Shape of Glove/Former.** Designed with Anatomic Curved Formers for a more natural fit and less fatigue. Ansell's exclusive formers are the result of a staff dedicated to the development and continuous improvement of formers.

**Stretchability (Modulus).** Because these gloves are latex, wearers experience a more secure fit that provides greater comfort without bagginess.

## Tactile Sensitivity

**Dependability.** Encore® gloves provide different balances of thickness, protection and sensitivity which cater to the complex needs of a diverse range of surgical procedures.

## Grip

**Application Driven.** Each glove in the Encore® line has a unique finish that is application driven for superior gripping capability.

	Encore®	Encore® Acclaim®	Encore® HydraSoft®	Encore® Underglove	Encore® MicrOptic®	Encore® Orthopaedic
<b>Material</b>	Natural rubber latex*	Natural rubber latex*	Natural rubber latex*	Natural rubber latex*	Natural rubber latex*	Natural rubber latex*
<b>Color</b>	Cream	Cream	Cream	Teal	Brown	Brown
<b>Donning Agent</b>	Polymer Coating	Polymer Coating	Polymer Coating	Polymer Coating	Polymer Coating	Polymer Coating
<b>Cuff Style</b>	SureFit Technology™	SureFit Technology™	SureFit Technology™	SureFit Technology™	SureFit Technology™	SureFit Technology™
<b>Finish</b>	Textured	Smooth	Smooth	Smooth	Smooth	Textured
<b>Minimum Length</b>	11.0 in (280 mm)	11.4 in (290 mm)	11.4 in (290 mm)	11.2 in (285 mm)	11.6 in (295 mm)	11.6 in (285 mm)
<b>Fingertip Thickness (single)</b>	9.0 mils 0.23 mm	8.5 mils 0.21 mm	8.5 mils 0.21 mm	7.1 mils 0.18 mm	7.1 mils 0.18 mm	13.0 mils 0.33 mm
<b>Cuff Thickness (single)</b>	7.0 mils 0.17 mm	6.2 mils 0.16 mm	6.2 mils 0.16 mm	5.0 mils 0.11 mm	5.0 mils 0.11 mm	9.6 mils 0.25 mm

CAUTION: This Product Contains Natural Rubber Latex, Which May Cause Allergic Reactions.

\*Caution: Safe use of these gloves by latex-sensitized individuals has not been established.

In addition to the overall benefits that Ansell Surgical gloves offer, each glove in the Encore® Family provides unique features for specific needs.

### IDEAL FOR GENERAL SURGERY-TEXTURED

Encore®

#### Addresses Powder-Related Problems

These gloves offer all the fit, tensile strength, protection and puncture resistance and tear resistance needed for most surgical procedures.

#### Unique features:

- Textured finish for reliable grip
- Recommended for the majority of procedures
- Produced on ergonomically designed formers to enhance comfort and fit



### IDEAL FOR GENERAL SURGERY-SMOOTH


Encore® Acclaim®

#### Smooth for Improved Sensitivity

These gloves offer all the fit, tensile strength, protection and puncture resistance and tear resistance needed for most surgical procedures.

#### Unique features:

- Anti-slip smooth finish provides worry-free instrument and suture control
- Recommended for the majority of procedures
- Produced on ergonomically designed formers to enhance comfort and fit



### WITH HYDRATING TECHNOLOGY

Encore® HydraSoft®


#### HydraSoft Technology™...for Hand Health

Ansell's moisturizing technology keeps hands hydrated to reduce dryness, roughness and cracked skin.

HydraSoft® utilizes two moisturizers: skin-friendly glycerine, which attracts and retains water, and dimethicone, which prevents the skin from drying out.

#### Unique features:

- Anti-slip smooth finish
- Recommended for the majority of procedures
- Produced on ergonomically designed formers to enhance comfort and fit
- Recommended to be used as an under glove because of the HydraSoft Technology™ incorporated in this product



### WITH HYDRATING TECHNOLOGY


Encore® Underglove

#### HydraSoft Technology™...for Hand Health

Ansell's moisturizing technology is designed to retain moisture and rehydrate your skin. It also helps to reduce dryness, roughness and cracked skin.

#### Unique features:

- Glove design allows for double gloving or single use
- Micro-thin design provides soft comfort and allows for finger tip sensitivity
- Recommended when there is a need for extra protection



### IDEAL FOR DELICATE PROCEDURES

Encore® MicrOptic®

#### Superb Tactile Sensitivity... Exceptional Strength

Encore® MicrOptic® gloves offer exceptional dexterity and tactile sensitivity while at the same time, outstanding strength and barrier protection. These gloves are ideal for delicate procedures in ophthalmic as well as micro- and cardiothoracic surgery.

#### Unique features:

- Anti-slip smooth finish
- Recommended for times when there is a need for more sensitivity to very low levels of tactile feedback
- Produced on ergonomically designed formers to enhance comfort and fit
- Micro-thickness allows for comfortable double-gloving



### ADDED PROTECTION


Encore® Orthopaedic

#### Extra Protection without the Powder

When we designed our Encore® Orthopaedic gloves, we made sure they provided the added strength and thickness for tear and puncture resistance that's vital to withstanding the rigors of high-risk, trauma and reconstructive procedures.

#### Unique features:

- Textured finish offers superior grip for handling instruments, sutures and situations when the application of high pressure is needed
- 30% thicker than our standard powder-free glove
- Recommended when there is a need for extra protection and puncture resistance



# Ordering Information

Encore®		Encore® Acclaim®		Encore® HydraSoft®		Encore® Underglove		Encore® MicroOptic®		Encore® Orthopaedic	
Reorder #	Size	Reorder #	Size	Reorder #	Size	Reorder #	Size	Reorder #	Size	Reorder #	Size
5785000	5 1/2	5795000	5 1/2	2018655	5 1/2	2018455	5 1/2	5787000	5 1/2		
5785001	6	5795001	6	2018660	6	2018460	6	5787001	6	5788001	6
5785002	6 1/2	5795002	6 1/2	2018665	6 1/2	2018465	6 1/2	5787002	6 1/2	5788002	6 1/2
5785003	7	5795003	7	2018670	7	2018470	7	5787003	7	5788003	7
5785004	7 1/2	5795004	7 1/2	2018675	7 1/2	2018475	7 1/2	5787004	7 1/2	5788004	7 1/2
5785005	8	5795005	8	2018680	8	2018480	8	5787005	8	5788005	8
5785006	8 1/2	5795006	8 1/2	2018685	8 1/2	2018485	8 1/2	5787006	8 1/2	5788006	8 1/2
5785007	9	5795007	9	2018690	9	2018490	9	5787007	9	5788007	9

1 Pair/Package - 50 Pairs/Dispenser Box - 4 Boxes/Case



**For information and samples, call Ansell at 800-952-9916 or visit us at [www.ansellhealthcare.com](http://www.ansellhealthcare.com)**

Ansell Healthcare Products LLC | [ansellhealthcare.com](http://ansellhealthcare.com)  
 200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)

Ansell, Encore, MicroOptic, Acclaim and Hydrasoft are trademarks owned by Ansell Limited or one of its affiliates.

©2011 Ansell Healthcare Products LLC. All Right Reserved.

SS200 5/2011

Printed on recycled paper using soy inks.



Everything you touch...we touch.

